

AMSCO Maintenance Manual



CLINICAL AUTOCLAVE

10x10x22" Series 7018

(12/87)

P-753817-091

INDEX

CONTENTS


	SHEET	Grid
Equipment Drawing	60232	A-3
Operating Instructions	83100-001	A-4
Maintenance Instructions	83100-001	A-6
Assembly Drawing — Sheet 1 of 3	SM-131	A-9
Assembly Drawing — Sheet 2 of 3	SM-131	A-11
Assembly Drawing — Sheet 3 of 3	SM-131	A-13
Door Assembly	SM-132	B-1
Measuring Chamber	SM-133	B-3
Operating Valve	SM-134	B-4
Principle of Operating Valve	SM-135	B-5
Pan/Water Assembly	30476-A	B-7
Timer Assembly (For Units Built Between Oct. 7, 1970 and Sept. 30, 1974)	134059-001	B-9
Timer Assembly (For Units Built Between Oct. 1, 1974 and May 31, 1983)	SM-136A	B-11
Timer Assembly (For Units Built After June 1, 1983)	136009-002	B-13
Timer Wiring Diagram (For Units Built Before July 31, 1965)	27474	C-2
Timer Wiring Diagram (For Units Built Between Aug. 1, 1965 and Sept. 30, 1974)	75091-091	C-3
Timer Wiring Diagram (For Units Built Between Oct. 1, 1974 and Feb. 6, 1975)	83030-001	C-4
Timer Wiring Diagram (For Units Built Between Feb. 7, 1975 and July 20, 1976)	83030-002	C-5
Timer Wiring Diagram (For Units Built Between July 21, 1976 and May 31, 1983)	83030-003	C-6
Timer Wiring Diagram (For Units Built After June 1, 1983)	83030-005	C-7
Timer/Sterilizer Schematic (For Units Built Between July 21, 1976 and May 31, 1983)	56194-001	C-9
Timer/Sterilizer Schematic (For Units Built After June 1, 1983)	56396-021	C-11
Sterilizer Wiring Diagram (For Units Built Before July 31, 1965)	30500	C-14
Sterilizer Wiring Diagram (For Units Built Between Aug. 1, 1965 and Sept. 30, 1974)	SM-137	D-1
Sterilizer Wiring Diagram (For Units Built Between Oct. 1, 1974 and May 31, 1983)	SM-137A	D-2
Sterilizer Wiring Diagram (For Units Built After June 1, 1983)	83040-003	D-3
Element Heating for 10x10x22 Clinical Autoclave	27331	D-4
Instructions for Thermostat Setting	SM-181	D-5
Instructions for Replacing Thermostat (For Units Built Between 10-74 and 4-76)	SM-181A	D-6

This document is intended for the exclusive use of AMSCO Service Technicians and customers. Reproduction in whole or in part is prohibited.


©AMSCO — 1979-1987

Printed in U.S.A.

REV. 12/87



PLAN VIEW



REAR VIEW

SIDE VIEW

NOTES:

1. 1/2" CONDUIT BOX. OTHER CONTRACTORS TO PROVIDE 240 VOLTS, 60 CYCLE, 20 AMP A.C. SERVICE AND MAKE CONNECTIONS. (FUSED FOR 4000 WATTS) FOR PROPER GROUNDING, 3 WIRE CONDUCTOR MUST BE USED FOR ELECTRICAL SERVICE.
2. AMSCO ASSUMES NO RESPONSIBILITY FOR CHANGES MADE NECESSARY THROUGH FAILURE TO OBSERVE THESE INSTRUCTIONS.
3. UNIT WEIGHT 280 LBS.
4. CLEARANCE INDICATED IS TO PROVIDE WORK SPACE FOR SERVICING WHERE ADJACENT EQUIPMENT IS PORTABLE, THIS DISTANCE MAY BE REDUCED.

AMERICAN STERILIZER CO., ERIE, PA.

10"x10"x22" CLINICAL MODEL
AUTOCLAVE
OPEN MOUNTED, ELECTRIC HEAT

DOOR SWING RH ☐ LH ☐

ITEM NO. _____ ROOM NO. _____

DWG. NO. B-60232

M100

REV. 0

SUPERSEDES F-24526

A-3



"AMERICAN"
10x10x22 CLINICAL AUTOCLAVE

OPERATING
INSTRUCTIONS

P-83100-001
(9/28/78)

WARNING
LIQUID STERILIZATION

To prevent possible personal injury or property damage resulting from bursting bottles and hot fluid, you must follow the recommended procedure listed below:

RECOMMENDED PROCEDURE:

- Use only vented closures — do not use screw caps or rubber stoppers with crimped seal.
- Use only Type 1 borosilicate (Pyrex) glass bottles — do not use ordinary glass jugs or any container not designated for sterilization.
- At end of cycle, open sterilizer door, no more than 1/2 inch. Wait 10 minutes before unloading sterilizer.
- Do not allow hot bottles to be jolted. This can cause hot-bottle explosions! Do not move bottles if any boiling or bubbling is present.
- Bottles should be cool to touch before attempting to move them from sterilizer loading car or shelves to the storage area.

OPERATING INSTRUCTIONS

Preliminary Preparation — Turn timer, thermostat and operating handle to "OFF" position. Pull sliding top panel forward until it stops and open lid on cover of reservoir. Pour distilled or softened water into reservoir until gauge on front of panel indicates "FULL" (Reservoir Capacity 12 Quarts). Allow 10 minutes for water to reach the measuring chamber.

Time required to bring the chamber to sterilizing temperature. The timer automatically engages only when the chamber is up to the selected temperature as indicated on the thermometer. Completion of the exposure period will be signalled by a buzzer, which will continue buzzing until unit is reset for the next cycle.

NOTE: CHECK THE WATER LEVEL EACH DAY.

To speed the first sterilizing cycle when starting with a cold autoclave, it is suggested that the following procedure be followed. Approximately 30 minutes before using autoclave turn thermostat knob to the desired temperature setting and push red reset button. The red signal light will glow, indicating that the current to the heaters is on. The autoclave door should be kept closed to conserve heat. The autoclave will continue to heat for about 30 to 35 minutes at this setting and will automatically turn off if the autoclave is not used during this period of time. When the red pilot light goes out, this indicates that the chamber has been sufficiently heated. The autoclave cannot be reset until it is used for a sterilizing cycle or until the chamber has cooled.

Timer — The knob may be turned clockwise from "OFF" position to any time setting. The total time involved in the sterilizing period, in any case, will be longer than the period for which the timer is set. This is due to the additional

RECOMMENDED EXPOSURE PERIOD

	250°F.	270°F.
Instruments without towel covers	15 Min.	3 Min.
Instruments with towel covers	20 Min.	7 Min.
Dressings and small packs in muslin covers	20 Min.	
Utensils, syringes and glassware, covered	15 Min.	3 Min.
Rubber gloves and tubing in muslin covers	15 Min.	
Solutions (125-200 ML Pyrex Flasks)	20 Min.	
Solutions (50 ML Pyrex Flasks, Test Tubes up to 100 MM)	15 Min.	

AMSCO | AMERICAN STERILIZER COMPANY • 2424 WEST 23rd STREET • ERIE • PENNSYLVANIA 16514

A-4

WARNING

TO PREVENT POSSIBLE PERSONAL INJURY RESULTING FROM BURSTING BOTTLES AND HOT FLUID, USE ONLY BOROSILICATE (PYREX) FLASKS WITH VENTED CLOSURES FOR STERILIZING LIQUIDS.

• SEE PAGE 4 FOR FURTHER INFORMATION

4 Operational Technique — Load Autoclave, close and tighten door, but not excessively. Turn "Operating Handle" to "Ster.", and recheck the thermostat knob setting to coincide with the above temperature chart and set timer knob at recommended exposure period. If the red signal is not on, press reset button to turn on current to heaters. No other mechanical operation is necessary until the end of the exposure period, since the unit now functions automatically.

Information — When the operating handle is turned to "Sterilize," the main heaters come on and a measured amount of water is released into the chamber well. Within a short period of time, steam is generated and the temperature begins to rise. The locked-in air settles to the bottom of the chamber and gradually is forced out through an air-eliminating system which is controlled by a thermostatic trap which locks only after all the air has passed from the chamber. When the thermometer reaches the desired temperature, this is the beginning of the exposure period. Should the thermometer reading indicate a temperature above or below the thermostat knob setting, this does not mean that the Autoclave is functioning improperly. Some variances may prevail, but always work from the thermometer reading, because the thermostat numerical setting is only an approximation.

This can be adjusted by turning the thermostat knob clockwise or counter-clockwise depending upon the required thermometer reading adjustment. If the thermostat knob has been turned to its full extent and thermometer does not reach 270° minimum, this is an indication that thermostat will have to be adjusted.

The unit is equipped with a safety valve which will protect the chamber from over-pressure conditions. When the safety valve opens, it is merely signalling the release of excess pressure from within the chamber. The Autoclave can still be operated with safety, until the arrival of your dealer's serviceman, by reducing the setting on the thermostat knob. In the event the Autoclave shuts off during a cycle, it means there is insufficient water in the chamber well. This can be corrected by turning the operating handle to "OFF" and allow 10 minutes for water from the reservoir to fill the measuring cylinder. Then turn the operating handle to "Sterilize" and press the red reset button to start the sterilizer.

5 To Exhaust and Dry Load — When the buzzer signals the end of exposure period, turn timer knob and operating handle to "OFF." A gurgling sound will indicate that the steam and unused water in the chamber is being exhausted into the reservoir. Leave door closed until thermometer indicates 212°F; then open door about 1/4" to allow cool air from room to enter chamber gradually. The heaters will stay on low heat after door is opened to accelerate drying. After 15 minutes turn off heaters by turning thermostat setting to zero, and remove trays from Autoclave. Close door (leave slightly ajar) so that radiated heat is contained within the chamber. This will benefit the following cycle.

6 Reloading — Wait at least 5 minutes after removing trays from previous cycle before starting the next one. After that, the process is repeated as outlined above.

ARRANGEMENT OF LOADS

- (a) **Instruments**: For routine sterilization of washed instruments use only trays with perforated bottoms as supplied with the Autoclave. Place a layer of muslin or towel in bottom of tray, arrange instruments in any desired order, and then cover instruments with a final layer of muslin or towel to facilitate drying and to prevent contamination in transit.
- (b) **Surgical Dressings, small packs and rubber gloves**: These supplies should always be wrapped in muslin, and placed on edge (never flat) in the Autoclave tray to assure ample steam circulation. Under no conditions permit crowding of packs into large or dense masses.
- (c) **Empty Glassware and Utensils** — Place all containers, items of glassware and utensils (wrapped or not) on their sides or in the inverted position in the Autoclave tray.
- (d) **Solutions**: Remove trays from chamber of Autoclave to accommodate flasks of 1000 ml or less capacity or racks of test tubes. Flasks of smaller capacity may be sterilized directly in the trays.

MAINTENANCE INSTRUCTIONS

1 The Door - To assure ease of operation and long life, place a small amount of high temperature grease periodically on screw threads of locking handle.

2 The Door Gasket - When door fails to close steam-tight under normal closing pressure, without straining handwheel, renew gasket. Gasket may be secured through the dealer from whom you purchased this autoclave.

When removing old gasket, scrape groove in door clean. Gasket is cut to a snug fit in groove and must be forced in, a short section at a time, without stretching. Should gasket appear to be too long, do not cut it but start over again, compressing short sections as inserted in groove, to take up full length.

3 Signal Lamp - To replace a burned out lamp, unscrew lamp lens, push in on lamp, and turn counter clockwise to

remove from socket. Insert lamp by placing in socket and turning clockwise. Pilot Lamp NE-51, Miniature Bayonet Type Base. Now screw lamp lens back into place.

4 Care of Chamber - The chamber is constructed with a removable inner liner. This liner must be removed at least once each month for inspection and cleaning of the outer chamber or steam jacket to remove lime and scale deposits on the heating surface. This is necessary to maintain maximum heating efficiency. To withdraw the liner from the Autoclave, remove the knurled strainer screw at the bottom center and lift out the drain tube fitting at the bottom rear of the liner.

For repairs or the replacement of parts not covered by the above routine maintenance suggestions, contact the dealer through which this Autoclave was purchased. Always mention the Serial and Model Number of the unit when writing in reference to or when ordering replacement parts.

REPLACEMENT PARTS

Name	Part No.	Req'd.
Light — Pilot	P-29517-091	1
Valve — Safety	P-150498-091	1
Lamp NE-51 (Box of 10)	P-764317-708	1
Thermometer	P-13638-091	1
Gasket — Heater	P-12854-091	2
Heater	P-27332-A-091	2
Thermostat	P-22446-091	1
Switch	P-32959-091	1
Gasket — Door	P-74375-091	1
Bellows — Measuring Chamber	P-12258-091	1
Gasket — Measuring Chamber	P-28198-091	1
Timer Assembly	P-136009-092	1

LIQUID STERILIZATION

Your AMSCO Sterilizer is designed to process liquids when borosilicate (Pyrex) flasks with vented closures are used.

Borosilicate (Pyrex) glass is recommended because it is a superior glass capable of containing higher pressures, of resisting thermal shock (such as cold air striking the hot glass), and of withstanding repeated handling.

Vented closures are recommended because, by design, they will prevent excess pressure by automatically venting a flask!

If other types of glass (such as flint glass) and non-venting (sealed) closures are used to sterilize liquids in your AMSCO Sterilizer, a potential dangerous condition, capable of causing personal injury and property damage, is created. As the liquid and residual air in a sealed flask are heated, they expand and create an internal pressure greater than the external pressure of the steam. With the weaker glass, a greater potential for bursting exists.

After the sterilization exposure, the chamber is exhausted slowly but it still exhausts more rapidly than the pressure within a sealed flask.

This pressure within the flask will exist until the residual air and the liquid have cooled (unlike a flask with a vented closure that prevents this excess pressure). Thus, the potential exists for the flask to burst and cause personal injury or property damage.

THIS PAGE
INTENTIONALLY LEFT
BLANK


LIMITATION OF LIABILITY AND INDEMNITY

In no event, whether as a result of breach of contract, warranty or tort (including negligence and strict liability) shall American Sterilizer Company or its suppliers be liable for any consequential or incidental damages including, but not limited to loss of profits or revenues, loss of use of the Products or any associated equipment, loss of the Buyer's Products, damage to associated equipment, cost of capital, cost of substitute products, facilities, service or replacement power, downtime cost, caused by such Products, or claims of the users for such damages. Buyer and ultimate user hereby agree to indemnify the American Sterilizer Company and to hold the American Sterilizer Company harmless from any and all liability for such consequential or incidental damages. The responsibility of the American Sterilizer Company for damages due to injuries or death or the death of employees of the Buyer or ultimate user of the Product, caused by the Product, shall be limited to that portion of such damages as might be attributable to the negligence or strict liability or other tortious conduct of the American Sterilizer Company. The Buyer and ultimate user agree to indemnify the American Sterilizer Company and hold the American Sterilizer Company harmless from any further damages, indemnity or contribution, if Buyer transfers title to or leases the Products sold hereunder to any third party. Buyer shall obtain from such third party a provision affording American Sterilizer Company and its suppliers the protection of this article relating to Limitations of Liability and Indemnities.

The American Sterilizer Company's liability for any claim of any kind (including negligence and strict liability) for any loss or damage arising out of, or resulting from this agreement, or from the performance or breach thereof, or from the Products or Services furnished hereunder, shall in no case exceed the price of the specified Product, system, component or service which gives rise to the claim. Except as to title, any such liability shall terminate one year from the date of installation of any Product or upon the expiration of the warranty period applicable to each type of Product covered hereby, whichever time period expires first.

*Note 1: Dial Plate, 31264, used on Units built before 10-1-74. Dial Plate, 134171, and Decal, 83073, (item 105) used on units built after 10-1-74.

A-9




NO		NAME		PT. NO. RD.		NO		NAME		PT. NO. RD.		NO		NAME		PT. NO. RD.	
96		NUT - HEX - JAM 5/8"		33961-045 4		48		NIPPLE		28905-091 1		6		1 SHELL & FRAME		27327-010 1	
97		SPACER		31334-042 2		49		SPACER		42284-091 6		2		2 DOOR ASSY.		39135-091 1	
98		PLATE - BOTTOM		31333-010 1		50		SCREW - TRUSS HD		10859-041 2		8		3 PANEL - SIDE ASSY		30469-010 2	
99		PLATE - BACK		31332-010 1		51		LOCKWASHER		19685-061 8		10		4 PAN - WATER ASSY		30476-NLA 1	
100		LABEL - CAUTION		50369-001 2		52		NUT - HEX. #10-32		2960-042 10		1		5 HEATER ASSY.		27332-091 2	
101		PLATE - DIAL *		134171-001 1		53		SWITCH		32959-091 1		2		6 CAP ASSY.		27338-042 1	
102		CONNECTOR TERMINAL		14166-091 1		54		COVER-CONDUIT-B.		22614-091 1		2		7 SCREW-RD. HD.		4003-051 2	
103		TURNBUCKLE		42266-061 1		55		SCREW - RD HD		9374-041 2		1		8 COVER-ASSEMBLY		33957-061 1	
104		SCREW - SELF TAP		23437-041 6		56		CHAMBER-MEAS.		30304-091 1		2		9 THERMOSTAT		22446-091 1	
105		DECAL - DIAL		83073-002 1		57		WASHER		27324-091 2		2		10 PANEL-LOWER ASSY		75090-061 1	
						58		BUSHING		27325-042 2		3		11 PANEL-UPPER		27357-061 1	
						59		LINER ASSY.		27300-061 1		3		12 BOARD-INSUL.		29526-091 2	
						60		TRAY ASSY.		27366-063 3		1		13 HOLDER-GLASS		22443-051 1	
						61		GAUGE-INDICATOR		27360-091 1		1		14 SPACER		27298-091 1	
						62		FLOAT ASSY.		30763-091 1		1		15 PLATE-HEATER		27333-061 2	
						63		PACKING-LINER		37274-002 1		1		16 PLATE-DIRECTION		41027-091 1	
						64		TUBE ASSY.		27320-042 1		1		17 PLATE-NAME		33114-091 1	
						65		PANEL - TOP FRONT		30466-010 1		1		18 TUBE - DRAIN		27336-042 1	
						66		SHELF		27361-061 1		2		19 COVER-HEATER		27341-010 1	
						67		CLAMP-HEATER		12823-091 2		1		20 TUBE-EXHAUST		30450-091 1	
						68		PLATE - NAME		33944-091 1		1		21 MOLDING 42" GMC		3000-965 1	
						69		BOX - TERMINAL ASSY		31331-091 1		1		22 NIPPLE		29302-091 1	
						70		SCREW - OVAL HD.		18976-048 1		1		23 LIGHT-PILOT		29517-091 1	
						71		TUBE - FILL		30470-091 1		1		24 VALVE-SAFETY		763710-001 1	
						72		TUBE		30471-091 1		1		25 LAMP - NE-51		23883-091 1	
						73		PLATE-NAME		27043-091 1		1		26 ELL - CONDUIT 90°		25388-091 1	
						74		GASKET		22448-091 1		2		27 PANEL-TOP REAR		43416-010 1	
						75		JUMPER-LINK		14278-091 1		2		28 SCREW-BUTTON		37344-048 20	
						76		GASKET-HEATER		12854-091 2		1		29 ELBOW		1635-051 1	
						77		NUT - HEX. - JAM		15359-042 1		1		30 THERMOMETER		13638-091 1	
						78		LOCKWASHER 5/8"		19692-061 1		2		31 SCREW-BINDER		20823-061 4	
						79		SCREW-SOCKET		15339-045 2		1		32 SCREW-RETAINER		13841-062 1	
						80		FITTING-COMPR.		6695-044 1		1		33 STUD		13198-061 1	
						81		ROD		27329-061 1		2		34 SCREW - RD. HD.		9315-041 2	
						82		SCREW-DRIVE		13289-045 2		2		35 FITTING-COMPR.		6750-044 1	
						83		LOCKWASHER 3/8"		19687-061 2		3		36 SCREW-SET		34518-061 3	
						84		HANDLE		33137-091 1		2		37 SCREW-RD. HD.		12451-041 2	
						85		SCREW-RD. HD.		33742-041 2		1		38 NIPPLE 1/2"x1 3/4"		29165-051 1	
						86		NIPPLE		29334-091 1		6		39 LOCKWASHER 5/16"		19691-061 6	
						87		KNOB		14378-091 1		4		40 ELL		1633-051 1	
						88		SCREW-TRUSS HD		8129-042 4		1		41 COCK-DRAIN ASSY		18543-051 1	
						89		RECEPTACLE		75077-091 1		4		42 ELBOW		1614-091 1	
						90		LOCKWASHER #6		19675-041 4		1		43 FLANGE-FLOOR		27570-045 4	
						91		CHASSIS-TIMER		136009-001 1		2		44 RING-RETAINING		30580-091 1	
						92		CLAMP		33943-061 2		6		45 ELBOW 90°		25389-091 1	
						93		SCREW-HEX. CAP		12552-061 6		1		46 STOP		33939-091 1	
						94		LID		33940-061 1		4		47 HANDLE		30460-051 1	
						95		NUT - HEX. #6-32		3037-041 4							


SEE SHEET 2 OF 3
SM-131 FOR
SECTION C-C

A-10

10 X 10 X 22 CLINICAL AUTOCLAVE


REAR





A-13

A-14



BILL OF MATERIAL			
NO.	NAME	PART NO.	QTY.
1	PLATE, DOOR	39136-034	1
2	RING, DOOR POST	39132-034	1
3	HUB	13195-057	1
4	HANDLE	14108-091	5
5	OILER	14178-091	1
6	GASKET	74375-091	1
7	NIPPLE, THERM.	41552-056	1
8	PIN, GROOVE	29848-045	1
9	PLATE	150312-001	1
10	POST, DOOR	13196-091	1
11	BEARING (4 PIECE SET)	13700-091	1
12	SCREW	11249-041	4
13	PIN, LOCATING	39133-061	2
14	BAR, LOCK R.H. BAR, LOCK — L.H.	13698-056 14665-056	1 1
15	SCREW, SOC. HD.	12176-041	3
16	LOCKWASHER	10436-041	3
17	NEEDLE THRUST BEARING	P-48267-091	1
18	THRUST BEARING RACE	P-78954-091	2

NOTE 1


REMOVE OLD GASKET, SCRAPE GROOVE CLEAN. GASKET IS CUT TO A CLOSE FIT IN GROOVE AND MUST BE FORCED IN, A SHORT SECTION AT A TIME, WITHOUT STRETCHING. SHOULD GASKET APPEAR TO BE TOO LONG, DO NOT CUT IT, BUT START OVER AGAIN, COMpressing SHORT SECTIONS AS INSERTED IN GROOVE, TO TAKE UP FULL LENGTH.

NOTE 2

AT INTERVALS OF THREE MONTHS OR SO PUT A SMALL AMOUNT OF HIGH TEMPERATURE GREASE ON THREAD OF STUD ON WHICH DOOR CENTERS AND A FEW DROPS OF OIL IN OIL HOLE BACK OF HANDWHEEL.

DOOR ASSEMBLY

BILL OF MATERIAL			
NO	NAME	PT. NO.	R'D
1	CHAMBER	27346-091	1
2	COVER	27304-091	1
3	BELLOWS	12258-091	1
4	GASKET	28198-091	1
5	LOCKWASHER	19678-045	4
6	SCREW	3999-051	4
7	FITTING	6750-091	2
8	PLUG	3439-051	1




MEASURING CHAMBER
P-30304-091

REV. 8/83

SM - 133

B-3

BILL OF MATERIAL			
NO	NAME	PT. NO.	R'D
1	BASE	14232-091	1
2	GASKET	14234-091	1
3	BONNET	14233-091	1
4	SCREW- FIL.	14257-061	3
5	PACKING	46355-091	2
6	NUT-PACK'G	14236-051	1
7	GLAND	7755-042	1
8	STEM ASSY.	21281-061	1
9	SPRING	9280-061	1
10	WASHER	10414-042	1
11	SEAT- VALVE	20512-091	1




OPERATING VALVE ASSEMBLY P-30541-091

REV. 8/83

SM-134

B-4




PRINCIPLE OF OPERATING VALVE

SM 135




THIS PAGE
INTENTIONALLY LEFT
BLANK

BILL OF MATERIAL			
NO	NAME	PT. NO.	R'O.
1	WATER PAN	27350-A	1
2	COIL - CONDENSER	30499-A	1
3	OPERATING VALVE	30541-091	1
4	GASKET	14245-091	1
5	SCREW - RD. HD.	9314	3
6	FITTING - COMPRESSION	19514-091	1
7	FITTING - COMPRESSION	6695-091	2




PAN-WATER ASSEMBLY 30476-A

Add. 8/83



BILL OF MATERIAL			
NO.	PT. NO.	NAME	R'Q.
	134059-001	CHASSIS (TIMER) ASSEMBLY	X
1	25842-091	SOCKET, OCTAL	1
2	8129-042	SCREW, TRUSS. HD. (6-32 x 3/8)	1
3	3037-041	NUT, HEX (6-32)	3
4	26720-091	CLAMP, TUBE (1/4 ALUM)	1
5	25839-091	RELAY, THERMAL	1
6	75086-063	CASE ASSEMBLY	1
7	24539-091	RING, RETAINING	1
8	55953-001	TIMER	1
9	20824-045	SCREW, RD. HD. (6-32 x 2-1/4)	2
10	20340-091	FUSE HOLDER	1
11	17796-091	WASHER, BRASS	2
12	19712-091	WASHER, MICA	8
13	18159-091	RESISTOR (2500 Ω -10 WATT)	2
14	75082-091	BUSHING	1
15	75089-091	WIRE AND PLUG ASSEMBLY	1
16	25853-045	NUT, HEX (3/8-32)	1
17	25844-091	BUSHING	1
18	26710-091	WIRE, BLACK	2
19	25847-061	COVER (NOT SHOWN)	1
20	12800-042	SCREW (6-32 x 3/16) (NOT SHOWN)	2
21	75091-091	WIRING DIAGRAM	1
22	27429-091	SCREW, SET (8-32 x 3/16)	1
23	33810-091	FUSE	1
24	77798-091	BUSHING, SNAP	1
25	82674-001	DECAL, CAUTION	1
26	82645-001	DECAL, WARNING	1

TIMER ASSEMBLY
(For Units Built Between Oct. 7, 1970 and Sept. 30, 1974)
134059-001
Rev. 4/85






*NOTE 1: Place caution label on the outside of cover (20).

‡NOTE 2: Place wiring diagram on the inside of cover (20).

‡‡NOTE 3: No. 18 stranded copper wire, thermoplastic insulation, 16/30 str., 300 V min., 105 C min., UL label required.

**NOTE 4: 92643-001 is a new plug assembly with a round post in the center. The part number of the female receptacle for the plug is P-75077-091.




BILL OF MATERIAL			
NO	PT. NO.	NAME	R/O
1	136009-001	TIMER ASSEMBLY	X
2	136010-001	CHASSIS ASSEMBLY	1
3	83035-001	TIMER WITH KNOB, CRAMER 271 — 60 MIN.	1
4	80518-091	RELAY, TIME DELAY (180 SEC.)	1
5	79739-001	BUZZER	1
6	12529-061	SCREW, R.H. (6-32 X 1/4)	2
7	41295-091	RELAY	1
8	83036-001	RESISTOR (500 Ω -50 WATT)	1
9	83037-001	RESISTOR (1500 Ω -25 WATT)	1
10	83041-001	BRACKET, MTG.	4
11	82674-001	LABEL, CAUTION (NOT SHOWN)*	1
12	3037-041	NUT, HEX (6-32)	6
13	32118-091	ADAPTOR, TERM.	2
14	24539-091	RING, RETAINING	1
15	75082-091	BUSHING	1
16	92643-001	PLUG ASSEMBLY**	1
17	12531-061	SCREW R.H. (6-32 X 3/8)	8
18	19675-041	WASHER, LOCK (#6)	8
19	83030-001	STICKER, W.D. For Units Built Between 10/74 and 3/75 (NOT SHOWN)‡	1
	83030-002	STICKER, W.D. For Units Built Between 3/75 and 4/76 (NOT SHOWN)‡	1
	83030-003	STICKER, W.D. For Units Built After 4/76 (NOT SHOWN)‡	1
20	25847-061	COVER	1
21		WIRE, BLACK (87 IN. LG)‡‡	1
22	81400-001	TERMINAL Q.C.	7
23	77299-091	WIRE TIES (NOT SHOWN)	3

TIMER ASSEMBLY (For Units Built Between

Oct. 1, 1974 and May 31, 1983)

REV. 8/83

SM-136A



BILL OF MATERIAL			
NO.	PT. NO.	NAME	R'Q.
1	136009-002	TIMER ASSEMBLY	X
2	136010-002	CHASSIS ASSEMBLY	1
3	83035-001	TIMER, CRAMER 271-60 MIN.	1
4	150822-044	RELAY, TIME DELAY (180 SEC.)	1
5	79739-001	BUZZER ASSEMBLY	1
6	12800-042	SCREW (6-32 x 1/4)	2
7	41295-091	RELAY	1
8	150822-045	RESISTOR (3900 Ω -12 WATT)	2
9	150822-046	RESISTOR (4700 Ω -12 WATT)	2
10	150822-157	SOCKET	1
11	82674-001	LABEL, CAUTION (NOT SHOWN)*	1
12	3037-041	NUT, HEX (6-32)	8
13	32118-091	ADAPTOR, TERM.	1
14	24539-091	RING, RETAINING	1
15	75082-091	BUSHING	1
16	56396-019	PLUG ASSEMBLY	1
17	12531-061	SCREW, RD. HD. (6-32 x 3/8)	10
18	19675-041	WASHER, LOCK (#6)	10
19	83030-005	STICKER, W.D. (NOT SHOWN)**	1
20	25847-061	COVER	1
21		WIRE, BLACK (#3, 99 IN. LG.)†	1
22	81400-001	TERMINAL, Q.C.	5
23	77299-091	WIRE TIES	3
24	78629-091	SHRINK TUBE††	8
25	3038-041	NUT, HEX (8-32)	2
26	11241-041	SCREW, RD. HD. (8-32 x 5/8)	2
27	19676-041	WASHER, LOCK (#8)	2
28	118177-091	TERMINAL, SPADE	8

*NOTE 1: Place caution label (11) on the outside of cover (20).

**NOTE 2: Place wiring diagram sticker (19) on the inside of cover (20).

†NOTE 3: Wire specs: #18 stranded copper wire, thermoplastic insulation, 16/30 str., 300 V. min., 105°C min., UL, CSA label required.


††NOTE 4: 3/16" I.D. black polyolefin tubing. Use 1/2" long piece of tubing for each wire at all solder connections at resistors, items 8 and 9.

TIMER ASSEMBLY
(For Units Built After June 1, 1983)


136009-002

ADD. 8/83


THIS PAGE
INTENTIONALLY LEFT
BLANK




TIMER WIRING DIAGRAM
(For Units Built Before July 31, 1965)



TIMER WIRING DIAGRAM
(For Units Built Between Feb. 7, 1975 & July 20, 1976)



TIMER WIRING DIAGRAM
(For Units Built Between July 21, 1976 & May 31, 1983)




TIMER WIRING DIAGRAM
(For Units Built After June 1, 1983)

ADD. 8/83

83030-005


THIS PAGE
INTENTIONALLY LEFT
BLANK







56396-021

THIS PAGE
INTENTIONALLY LEFT
BLANK







D-1



D-2




STERILIZER WIRING DIAGRAM
(For Units Built After June 1, 1983)

ADD. 8/83

83040-003

D-3



NOTE : COLD RESISTANCE - LOW 23.7 OHMS,
NOMINAL - 26.3 OHMS, HIGH 27.6 OHMS


ELEMENT HEATING
10 X 10 22 CLINICAL AUTOCLAVE

ADD. 12/87

27331

D-4

INSTRUCTIONS FOR THERMOSTAT SETTING



DESCRIPTION

1. Attach test lamp

2. Turn control knob clockwise to extreme high position. Remove knob and operate autoclave until temperature reaches 272°F. If temperature continues to rise, loosen clamp screws and hold post firmly against stop, then turn graduated dial in direction of arrow L until test light goes on. **

3. If temperature does not reach 272°F, test light goes on repeat above procedure, but turn dial in direction of arrow R to increase temperature. **

4. Remove test light and replace knob.

** **NOTE:** This adjusting procedure should be repeated for at least 3 cycles - moving dial as necessary to obtain proper setting.

SM-181

D-5

INSTRUCTIONS FOR REPLACING THERMOSTAT (For Units Built Between 6-76 and 4-15-76)

When replacing thermostat, item 9, SM-131, proceed as follows:

NOTE: The old thermostat bulb was secured by one clamp; therefore, an additional clamp, P-33943-061, will be required for installing the new thermostat bulb.

1. Disconnect electric power to sterilizer.
2. Remove button screws (28, SM-131) securing side panels (3, SM-131), then remove panels.
3. Identify wires connected to thermostat (9, SM-131), then disconnect wires.
4. Remove thermostat bulb from clamp (92, SM-131) that secures it near heater assembly (5, SM-131).
5. Unscrew plug on thermostat capillary tube from stuffing box. Withdraw thermostat bulb from sterilizer.
6. Remove knob from thermostat, then remove two round head screws (37, SM-131) that secure thermostat to upper panel (11, SM-131).
7. Remove thermostat including attached bulb.
8. Loosen two 6-32 x 1/4 round head screws (6, SM-136A) that secure cover (20, SM-136A) to timer chassis assembly (2, SM-136A). Also refer to item 19, SM-131.
9. (Refer to 136009 and 136009A.) Remove wire number 9 that is soldered to TD1-8 and TD1-7.


NOTE: For steps 10 and 11, use number 8 stranded copper wire that meets the requirements in note 3 of SM-136A.

10. (Refer to 136009 and 136009A.) Add a 9-inch long black wire (number 3) by soldering one end to CR1-C (common) and soldering the other end to R1-A.
11. (Refer to 136009 and 136009A.) Add a 9-inch long black wire (number 11) by soldering one end to CR1-C (normally open) and soldering the other end to TD1-7.
12. Remove timer wiring diagram sticker (19, SM-136A) from cover (20, SM-136A), and replace it with new timer wiring diagram, P-83030-003.
13. Install new thermostat, P-22446-091 (9, SM-131) in upper panel (11, SM-131) using the two round head screws (37, SM-131).
14. Connect wires to new thermostat.
15. Insert thermostat bulb through hole in stuffing box, then screw plug into stuffing box.
16. Install thermostat bulb as shown using two clamps.
17. Replace timer assembly cover (20, SM-136A) and tighten round head screws (6, SM-136A).
18. Secure side panels (3, SM-131) to sterilizer with button screws (28, SM-131).
19. Reestablish electric power to sterilizer.

Sheet 1 of 2
REV. 2/79

SM-181A

D-6



*Note: Bulb must
not touch heater.



**AMSCO
SERVICE**

**CLINICAL AUTOCLAVE
SERIES 7018 10x10x22"
P-753817-002**

12/87

1 of 1

