

AMSCO  
Maintenance  
Manual

EO DISPOSER™

(9/30/92)

P-764323-253

## Table of Contents

|                                                                                                                        | Grid       |
|------------------------------------------------------------------------------------------------------------------------|------------|
| <b>SECTION 1: GENERAL INFORMATION</b>                                                                                  | <b>1-1</b> |
| 1.1 Application and Design                                                                                             | 1A-10      |
| 1.2 Major Components—Disposer System                                                                                   | 1A-10      |
| 1.2.1 Sequencer Kit                                                                                                    | 1A-10      |
| 1.2.2 Additional Materials                                                                                             | 1A-10      |
| <b>SECTION 2: OPERATING INSTRUCTIONS</b>                                                                               | <b>2-1</b> |
| 2.1 Sterilizer/Disposer Communications                                                                                 | 1A-12      |
| 2.2 General Operation of Sterilizer with EO Disposer                                                                   | 1A-12      |
| 2.2.1 2-1                                                                                                              |            |
| 2.2.2 2-1                                                                                                              |            |
| 2.3 Sterilizer Operation with Disposer                                                                                 | 1A-14      |
| 2.3.1 Disposer Operation with Eagle 3000 Series Sterilizers                                                            | 1A-14      |
| 2.3.2 Disposer Operation with Eagle 2000 Printcon Sterilizers                                                          | 1B-1       |
| 2.3.3 Disposer Operation with Eagle 2000 IRC Sterilizer                                                                | 1B-3       |
| 2.3.4 Disposer Operation with Medallion Sterilizers                                                                    | 1B-4       |
| 2.4 Sequencer — General                                                                                                | 1B-6       |
| 2.4.1 Indicator Lights — Sequencer                                                                                     | 1B-6       |
| 2.5 Aerators — General                                                                                                 | 1B-7       |
| 2.6 Disposer                                                                                                           | 1B-7       |
| <b>SECTION 3 — PRINCIPLES OF OPERATION</b>                                                                             | <b>3-1</b> |
| 3.1 General                                                                                                            | 1B-10      |
| 3.2 Cycle Description — Disposer                                                                                       | 1B-10      |
| 3.2.1 Disposer READY Mode Drop-out                                                                                     | 1B-11      |
| 3.2.2 Disposer Outputs                                                                                                 | 1B-11      |
| 3.2.3 Other Disposer Communications                                                                                    | 1B-11      |
| 3.3 Cycle Description — Eagle 3000 Sterilizer Operation with Disposer                                                  | 1B-11      |
| 3.3.1 General                                                                                                          | 1B-11      |
| 3.3.2 Disposer On Signal                                                                                               | 1B-11      |
| 3.3.3 Exhaust phase                                                                                                    | 1B-12      |
| 3.3.4 Air Break                                                                                                        | 1B-12      |
| 3.3.5 Second post vac and air break                                                                                    | 1B-12      |
| 3.3.6 Air flush                                                                                                        | 1B-12      |
| 3.3.7 Water dump                                                                                                       | 1B-12      |
| 3.3.8 Aeration                                                                                                         | 1B-12      |
| 3.4 Cycle Description - Eagle 2000 Printcon Sterilizer and Eagle 2000 Indicator Round Chart Recorder (IRC) Sterilizers | 1B-12      |
| 3.4.1 Cycle start                                                                                                      | 1B-12      |
| 3.4.2 Exhaust                                                                                                          | 1B-12      |
| 3.4.3 First Post Vacuum and Air Break                                                                                  | 1B-13      |
| 3.4.4 Second Post Vacuum and Air Break                                                                                 | 1B-13      |
| 3.4.5 Air Flush                                                                                                        | 1B-13      |
| 3.4.6 Water Dump                                                                                                       | 1B-13      |
| 3.4.7 Aeration                                                                                                         | 1B-14      |

|                                                                            | Grid            |
|----------------------------------------------------------------------------|-----------------|
| 3.5 Cycle Description: Medallion Sterilizers 16X16, 20X20, and 24X36 ..... | 3-5 1B-14       |
| 3.5.1 General .....                                                        | 3-5 1B-14       |
| 3.5.2 Cycle start .....                                                    | 3-5 1B-14       |
| 3.5.3 Exhaust .....                                                        | 3-5 1B-14       |
| 3.5.4 Air break .....                                                      | 3-6 1C-1        |
| 3.5.5 Second post vacuum and air break .....                               | 3-6 1C-1        |
| 3.5.6 Air flush .....                                                      | 3-6 1C-1        |
| 3.5.7 Water dump .....                                                     | 3-6 1C-1        |
| 3.5.8 Aerate .....                                                         | 3-6 1C-1        |
| 3.5.9 Complete .....                                                       | 3-6 1C-1        |
| 3.5.10 Operator Unloading .....                                            | 3-6 1C-1        |
| 3.5.11 Short cycle .....                                                   | 3-6 1C-1        |
| 3.6 Cycle Description — Sequencer Control .....                            | 3-6 1C-1        |
| 3.6.1 General .....                                                        | 3-6 1C-1        |
| 3.6.2 Cycle Operation .....                                                | 3-7 1C-2        |
| 3.6.3 Alarms and abnormal conditions .....                                 | 3-7 1C-2        |
| 3.7 Cycle Description — Aerators .....                                     | 3-7 1C-2        |
| <b>SECTION 4: GENERAL INSPECTION AND MAINTENANCE .....</b>                 | <b>4-1 1D-6</b> |
| 4.1 General Inspection and Maintenance .....                               | 4-1 1D-6        |
| 4.2 Recommended Preventive Maintenance - AMSCO EO Disposer System .....    | 4-2 1D-7        |
| <b>SECTION 5: TROUBLESHOOTING .....</b>                                    | <b>5-1 1D-8</b> |
| 5.1 Troubleshooting Chart .....                                            | 5-1 1D-8        |
| <b>SECTION 6: COMPONENT REPLACEMENT .....</b>                              | <b>6-1 1D-8</b> |
| 6.1 Adjustment Procedures .....                                            | 6-1 1D-8        |
| 6.1.1 Sterilizer Components .....                                          | 6-1 1D-8        |
| 6.1.2 Disposer Components .....                                            | 6-1 1D-8        |
| 6.2 Component Replacement Procedures .....                                 | 6-1 1D-8        |
| 6.2.1 Prefilter .....                                                      | 6-1 1D-8        |
| 6.2.2 Catalytic Cell .....                                                 | 6-1 1D-8        |
| 6.2.3 Power Switch Cable Adjustment .....                                  | 6-2 1E-1        |

## SUMMARY OF WARNINGS AND CAUTIONS

Following is a summary of the safety warnings and cautions which must be observed when operating or servicing this equipment. WARNINGS indicate the potential for danger to personnel and CAUTIONS indicate the potential for damage to equipment. These warnings and cautions are repeated, where applicable, throughout the manual.

### WARNINGS:

- ▲ **WARNING** — EO is released from the drain funnel and from the door by this procedure. Avoid the area around the sterilizer. Insure the Envirogard fan is operational. Use of protective equipment (respirator, eye protection, gloves, jumpsuit) may be necessary.
- ▲ **WARNING** — Activating the try lever will require that the access room behind the sterilizer be entered. Protective equipment must be worn. EO will be present in the access room. Do not perform this procedure unless the protective gear is capable of protecting against expected EO levels.

### CAUTIONS:

- ▲ **CAUTION** — Actuation of safety valve in the manner described above may cause the valve to leak when the sterilizer is put back into service. It will be necessary to leak check the safety valve prior to turning the unit back over to the customer.

THIS PAGE  
INTENTIONALLY LEFT  
BLANK

# AMSCO EO Disposer System

## DESCRIPTION

The EO Disposer System converts ethylene oxide gas into carbon dioxide and water vapor by means of catalytic combustion. An efficiency rate of 99.9% minimizes environmental EO with correct installation.

The system is designed to interface with AMSCO Eagle and Medallion Series gas sterilizers and aerators to remove EO from the equipment exhaust. The Disposer makes any of these straight gas sterilizers function as a sterilizer/aerator. The system consists of a Disposer and an Interconnect to the sterilizer and/or aerator. In multiple-sterilizer installations, a Sequencer allows the Disposer to handle up to three sterilizers and three aerators (maximum 110 cfm) by controlling exhaust flow through the system. (Note: Only two small aerators may be connected to the disposer.)

The Disposer assembly consists of a catalytic cell, pre-filter, vacuum pump, heater and blower. The Disposer uses flameless oxidation to convert EO.

Interconnect kits allow the Disposer to be linked to any model Eagle Series or Medallion Series straight gas sterilizer or aerator. The Interconnect kit includes a printed circuit board which allows the operator to use the Disposer system simply by initiating normal sterilization or aeration procedures. No further operator involvement is needed.

The Sequencer control, if specified, allows one sterilizer at a time to vent

to the Disposer, with up to three small or medium sterilizers on the system. The use of a sequencer will lengthen cycle times and productivity may be affected. If aerators are also on the system, they will vent continuously to the Disposer through a separate connection.

## STANDARDS

Meets applicable requirements of the following listings and standards and carries the appropriate symbols:

ETL Testing Laboratories, Inc.

Canadian Standards Association (CSA)


Seismic Stress Calculation—Cal. Admin. Code—Title 24.

## SYSTEM COMPONENTS

Disposer cabinet houses:

- **Vacuum pump** to provide suction on sterilizer chamber during exhaust and aeration
- **Prefilter** to remove particulate matter from the airstream
- **Blower** to feed air to system
- **Heater** to raise airstream and catalyst temperature to catalytic reaction temperature (280° F [138° C] minimum)
- **Catalytic cell** to remove EO from the airstream through an exothermic (heat-producing) reaction.

**Interconnect Kits**—Designed to be field-installed on Eagle or Medallion sterilizers. Kit includes piping



Typical only—some details may vary

THE SELECTIONS BELOW APPLY TO THIS EQUIPMENT

### INTERCONNECT KIT (ONE PER STERILIZER OR AERATOR):

- ☐ Eagle 3000 EO Sterilizer
  - Small
  - Medium
- ☐ Eagle 2000 EO Sterilizer with Printcon
  - Small
  - Medium
- ☐ Eagle 2000 EO Sterilizer with Round Chart Recorder
  - Small
  - Medium
- ☐ Medallion EO Sterilizer
  - Small
  - Medium
- ☐ Aerator Kit

### OPTIONS:

- ☐ Sequencer (if multiple sterilizers)

Item \_\_\_\_\_

Location(s) \_\_\_\_\_

modifications and a box with control boards, program chips, and wiring. Aerator interconnect kit is to be field-installed on Eagle or Medallion aerators.

**Sequencer** (optional)—Sends and receives signals to control access to Disposer system in multi-sterilizer installations.

## TECHNICAL DATA

**Control**—System control is designed as part of the Interconnect kit so that the operator can use the Disposer simply by initiating normal sterilization or aeration procedures. No further operator involvement is needed.

**Catalyst**—The catalyst is proprietary, a granular material composed of manganese dioxide and copper oxide housed in a cell made of type 304 stainless steel approximately 20x20x6 inches (508x508x152 mm) with a 4-inch (102 mm) bed depth. The catalyst removes EO from the airstream through an exothermic (heat-producing) reaction. The amount of temperature increase within the catalytic cell depends upon the concentration of ethylene oxide in the airstream, the flow rate, heat loss and amount of catalyst. The expected life of the catalyst is one to three years, dependent on usage. Spent catalyst is not a hazardous material, but must be disposed of according to local, state and federal regulations.

**Vacuum Pump**—0.5 hp, nominal 25 inch Hg.

**Prefilter**—Resistance 1 inch w.g. at 110 scfm.

**Blower**—Single speed. 0.6 horsepower blower feeds air to system at rate of 110 scfm.

**Heater**—10 Kw electric, finned tube duct heater heats airstream and catalyst to approximately 280° F (138° C) minimum, the temperature at which catalytic reaction is 99.9% efficient.

**Exhaust**—Exhaust port from Disposer unit is 5 inch diameter with 5 inch to 6 inch adaptor. Duct for venting to atmosphere is to be provided by others. Idling exhaust temperature is 280° F (138° C). Normal operating exhaust temperature is 450° F (233° C); maximum exhaust temperature is 500° F (260° C). If overtemperature is reached, system will prevent introduction of additional EO to the catalyst until proper operating temperature resumes.

## MATERIAL

**Disposer Cabinet**—NEMA 3 Enclosure; Painted Carbon Steel

**Plenum**—Aluminized Steel and Painted Carbon Steel

**Catalyst Cell**—Stainless Steel Type 304

## MOUNTING

**Electricity**—All internal wiring within the Disposer cabinet has been completed at the factory. Customer is responsible for providing 208 V, 60 Hz, 3 Phase, 40 Amp power to the main disconnect located within the Disposer cabinet.

**Piping**—Piping modifications to existing sterilizers or aerator are provided as part of the Interconnect kit, which is designed to fit specified models. Piping from the sterilizer to the Disposer and exhaust piping from the Disposer are to be supplied by others.

**Installation**—Disposer has lifting lugs and can be moved with a forklift. Cabinet can be tilted or laid on side (not on front or back) during installation, and meets seismic requirements when properly bolted to building floor or roof.

## CYCLES, READOUTS AND ALARMS

All indicators, controls and safety devices are included to assure proper system operation and safety. If used with Printcon and Eagle 3000 sterilizers, the printer tape reflects the elapsed Disposer exhaust and aeration time as part of the cycle summary. Displays for Eagle 3000 and Printcon sterilizers show countdown aeration time remaining. Disposer installation with Printcon and Eagle 3000 sterilizers will add printed alarm messages indicating Disposer readiness, air flow, exhaust and catalyst status. For Eagle Round Chart and Medallion units, a remote alarm box provides this information. The separate interface box (all units) has LED readouts indicating cycle phases for service personnel.

## WARRANTY

AMSCO warrants that each Disposer is carefully tested, inspected and leaves the factory in proper working condition, free of visible defects. Standard coverage includes one (1) year on parts (except expendables) and labor from date of installation. AMSCO representatives can provide full details of the warranty program upon request.

## OPTIONAL PREVENTIVE MAINTENANCE AGREEMENT


A coast-to-coast network of skilled and competent specialists can provide periodic PMA inspection and adjustment to assure low-cost peak performance.

## ENGINEERING DATA

**Size:** Nominal 29" wide x 29" deep x 82" high (737x737x2083mm)  
**Disposer Weight:** Nominal 830 lbs (376 kg)  
**Air Flow:** 110 scfm  
**Processing Capability:** 0.034 lbs/min of EO in a continuous flow of 12% ethylene oxide 88% refrigerant mixture

### Temperature Specifications:

- Ambient storage/transport: Minus 20° F to 120° F (Minus 29° C to 49° C)
- Ambient Operating Temp. Range: 50° F min - 120° F max (10° C min - 304° C max)
- Operating temperature:
  - Catalyst: 280° F to 450° F (138° C to 233° C)
  - Exhaust: 450° F (233° C)
  - Overtemperature: 500° F (260° C)
- Heat Loss to Surroundings: 3000 BTU/hr at 70° F (21° C) (Disposer only; no exhaust pipe considered)



DIMENSIONS ARE INCHES (MILLIMETERS)—DRAWING IS NOT TO SCALE

This print is for guidance when planning space and utility services. Actual installation prints may be obtained from any AMSCO office representative.

## OPERATING REQUIREMENTS

- Ethylene Oxide inlet—3/4" NPT
- Exhaust 7-1/8" dia. flange (removable), 5" dia. OD outlet Temperature Range, Ambient: 500° F
- Electrical—208 VAC, 3 Phase, 60 Hz, 40 Amp
- Interface with sterilizer 22/24 AWG, 24 VAC, 12 conductor cable (not by AMSCO). House in 1" electrical conduit.
- Optional Catalyst Bed Temperature Recorder—Type J thermocouple connector
- Installation outside the facility must be limited to climates 50° F (10° C) and above. Otherwise disposer

must be enclosed in an insulated, heated structure and intake air must be supplied at 50° F (10° C) and above.

- Aerator inlet 4" ODT.
- Piping from sterilizer(s) to disposer 200' maximum. Minimum 5' vertical rise from 3/4" NPT inlet (see 1). Note: The 200' maximum is affected by vertical and horizontal runs and must also reflect losses due to pipe elbows and other fittings. See Equipment Drawing P-65435-332 for installation details.
- For installation details order Installation Guidelines Booklet P-129362-980




...CHECK LOCAL CODES...

## SPECIFICATION WORKSHEET

ITEM \_\_\_\_\_: Furnish an AMSCO EO Disposer System per Tech Data Sheet SD-433. Furnish Disposer Assembly and Interconnect Kit for: \_\_\_\_\_ small Eagle 3000 sterilizer with Stage 3 control; \_\_\_\_\_ medium Eagle 3000 sterilizer with Stage 3 control; \_\_\_\_\_ small Eagle EO sterilizer with Printcon; \_\_\_\_\_ medium Eagle EO sterilizer with Printcon; \_\_\_\_\_ small Eagle EO sterilizer with Round Chart Recorder; \_\_\_\_\_ medium Eagle EO sterilizer with Round Chart Recorder; \_\_\_\_\_ small Medallion EO sterilizer; \_\_\_\_\_ medium Medallion EO sterilizer; \_\_\_\_\_ Eagle gas aerator. Provide \_\_\_\_\_ Sequencer option for multi-sterilizer installation.



## TYPICAL CONFIGURATIONS (Details May Vary)



## SECTION 1: GENERAL INFORMATION

### 1.1 APPLICATION AND DESIGN

The AMSCO EO Disposer System converts ethylene oxide gas into carbon dioxide and water vapor by means of catalytic reaction. The system is designed to interface with AMSCO Eagle and Medallion gas sterilizers to remove EO from the equipment exhaust. The Disposer makes any of these straight gas sterilizers function as a sterilizer/aerator.

The rate of discharge to the Disposer is controlled by the sterilizer and the connecting plumbing. This interface is designed for optimum Disposer performance and maximum destruction of EO. All indicators, controls, and safety devices are included in the system to assure proper system operation and safety.

The AMSCO EO Disposer utilizes a separate vacuum system for exhaust that is not a part of the sterilizer. The sterilizer water ring vacuum pump is necessary for operation of the sterilizer condition but is no longer used during the exhaust phase.

### 1.2 MAJOR COMPONENTS—DISPOSER SYSTEM

The EO Disposer system consists of several major components that will be briefly described in this section. Details of cycle operation are listed in the Principles of Operation section of this manual. The overall Disposer system consists of the following:

1. Disposer unit
2. Sterilizer piping modifications
3. Sterilizer electrical modifications
4. Sequencer Control (necessary only when connecting multiple sterilizers to one Disposer)
5. Exhaust piping from sterilizer to Disposer
6. Communication wiring between sterilizer and Disposer
7. Exhaust stack for Disposer exhaust

The Disposer is sold as a separate product. To operate a sterilizer with a Disposer requires a modification kit ordered for that particular model sterilizer. The modification kit includes items 2 and 3 above, the piping and electrical modifications needed. Separate kits are used for Eagle series sterilizers, Medallion series sterilizers, and Aeration cabinets.

### 1.2.1 Sequencer Kit.

In the event more than one sterilizer (maximum of three) is to be connected to the same Disposer, a Sequencer kit will be needed. This Sequencer allows only one sterilizer at a time to exhaust to the Disposer, thus preventing a Disposer from being overloaded with EO.

### 1.2.2 Additional Materials.

Items 5, 6, and 7 listed in paragraph 1.2 are not supplied by AMSCO and must be obtained locally by the installer. Refer to the Disposer equipment drawing (P-65435-332) for details:

- The exhaust piping between the sterilizer and Disposer is to be as specified in the equipment drawing. 1/4" plugged ports, with shutoff valves, are supplied by AMSCO and are to be located at each end to allow for servicing.
- The communication wiring is to be 12 conductor cable, 22/24 AWG for low voltage (24 VAC) operation.
- The exhaust stack for the Disposer must be capable of handling temperatures up to 500° F. The stack must have a plugged test port for stack sampling. Local codes must be checked for ducting requirements. Refer to the equipment drawing for details.

THIS PAGE  
INTENTIONALLY LEFT  
BLANK

## SECTION 2: OPERATING INSTRUCTIONS

### 2.1 STERILIZER/DISPOSER COMMUNICATIONS

The disposer system consists of an AMSCO EO Disposer, one or more AMSCO sterilizers or Aerators, an interface box for each sterilizer and a sequencer assembly for installations with multiple sterilizers. Refer to Figure 2-1 for a block diagram of the disposer system.

### 2.2 GENERAL OPERATION OF STERILIZER WITH EO DISPOSER

#### 2.2.1

Operation of the sterilizer requires that the Disposer's power switch be ON, signified by "POWER" LED being lit on the sterilizer's interface box.

#### 2.2.2

Cycle selection is handled the same way as with standard sterilizer operation. Sterilizer control automatically turns on Disposer, as signified by "DISPOSER ON" LED lighting up. Note these differences, however:

1. Eagle 2000 and Medallion style sterilizers turn Disposer on when cycle is selected.
2. Eagle 3000 Series sterilizers turn Disposer on upon entering charge phase.
3. Sterilizers equipped with a Disposer run sterilize/aerate cycle only. These units do not run sterilize only cycles. To run a sterilize only cycle it is necessary to perform an AERATE PAUSE or an AERATE ABORT sequence.

Two situations may occur involving interruption of the aeration process, a PAUSE and a PAUSE/ABORT. (See Figures 2-5, 2-6 and 2-7 for sample printouts and round chart record.)

AERATE PAUSE is used when it is necessary to remove sterilized goods or biological indicators from the chamber before aeration is complete. This procedure does not change the requirements for aeration time. In order to assure maximum operator safety and minimize the possibility of exposure to ethylene oxide, the chamber is always swept by an air evacuation pulse just prior to allowing the door to be opened.

AERATE PAUSE/ABORT is used when it is necessary to remove the sterilized load from the sterilizer and transfer it to a separate aerator. The aeration phase of the sterilizer cycle must be aborted to terminate the cycle and initialize the sterilizer control for another cycle start. In order to assure maximum operator safety and minimize the possibility of exposure to ethylene oxide, the same procedure as above is used to initialize the Aerate Pause and then end with Abort.

Refer to section 2.3 for specific instruction on the aerate pause procedures.

4. The cycle times for a sterilizer equipped with Disposer are extended by the operation of the Disposer. Thus, cycles run longer than for units without a Disposer.
  - a. Cycle time is extended approximately 45 minutes for a 16x16-inch or 20x20-inch unit; and approximately 90 minutes for a 24x36-inch unit.
  - b. If multiple sterilizers are connected to a Disposer through a sequencer, cycle times may be extended further, but only if sterilizer cycles are started at the same time. The cause of this is the necessary sequencing of the sterilizer exhaust to prevent overloading the Disposer with EO. See chart 2-1 below.

Chart 2-1

Multiple (MAX. 3) EO Sterilizers(s) cycle times with one Disposer when all sterilizers are started simultaneously. (NOTE: Time shown is additional to standard exhaust phase time.)

| EO Gas Sterilizer | (Scenario) Two Sterilizers | (Scenario) Three Sterilizers |
|-------------------|----------------------------|------------------------------|
| 20x20x38          | Add 95 min to 2nd unit     | Add 150 min to 3rd unit      |
| 24x36x48          | Add 164 min to 2nd unit    | Add 256 min to 3rd unit      |
| 24x36x60          | Add 170 min to 2nd unit    | Add 265 min to 3rd unit      |

For example: If three 20x20x38-inch sterilizers are connected to one Disposer and all machines were started on 130° F Cycle at the same time, cycle times would be:

| Sterilizer #1 | Sterilizer #2      | Sterilizer #3       |
|---------------|--------------------|---------------------|
| 220 min       | 315 min (220 + 95) | 370 min (220 + 150) |




Figure 2-1. Interface Box — Eagle 3000 and 2000 Series Sterilizers

2-2  
764323-253

A- 13

#### 5. Disposer Failure — Goods Locked in Chamber

If the Disposer should fail to become ready, for any reason, the condition may not be detected until a cycle has begun. The sterilized goods cannot safely be removed from the chamber until the gas has been removed from the chamber. Every attempt should be made to repair the Disposer so that the sterilizer may be exhausted normally. The most probable reasons the Disposer has failed to become ready should be investigated first. See Troubleshooting, Section 5.

In the event the Disposer cannot be repaired in a reasonable period of time, it is recommended to allow the sterilizer to remain with gas in the chamber until a qualified service representative can repair the Disposer. The sterilizer will not vent gas through any of the valves. As long as the Envirogard fan is operating, there will not be an appreciable amount of EO lost to the atmosphere.

### 2.3. STERILIZER OPERATION WITH DISPOSER

The AMSCO EO Disposer is completely automatic. The power disconnect switch on the Disposer cabinet door (Figure 2-1) should be left ON at all times except when a qualified technician is servicing the unit. No operator actions are necessary to activate or operate the Disposer. The operator need only initiate the sterilization or aeration procedures in the normal manner. Refer to the Troubleshooting section for more information on causes of alarms.

#### 2.3.1 Disposer Operation with Eagle 3000 Series Sterilizers

1. Interface Box — The sterilizer communicates with the disposer through the interface box. This box is mounted on the sterilizer framework. Eight interface signal lights are located on the Interface Box (see Figure 2-1). The lights signify:

**"POWER"** — This indicates that the Disposer on/off switch on the door is in the ON position. The POWER light will be ON even if the sterilizer is off.

**"DISPOSER ON"** — This indicates that the sterilizer has sent a signal to the Disposer to turn on.

**"VACUUM"** — This indicates that the sterilizer is sending a signal to the Disposer to turn on the vacuum pump.

**"DISPOSER READY"** — This indicates that the Disposer has reached the minimum operating temperature and is ready to accept the sterilizer exhaust.

**"READY OUT"** — This indicates that the interface box has received the DISPOSER READY signal from the Disposer and has communicated this information to the sterilizer control.

**"OVERTEMP"** — This indicates that the Disposer has exceeded 500° F (260° C) and now has shut off. Qualified service personnel will be required to manually reset the Disposer. There are two locations for the manual reset; one is located in the interface box and the other is located in the Disposer. Pressing either resets the Disposer.

**"AIR FLOW"** — This indicates that there is insufficient air flow through the Disposer.

**"CATALYST"** — This indicates that the temperature achieved in the catalyst bed during the sterilizer exhaust phase did not reach a normal minimum operating temperature (320° F [160° C]).

2. Alarms and Signals — Four alarms are added when the Eagle 3000 sterilizer is interfaced with the Disposer. The following will print at the sterilizer as appropriate:

**"ALARM: DISPOSER NOT READY"** — A tone sounds and this message is printed if the Disposer is not ready within 60 minutes after the sterilizer exposure phase is completed. Pressing the #1, #2, #3 or #4 touchpad on the sterilizer column silences the tone while the sterilizer tries to complete the cycle. If the alarm sounds again, call qualified service personnel.

**"ALARM: NO DISPOSER VACUUM"** — A tone sounds and this message is printed if insufficient vacuum is drawn in the Disposer inlet piping. This alarm indicates a failure in the vacuum system. Pressing the #1, #2, #3 or #4 touchpad on the sterilizer silences the tone while the sterilizer tries to complete the cycle. If the alarm sounds again, call qualified service personnel.

**"ALARM: DISPOSER LINE RESTRICTED"** — A tone sounds and this message is printed if the chamber pressure fails to drop at the beginning of the exhaust phase. This alarm indicates a failure in the exhaust plumbing. Pressing the #1, #2, #3 or #4 touchpad on the sterilizer silences the tone while the sterilizer tries to complete the cycle. If the alarm sounds again, call qualified service personnel.

2-3  
764323-253

A- 14

**"CHECK DISPOSER CATALYST"** — Printed if the catalytic cell fails to minimum operating temperature. Repeated occurrence of this alarm indicates that the catalytic cell should be replaced. Call qualified service personnel.

### 3. Aerate Pause and Pause/Abort — Eagle 3000 Sterilizers with Stage 3 Control

The sterilizer must have completed the exhaust phase and entered the aeration phase before an aerate pause or pause/abort can be initiated.

To Pause an Eagle 3000 unit, the operator must:

- Press the #1 cycle selection touch pad. This causes the control to initiate the air evacuation pulse: "Pause Initiated" and "Wait to open door" messages are displayed on the control.
- After the pulse is completed, an alarm tone sounds. The operator then has 5 minutes to open the door, remove the object and re-close and re-lock the door.
- When the door is opened, "Close Door to Resume" and "Press Reset to Abort" are displayed alternately on the control and goods may be removed from the chamber.
- The aeration cycle will resume if the door is reclosed and relocked. Aeration will then continue to completion.

To Pause/Abort an Eagle 3000 unit, the operator must:

- Press the #1 cycle selection touch pad. This causes the control to initiate the air evacuation pulse: "Pause Initiated" and "Wait to open door" messages are displayed on the control.
- After the pulse is completed, an alarm tone sounds. The operator then has 5 minutes to unlock and open the door.
- When the door is opened, "Close Door to Resume" and "Press Reset to Abort" are displayed alternately on the control.
- Press the RESET button while the door is open. The cycle aborts and the display changes to the main cycle selection menu.

### 2.3.2. Disposer Operation with Eagle 2000 Printcon Sterilizers

- Interface Box** — The sterilizer communicates with the disposer through the interface box. This box is mounted on the sterilizer framework. Eight interface signal lights are located on the Interface Box (see Figure 2-1). The lights signify:

**"POWER"** — This indicates that the Disposer on/off switch on the door is in the ON position. The POWER light will be ON even if the sterilizer is off.

**"DISPOSER ON"** — This indicates that the sterilizer has sent a signal to the Disposer to turn on.

**"VACUUM"** — This indicates that the sterilizer is sending a signal to the Disposer to turn on the vacuum pump.

**"DISPOSER READY"** — This indicates that the Disposer has reached the minimum operating temperature and is ready to accept the sterilizer exhaust.

**"READY OUT"** — This indicates that the interface box has received the DISPOSER READY signal from the Disposer and has communicated this information to the sterilizer control.

**"OVERTEMP"** — This indicates that the Disposer has exceeded 500°F (260°C) and now has shut off. Qualified service personnel will be required to manually reset the Disposer. There are two locations for the manual reset; one is located in the interface box and the other is located in the Disposer. Pressing either resets the Disposer.

**"AIR FLOW"** — This indicates that there is insufficient air flow through the Disposer.

**"CATALYST"** — This indicates that the temperature achieved in the catalyst bed during the sterilizer exhaust phase did not reach a normal minimum operating temperature (320°F [160°C]).

- ALARMS and SIGNALS** — Eagle 2000 Series with PRINTCON Control

Four printed alarms are added to the Printcon by the Disposer controls. The following will print at the sterilizer as appropriate:

**"ALARM: DISPNRDY" (Not Ready)** — A buzzer sounds and this message is printed if the Disposer has not reached operating temperature within 60 minutes after the sterilizer exposure phase is completed. Pressing the Gas Cycle selector on the sterilizer will silence the buzzer and allow the Disposer to try and reach operating temperature

```

-----
--- HIGH TEMP ---
-----
CYCLE START AT 9:53:31A
ON 10/05/90

CYCLE COUNT 00303
OPERATOR
STERILIZER

STER TEMP = 130.0 F
PREHEAT = 20 M
STER TIME = 0:05 H:M
AERATE TIME = 2:00 H:M

P=psig
V=inHg
- TIME T=F V=
H 9:53:32A 127.8 0.0P
C 10:13:30A 131.1 0.0P
C 10:24:55A 130.4 26.6V
C 10:28:07A 129.8 26.6V
C 10:31:22A 128.7 26.6V
C 10:34:29A 127.2 26.6V
S 10:43:06A 130.2 8.5P
E 10:48:05A 128.5 8.7P
E 11:03:44A 127.5 22.0V
E 11:04:20A 130.5 1.9V
E 11:11:54A 125.9 22.0V
E 11:12:31A 132.0 1.9V
A 11:34:18A 132.1 0.0P
* SET AERATION TIME
COMPLETED AT 1:44:05P
E 1:50:07P 139.6 0.7V

LOAD 100501

TEMP MAX=130.2 F
TEMP MIN=128.3 F

PREHEAT = 0:20:00
CONDITION = 0:29:36
STERILIZE = 0:05:00
EXHAUST = 0:46:14
AERATE = 2:15:51
TOTAL CYCLE = 3:56:41

COMPLETE CYCLE 1:50:07P
ON 10/05/90

- READY TO UNLOAD -

* DOOR 1:50:07P
UNLOCKED

```

```

-----
--- HIGH TEMP ---
-----
CYCLE START AT 9:53:31A
ON 10/05/90

CYCLE COUNT 00303
OPERATOR
STERILIZER

STER TEMP = 130.0 F
PREHEAT = 20 M
STER TIME = 0:05 H:M
AERATE TIME = 2:00 H:M

P=psig
V=inHg
- TIME T=F V=
H 9:53:32A 127.8 0.0P
C 10:13:30A 131.1 0.0P
C 10:24:55A 130.4 26.6V
C 10:28:07A 129.8 26.6V
C 10:31:22A 128.7 26.6V
C 10:34:29A 127.2 26.6V
S 10:43:06A 130.2 8.5P
E 10:48:05A 128.5 8.7P
E 11:03:44A 127.5 22.0V
E 11:04:20A 130.5 1.9V
E 11:11:54A 125.9 22.0V
E 11:12:31A 132.0 1.9V
A 11:34:18A 132.1 0.0P
P 11:56:18A 131.9 6.7P
X 12:03:43P 131.9 1.0V
D 12:04:01P 131.1 3.6V
* SET AERATION TIME
COMPLETED AT 1:44:05P
E 1:50:07P 139.6 0.7V

LOAD 100501

TEMP MAX=130.2 F
TEMP MIN=128.3 F

PREHEAT = 0:20:00
CONDITION = 0:29:36
STERILIZE = 0:05:00
EXHAUST = 0:46:14
AERATE = 2:15:51
TOTAL CYCLE = 3:56:41

COMPLETE CYCLE 1:50:07P
ON 10/05/90

- READY TO UNLOAD -

* DOOR 1:50:07P
UNLOCKED

```

Figure 2-2. Eagle 3000 with Disposer: Sample Tape

for an additional 60 minutes. If the alarm sounds again, the condition must be checked by qualified service personnel.

**"ALARM: DISP FLOW"** — A buzzer sounds and this message is printed if the Disposer piping vacuum switch fails to close when sensed. This switch is sensed at the beginning of each vacuum pulse to check the integrity of the connecting piping and the Disposer vacuum pump. If the VS-23 vacuum switch is not activated within 60 seconds, the alarm sounds. This alarm indicates a failure in the Disposer vacuum system. Pressing the Gas Cycle selector on the sterilizer silences the buzzer and the control allows another 60 seconds for the unit to reach VS-23. If the alarm sounds again, the condition must be checked by qualified service personnel.

**"ALARM: DISP EXHST"** (Exhaust) — A buzzer sounds and this message is printed if the sterilizer has been in exhaust for too long. This alarm indicates a failure in the exhaust plumbing. Pressing the Gas Cycle selector on the sterilizer will silence the buzzer while the sterilizer tries to complete the cycle. If the alarm sounds again, the condition must be checked by qualified service personnel.

**"CHECK DISP CATALYST"** — Printed if the catalytic cell fails to achieve a minimum operating temperature during exhaust. Repeated occurrence of this alarm indicates that the catalytic cell should be replaced. See section 6.2.2, Catalytic Cell, for replacement procedure.

### 3. AERATE PAUSE and PAUSE/ABORT — Eagle 2000 Sterilizers with Printcon

To Pause a Printcon unit, the operator must:

- Wait until the sterilizer is in the Aeration phase ("Aerate" light on the primary panel is lit), or, until time on the primary panel is between 11:00 and 12:00 (if primary panel does not have an "AERATE" light.)
- Press the GAS cycle button. This will cause the control to initiate the air evacuation pulse. The sterilizer then pulls one post vacuum pulse.
- After the pulse is completed (approximately 20 minutes), the buzzer sounds for 20 seconds and the aeration time is blanked out. The operator then has 10 minutes to open the door, remove the object and reclose and relock the door.

**NOTE:** If the door is not opened within 10 minutes, the sterilizer will continue the cycle automatically.

- When the door is opened, "ALARM: DOOR OPEN" will print and the indicator will light. Goods may be removed from the chamber.
- The aeration cycle will resume if the door is reclosed and relocked, and will continue until the cycle is complete.
- If the door is left opened for more than 10 minutes, the buzzer will sound. If the cycle is to be continued, close and lock the door. The buzzer will cease, and the cycle will resume.

To Pause/Abort a Printcon unit, the operator must:

- Wait until the sterilizer is in the Aeration phase ("Aerate" light on the primary panel is lit), or until time on primary panel is between 11:00 and 12:00 (if primary panel does not have an "AERATE" light).
- Press the GAS cycle button. This will cause the control to initiate the air evacuation pulse. The sterilizer then pulls one post vacuum pulse.
- After the pulse is completed (approximately 20 minutes), the buzzer sounds for 20 seconds. The operator then has 10 minutes to press the gas cycle select button.
- When the gas cycle select button is pressed, the complete lamp and cycle complete buzzer will come on for 90 seconds or until the door is opened. The Cycle Summary will print.

**NOTE:** Do not press RESET.


- If the door is not opened within 10 minutes, the cycle will continue automatically.

### 2.3.3 Disposer Operation with Eagle 2000 IRC Sterilizer.

- Interface Box** — The sterilizer communicates with the disposer through the interface box. This box is mounted on the sterilizer framework. Eight interface signal lights are located on the Interface Box (see Figure 2-1). The lights signify:

**"POWER"** — This indicates that the Disposer on/off switch on the door is in the ON position. The POWER light will be ON even if the sterilizer is off.

**"DISPOSER ON"** — This indicates that the sterilizer has sent a signal to the Disposer to turn on.



**"VACUUM"** — This indicates that the sterilizer is sending a signal to the Disposer to turn on the vacuum pump.

**"DISPOSER READY"** — This indicates that the Disposer has reached the minimum operating temperature and is ready to accept the sterilizer exhaust.

**"READY OUT"** — This indicates that the interface box has received the DISPOSER READY signal from the Disposer and has communicated this information to the sterilizer control.

**"OVERTEMP"** — This indicates that the Disposer has exceeded 500°F (260°C) and now has shut off. Qualified service personnel will be required to manually reset the Disposer. There are two locations for the manual reset; one is located in the interface box and the other is located in the Disposer. Pressing either resets the Disposer.

**"AIR FLOW"** — This indicates that there is insufficient air flow through the Disposer.

**"CATALYST"** — This indicates that the temperature achieved in the catalyst bed during the sterilizer exhaust phase did not reach a normal minimum operating temperature (320°F [160°C]).

2. **ALARMS and SIGNALS** — A remote alarm box (see Figure 2-4) is mounted near the sterilizer to provide an audible and visual indication of the following alarm conditions:




Figure 2-4. Remote Alarm Box — Eagle 2000 IRC Sterilizers

**"NOT READY"** — A buzzer in the sterilizer sounds and this lamp lights if the Disposer has not reached operating temperature within 60 minutes after the sterilizer exposure phase is completed. Pressing the Gas Cycle selector on the sterilizer silences the buzzer, allowing the Disposer to try

and reach operating temperature for an additional 60 minutes. If the alarm sounds again, the condition must be checked by qualified service personnel.

**"NO FLOW"** — A buzzer sounds and this lamp lights if the Disposer piping vacuum switch fails to close when sensed. This switch is sensed at the beginning of each vacuum pulse to check the integrity of the connecting piping and the Disposer's vacuum pump. If the VS-23 vacuum switch is not activated within 60 seconds, the alarm sounds. This alarm indicates a failure in the Disposer vacuum system. Pressing the Gas Cycle selector on the sterilizer will silence the buzzer and the control allows another 60 seconds to reach the VS-23 set point. If the alarm sounds again, the condition must be checked by qualified service personnel.

**"CHECK CATALYTIC BED"** — This lamp lights if the catalytic cell fails to achieve a minimum operating temperature during exhaust. Repeated occurrence of this alarm indicates that the catalytic cell should be replaced. See section 6.2.2, Catalytic Cell, for replacement procedure.

**"TOO LONG IN EXHAUST"** — A buzzer sounds and this lamp lights if the sterilizer has been in exhaust too long, indicating a failure in the exhaust plumbing. Pressing the Gas Cycle selector on the sterilizer will silence the buzzer while the sterilizer tries to complete the cycle. If the alarm sounds again, the condition must be checked by qualified service personnel.

3. **AERATE PAUSE and PAUSE/ABORT** — Eagle 2000 with IRC

To Pause an IRC unit, the operator must:

- a. Wait until the sterilizer is in the Aeration phase ("AERATE" light on the primary panel is lit).
- b. Press the GAS cycle button. This causes the control to initiate the air evacuation pulse. The sterilizer then pulls one post vacuum pulse.
- c. After the pulse is completed (approximately 20 minutes), the buzzer sounds for 20 seconds and the aeration time is blanked out. The operator then has 10 minutes to open the door, remove the object and reclose and relock the door.
- d. When the door is opened, the "DOOR UNLOCKED" lamp lights and goods may be removed from the chamber.




Figure 2-5. IRC Sterilizer with Disposer: Sample Chart

- e. The aeration cycle resumes if the door is reclosed and relocked, and continues until the cycle is complete.
- f. If the door is left opened for more than 5 minutes, the buzzer sounds. If the cycle is to be continued, close and lock the door. The buzzer ceases, and the cycle resumes.

To Pause/Abort an IRC unit, the operator must:

- a. Wait until the sterilizer is in the Aeration phase ("AERATE" light on the primary panel is lit).
- b. Press the GAS cycle button. This causes the control to initiate the air evacuation pulse. The sterilizer then pulls one post vacuum pulse.
- c. After the pulse is completed, the buzzer sounds for 20 seconds. The operator then has 10 minutes to press the gas select button.
- d. When the GAS select button is pressed, the complete lamp and cycle complete buzzer sounds for 90 seconds, or until the door is opened.

### 2.3.4 Disposer Operation with Medallion Sterilizers

1. Interface box — The sterilizer communicates with the disposer through the interface box. This box is mounted on the sterilizer framework. Thirteen signal lights are located on the Medallion Interface Box (See Figure 2-7) to signify:

**"EXHAUST"** — This indicates the sterilizer is in the exhaust phase. The light will go out when the sterilizer is in the aeration phase.

**"WATER PURGE"** — This light will be on for 60 seconds at the end of the exhaust phase to indicate that the sterilizer vacuum pump is removing water from the drain piping.

**"AERATE"** — This light indicates the sterilizer is in the aeration phase.

**"COMPLETE"** — This light indicates the sterilizer has completed the 12-hour aeration and is ready to be unloaded.

**"OFF"** — This indicates the sterilizer rotary switch is in the OFF position.

**"POWER"** — This indicates that the Disposer on/off switch on the door is in the ON position. The POWER light will be ON even if the sterilizer is off.

**"DISPOSER ON"** — This indicates that the sterilizer has sent a signal to the Disposer to turn on.

**"VACUUM"** — This indicates that the sterilizer is sending a signal to the Disposer to turn on the vacuum pump.

**"DISPOSER READY"** — This indicates that the Disposer has reached the minimum operating temperature and is ready to accept the sterilizer exhaust.

**"READY OUT"** — This indicates that the interface box has received the DISPOSER READY signal from the Disposer and has communicated this information to the sterilizer control.

**"OVERTEMP"** — This indicates that the Disposer has exceeded 500 F (260 C) and now has shut off. Qualified service personnel will be required to manually reset the Disposer. There are two locations for the manual reset; one is located in the interface box and the other is located in the Disposer. Pressing either resets the Disposer.

**"AIR FLOW"** — This indicates insufficient air flow through the Disposer.

**"CATALYST"** — This indicates the temperature achieved in the catalyst bed during the sterilizer exhaust phase did not reach a normal minimum operating temperature (320° F [160° C]).

### 2. ALARMS and SIGNALS — Medallion Series Sterilizers

A remote alarm box (see Figure 2-6) is mounted near the sterilizer to provide an audible and visual indication of the following alarm conditions:

**"NOT READY"** — A buzzer sounds and this lamp lights if the Disposer has not reached minimum operating temperature within 60 minutes after the sterilizer exposure phase is completed. Pressing the Gas Cycle selector on the sterilizer silences the buzzer allowing the Disposer to try and reach operating temperature for an additional 60 minutes. If the alarm sounds again, the condition must be checked by qualified service personnel.

**"TOO LONG IN STEP/NO AIR FLOW"** — A buzzer sounds and this lamp either blinks or lights steadily to indicate:

- Blinking Light — Too long in exhaust or air break.
- Steady light — Disposer piping vacuum switch fails to close.

Either condition indicates a failure in the sterilizer exhaust piping or the disposer vacuum supply. Call qualified service personnel.

**"CHECK CATALYTIC BED"** — This lamp lights if the catalytic cell fails to achieve minimum operating temperature. Repeated occurrence of this alarm indicates that the catalytic cell should be replaced. See section 6.2.2, Catalytic Cell, for replacement procedure.




Figure 2-6. Remote Alarm Box — Medallion Sterilizer

**"SILENCE BUZZER/AERATE PAUSE"** — This button on the remote alarm box allows the operator to silence alarm buzzers or to initiate an aerate pause.

### 3. Aerate Pause and Pause/Abort — Medallion Gas Sterilizers

To Aerate/Pause a Medallion unit, the operator must first wait until unit is in AERATE, as indicated by the "AERATE" display on the rotary switch on the control panel. Then:

- a. Press the SILENCE BUZZER/AERATE PAUSE button. This causes the control to initiate the air evacuation pulse. Note that the AERATE PAUSE/DISP. NOT READY lamp lights steadily at this time.
- b. After the pulse is completed, the buzzer sounds for 20 seconds. The operator then has 5 minutes to open the door, remove the object and reclose and relock the door.
- c. The aeration cycle resumes if the door is reclosed and re-locked, and continues to cycle completion. If the door is not re-locked within five minutes, the BUZZER begins pulsating, reminding the operator to re-lock the door. The BUZZER silences automatically when the door is re-locked.

Aerate Pause/Abort is used when it is necessary to remove the sterilized load from the sterilizer and transfer it to a separate aerator. The aeration phase of the sterilizer cycle must be aborted to terminate the cycle and initialize the sterilizer control for another cycle start. In order to assure maximum operator safety and minimize the possibility of exposure to ethylene oxide, the same procedure as above is used to initialize the Aerate Pause and then end with Abort.

To Aerate Pause/Abort a Medallion unit, the operator must first wait until unit is in AERATE, as indicated by the "AERATE" display on the rotary switch on the control panel. Then:

- a. Press the SILENCE BUZZER/AERATE PAUSE button. This causes the control to initiate the air evacuation pulse. Note that the AERATE PAUSE/DISP. NOT READY lamp lights steadily at this time.
- b. After the pulse is completed, the buzzer sounds for 20 seconds. The operator then has 5 minutes to abort the cycle by pressing the SILENCE BUZZER/AERATE PAUSE button.
- c. When the SILENCE BUZZER/AERATE PAUSE button is pressed, the control advances to the COMPLETE phase. The STERILIZE lamp lights and the BUZZER sounds for 90 seconds, or until the door is opened.

When the door is opened (or STOP pushbutton pressed, if equipped), the control advances to the OFF/READY position.




Figure 2-7. Interface Box — Medallion Sterilizer

## 2.4 SEQUENCER — GENERAL

The POWER and AUTO/MANUAL switches on the Sequencer should be set to ON and AUTO at all times except when a qualified technician is servicing the unit. The Sequencer microprocessor control allows up to three sterilizers and up to three aerators to automatically vent EO to a single AMSCO Disposer. The Sequencer control allows one sterilizer at a time to vent to the Disposer. The aerator(s) continuously vents to the Disposer through a separate connection. Multiple aerators are wired and piped in parallel.

The three sterilizers connected to the Sequencer are designated #1, #2 and #3 for reference only because the Sequencer shows no preference or priority among the sterilizers. No modifications are required to the sterilizers when installing a Sequencer system.

### 2.4.1 Indicator Lights — Sequencer

LED lamps on the front panel of the Sequencer control indicate operational status as follows:

**"AERATOR ON"** — This indicates that an aerator has been turned on and is sending a signal to the Disposer to turn on.

**"STER 1 ENABLE"** — This light can be operating in one of three ways. It will be on solid if the sterilizer designated as #1 is exhausting. The light will sequence on for two seconds and then be off for 10 seconds if the sequencer is "scanning" the sterilizer to see if it is ready for exhaust. Finally, the light will blink on and off in one second intervals if sterilizer #1 is ready to exhaust but the Disposer is Not Ready.

**"STER 2 ENABLE"** — Operation of these lights is same as above, except designated sterilizer is #2.

**"STER 3 ENABLE"** — Operation of these lights is same as above, except designated sterilizer is #3.

**"POWER"** — This indicates that the Disposer on/off switch on the door is in the ON position. The POWER light will be ON even if the sterilizer is off.

**"DISPOSER ON"** — This indicates that a sterilizer has sent a signal to the Disposer to turn on.

**"VACUUM"** — This indicates that the sterilizer is sending a signal to the Disposer to turn on the vacuum pump.

**"DISPOSER READY"** — This indicates that the Disposer has reached the minimum operating temperature and is ready to accept the sterilizer exhaust.

**"READY OUT"** — This indicates that the interface box has received the DISPOSER READY signal from the Disposer and has communicated this information to the sterilizer control.

**"OVER TEMP"** — This indicates that the Disposer has exceeded 500° F (260° C) and now has shut off. Qualified service personnel will be required to manually reset the Disposer. There are two locations for the manual reset; one is located in the interface box and the other is located in the Disposer. Pressing either will reset the Disposer.

**"AIR FLOW"** — This indicates that there is insufficient air flow through the Disposer.

## 2.5 AERATORS — GENERAL

As with sterilizers, the EO Disposer is completely automatic in processing vent air from an aerator. The operator need only initiate normal aeration procedures. The aerator will not begin operation until the Disposer has reached a READY state. Once the Disposer has reached READY, the aerator's HEAT and OPERATING lamps will come on.

## 2.6 DISPOSER

1. General operation — unit is completely automatic and is controlled by the sterilizer. The only operator involvement is to check that the power switch is ON when operating the sterilizer. The switch should remain ON unless being serviced by a qualified technician.

2. Indicator lights

Five lights are located inside the Disposer cabinet to signify:

- a. Disposer On — indicates Disposer power is on.
- b. Ready — Disposer has reached minimum operating temp.
- c. Low Air Flow — insufficient air flow through the Disposer.
- d. Oper Limit — Disposer has reached 450° F (232° C).
- e. Overtemp — Disposer has reached 500° F (260° C). This remains latched on until ILPB-1 reset is pressed. The ILPB-1 reset is located next to the overtemp light.

## Section 3 — Principles of Operation

### 3.1 GENERAL

When a sterilization cycle is initiated, the sterilizer will activate the EO Disposer, which then begins heating the catalytic cell. It is necessary for the Disposer catalytic bed to be at operating temperature before processing EO. The time for the catalytic bed to reach operating temperature is dependent on the ambient temperature.

The sterilizer cycle is normal during the condition and exposure phase. After the exposure phase of the sterilization cycle (or if the cycle is aborted after the chamber is charged with gas), the gas mixture in the chamber is vented to the Disposer. Two vacuum pulses are drawn and a 20-minute air wash follows. After the air wash is complete, the sterilizer vacuum pump operates for 60 seconds to remove accumulated moisture in the drain plumbing. Then the sterilizer proceeds to the aeration phase whereby one post vacuum is pulled per hour for a 12 hour period.

As the gas is vented to the Disposer it is mixed with air (see Figure 3-1). The ethylene oxide is removed from the air by the equation:  $2\text{CH}_3\text{CH}_2\text{O} + 5\text{O}_2$ , catalyst  $4\text{CO}_2 + 4\text{H}_2\text{O} + \text{heat}$ . Ambient air is pulled into the Disposer through the air inlet (1) and the pre-filter (2) by the blower (3). The air then passes through the electric heater (4) where it is heated, and then into the transition plenum (5) where EO is fed into the airstream through a solenoid valve (6). EO is pulled from the sterilizer by vac pump (7). As EO enters the airstream it is diluted by the preheated inlet air and then directed into the catalytic cell (8) where the EO is catalyzed in accordance with the above reaction. The heat given off by this reaction raises the temperature of the catalytic cell and the effluent airstream in proportion to the amount of EO introduced into the stream to a maximum of 500° F (260° C).

From the catalytic cell, air, carbon dioxide and water vapor and refrigerant 12 are discharged through the air outlet (9) to the atmosphere.

The sterilizer regulates the amount of EO released to the Disposer to maintain normal operating temperature well below the maximum allowable temperature. Should an event occur in the Disposer that causes the temperature to rise to the maximum, the Disposer signals the sterilizer that the maximum temperature has been reached and that the Disposer is no longer in the READY status. The sterilizer immediately suspends release of EO to the Disposer until the

Disposer cools to approximately 440 F and the Disposer Overtemp Reset is pressed. After cooling, the Disposer signals the sterilizer that it is again READY and the process resumes.




Figure 3-1. Disposer Air Flow

### 3.2 CYCLE DESCRIPTION — DISPOSER

Refer to Disposer electrical schematic (Donaldson P19-3610, sheet 2.)

Incoming 208 VAC passes through a door safety interlock switch and feeds the blower contactor M1, heater contactor R1, a 208/220 VAC step up transformer for the vacuum pump contactor M2, and a 220/110 VAC step down transformer for the Disposer control. The Disposer is turned on by CR6 relay contacts. CR6 coil is powered by the sterilizer control "Disposer On signal" or by the LS-2 service by-pass switch (located by the door hinge).

When CR6 contacts close, the "Disposer On" light comes on and M1 blower contactor pulls in (causing the blower to operate). M1 is controlled by a 30-



minute timer. This timer keeps the blower running for 30 minutes after the Disposer is shut off. The purpose for this delay is to allow the heaters to cool down properly. With the blower operating, the PS-1 air flow switch contacts close. The air flow switch contacts power up CR1 relay. A red indicating light (indicating Low Air Flow) turns off. The PS-2 switch monitors the differential pressure between the inside and outside atmosphere of the Disposer. When the differential pressure exceeds 1.8" (indicating a plugged filter), PS-2 opens, de-energizing CR1.

The temperature controller in the Disposer (ITG-2600 on the schematic) is a microprocessor device which monitors the catalytic bed temperature using a type "J" thermocouple. The controller has four outputs which control relays CR2, CR3, CR4, and CR5. CR2 is powered when 280° F (138° C) is achieved, CR3 is powered when 320° F (160° C) is achieved, CR4 is de-energized when 450° F (232° C) is achieved, and CR5 is powered when 500° F (260° C) is achieved.

R1 heater control relay is energized through CR1 and CR4 contacts. The heaters energize through R1 relay. R1 relay is controlled by two thermostats. One thermostat, TC-2, is the control thermostat and is set at 320° F (160° C). The other thermostat, TC-1, is an overtemperature safety set at 550° F (270° C). Note the control thermostat has its RTD probe in the chamber just above the catalyst bed, whereas the high limit thermostat has its RTD brazed to the bottom heater element. This latter RTD location is required for electrical codes.

When the catalyst bed reaches 280° F (138° C) the disposer temperature controller powers CR2. CR7 relay is then powered up through CR5, CR2, CR4 and CR1 contacts. A Disposer Ready signal is then sent to the sterilizer through CR7 contacts, the green "READY" light in the Disposer is lit, and FSV-1 EO injection solenoid is energized. FSV-1 solenoid is located in the Disposer. This valve opens to allow the Disposer to accept EO from the sterilizer exhaust piping. If the catalyst bed goes to an excessive operating temperature condition (over 450° F [232° C]), the temperature controller will de-energize CR4 relay. CR4 contacts would then open, causing CR7 to lose power and remove the disposer ready signal to the sterilizer.

When the sterilizer is ready to exhaust, a signal from the sterilizer (Disposer Vacuum) turns on CR8 relay. CR8 contacts close and power up M2 vacuum pump contactor. This allows a vacuum to be pulled in the piping between the disposer and sterilizer.

### 3.2.1 Disposer READY Mode Drop-out

The disposer will drop its Ready mode if the following conditions occur:

1. Inadequate pressure drop across the catalyst bed (PS-1 opens, de-energizes CR1).
2. Plugged air filter (PS-2 opens, de-energizes CR1).
3. Catalyst bed temperature drops below 280° F (138° C) (CR2 will de-energize).
4. Catalyst bed temperature rises above 450° F (232° C) (CR4 will de-energize).
5. Catalyst bed temperature rises above 500° F (260° C) (CR5 will energize). This acts as a redundant heater shutoff to the 450° F (232° C) switch above.
6. Sterilizer signal to power up the disposer is removed (CR6 will de-energize).

### 3.2.2 Disposer Outputs

The outputs from the disposer to the sterilizer are as follows:

1. Disposer Ready (CR7 contact)
2. Low Air Flow condition (CR1 contact)
3. Operating temperature (320° F [160° C]) achieved (CR3 contact)
4. Disposer overtemp (500° F [260° C]) (CR5 contact)

### 3.2.3 Other Disposer Communications

The light panel in the Disposer has five indicating lights:

1. Disposer On - indicates Disposer is turned on by the sterilizer
2. Ready - Disposer has reached minimum operating temp
3. Low Air Flow - insufficient air flow
4. Oper. Limit - Disposer has reached 450° F (232° C).
5. Overtemp - Disposer has reached 500° F (260° C). This will remain latched on until ILPB-1 reset is pressed. The ILPB-1 reset is located next to the overtemp light.

## 3.3 CYCLE DESCRIPTION — EAGLE 3000 STERILIZER OPERATION WITH DISPOSER

### 3.3.1 General

The Eagle 3000 control EPROM contains cycle information for operation with the EO Disposer. To access the program it is necessary to actuate Dip switch 3 on the control board to ON. It is then necessary to perform a manufacture's burn-in to load Disposer related programming into the memory.

The 3000 sterilizer also comes with a wiring harness with solenoid coil connectors ready for Disposer plumbing.

To convert the unit to be Disposer compatible a drain line modification is installed and a electrical interface box is installed.

### 3.3.2 Disposer On Signal

The Eagle 3000 control does not turn on the EO Disposer until the sterilizer starts its charge phase. This is different from the Eagle 2000 and Medallion sterilizers, which start the Disposer as soon as a cycle is selected. The reason for this is that the 3000 will not abort to a Exhaust phase unless EO has been admitted to the chamber. If the reset button is pressed before EO is admitted to the chamber, the control will simply cause the chamber to air break. Once EO is admitted, pressing the reset button will cause the unit to progress through the two post vacuums and air wash before the door can be opened.

### 3.3.3 Exhaust phase

At the completion of the exposure phase the sterilizer will progress to exhaust. First the control will check for a Disposer Ready signal from wire 65 from the interface box. If present, the sterilizer control turns on the Disposer vacuum pump and the S35 gas purge line solenoid valve. On some early units the purge line transducer measures the vacuum level in the exhaust line. On current units a separate VS-23 vacuum switch is used. If 180 mm Hg (approx 23" Hg) is measured within 60 seconds, the control will then close S35, it will energize S41 exhaust valve, and the Disposer vacuum pump then pulls a vacuum on the chamber.

At 10" Hg as measured by the chamber transducer, on 24x36 units a second exhaust valve, S42, opens allowing a faster exhaust rate. The Disposer vacuum pump continues to pull down until 180 mm Hg (approx 23" Hg) as measured by the chamber transducer, is reached in the chamber. At that point the unit goes to Air Break.

### 3.3.4 Air Break

The Disposer vacuum pump turns off and S41 (and S42 on 24x36 units) turn off. The filtered air valve S1 turns on to let the chamber break to VS2 setting.

### 3.3.5 Second post vac and air break

The unit pulls a post vacuum and air breaks in the same manner as the first, except on 24x36 units where S42 opens at the same time as S41.

### 3.3.6 Air flush

At the end of the second air break the Disposer vacuum pump turns back on and the S1 valve remains on. This occurs for a preprogrammed 20 minute time period.

### 3.3.7 Water dump

At the end of the 20 minute air flush, the control turns off the Disposer vacuum pump, S41 (and S42 on 24x36 units) and the S1 valve. The sterilizer vacuum pump is turned on and S28 is turned on, pulling water from the exhaust piping. This occurs for 60 seconds.

### 3.3.8 Aeration

The control continues to aeration, pulling a post vacuum in the same manner as above once per hour.

## 3.4 CYCLE DESCRIPTION - EAGLE 2000 PRINTCON STERILIZER AND EAGLE 2000 INDICATOR ROUND CHART RECORDER (IRC) STERILIZERS

Refer to sterilizer electrical schematic (146653-873 for Printcon units, and 146655-121 for IRC units), disposer electrical schematic, disposer piping schematic, and cycle graph.

Note that all communication between the sterilizer control and the Disposer is made through the interconnect box mounted on the sterilizer. The Disposer Interface board in the Interconnect box provides optical signal isolation between the sterilizer and Disposer. The Disposer provides 24 VAC to the photocouplers. See sterilizer schematic sheet 3 of 3 for reference connection.

### 3.4.1 Cycle start

The Disposer will be off when the sterilizer is in the off or idle mode. When a cycle is initiated, the Disposer receives a signal from the sterilizer to start up (wire 828 from the sterilizer to the interconnect box). When the Disposer is at its minimum operating tem-

perature (280 F [138 C]) a signal is sent back to the sterilizer, through wire 823 (from the interconnect box to the sterilizer), that the Disposer is Ready to receive EO.

#### 3.4.2 Exhaust

When the exposure phase times out, the sterilizer control waits for the Disposer Ready signal (wire 823 goes to 0 volts DC when Disposer is ready; if Disposer is not ready, wire 823 is at 5 volts DC). If the Disposer is not ready within 60 minutes after the exposure phase is completed, a Disposer Not Ready alarm will occur. Printcon units will print "ALARM: DISPOSER NRDY"; IRC units will display the "DISPOSER NOT READY" light on the alarm signal box. The sterilizer buzzer sounds intermittently. The buzzer will sound until the Disposer is ready. Pressing the GAS button will silence the buzzer until another 60 minutes passes. If still not ready, the alarm occurs again. The buzzer will silence if the Disposer becomes ready.

If the Disposer ready signal is present, the sterilizer sends a signal to the Disposer, through wire 829, to turn on the Disposer vacuum pump. This causes a vacuum to be pulled in the piping between the sterilizer and Disposer. The control checks vacuum switch VS-23 (set at 23" Hg). If the switch is not closed within 60 seconds, indicating 23" Hg has not been reached, the Disposer flow alarm will occur. Printcon units will print "ALARM: DISPOSER FLOW"; IRC units will display the "DISPOSER FLOW" light on the alarm box. The sterilizer will wait until the 23" Hg switch is made before proceeding. This is referred to as the "Vacuum Test". The flow alarm will also sound if VS-23 does not open within 5 seconds after S41/S42, which open during the "Vacuum Test".

With the 23" Hg switch made, the sterilizer turns off the chamber heat (S18), the solenoid valve S41 (sterilizer exhaust) opens, and the sterilizer exhausts to the Disposer. If the chamber fails to reach a specific vacuum level (2" Hg on 20x20 sterilizers, 10" Hg on 24x36 sterilizers) within 60 minutes, a too long in exhaust alarm occurs. Printcon units will print "ALARM: DISP. EXHST"; IRC units will display the "TOO LONG IN EXHAUST" light on the alarm box. An intermittent alarm sounds. The alarm can be silenced by pressing the GAS cycle select button. If the Disposer Ready signal drops out at any time while the disposer vacuum pump is on, the pump is shut off and S41/S42 close until the Disposer Ready signal is restored.

#### 3.4.3 First Post Vacuum and Air Break (Figures 3.2 and 3.3)

On 24x36 sterilizers a second exhaust solenoid valve, S42, opens at 10" Hg to allow a faster exhaust rate. On Printcon units this 10" Hg set point is programmed into the Printcon control and is not adjustable. On IRC units the 10" Hg set point is controlled by the VS-10 vacuum switch.

The chamber exhausts to the VS-23 switch setting (23" Hg). The control allows 40 minutes for the chamber to go from 2" Hg to 23" Hg on 20x20 sterilizers or from 10" Hg to 23" Hg on 24x36 sterilizers; if 23" Hg is not reached within 40 minutes a Disposer flow alarm will occur. The flow alarm will also sound if VS-2 is not activated when 23" pressure is reached during the vacuum pulse. Once 23" is reached, solenoid S41 (and S42 on 24x36 sterilizers) closes, solenoid S1 opens, and the sterilizer air breaks to VS-2 setting (2" Hg).

#### 3.4.4 Second Post Vacuum and Air Break (Figures 3.2 and 3.3)

After reaching VS-2 the second post vacuum begins. The sterilizer looks for the Disposer ready signal. If not present, the control waits for 60 minutes (during this time the sterilizer is waiting for the Disposer ready signal, the chamber heaters will come on to maintain temperature). If still not present at the end of 60 minutes, the Disposer not ready alarm sequence occurs. When the Disposer ready signal is received, the control turns on the Disposer vacuum and looks for the VS-23 switch. If switch is not made within 60 seconds, the Disposer flow alarm sequence occurs. Once the 23" switch is made the chamber heat is turned off and S41 opens, allowing a vacuum to be pulled on the chamber. The 23" switch must deactivate within 5 seconds or the Disposer flow alarm sequence will occur. The sterilizer control allows 40 minutes to pull from the VS-2 setting to VS-23.

After reaching VS-23 switch, solenoid S41 (and S42) closes and S1 opens, allowing the chamber to air break.

At the beginning of the exhaust phase, the control activates a "service pulse", wire 826 to the interconnect box as an input to a flip-flop circuit on the Interface board. If the catalyst bed did not reach an efficient operating temperature (320 F [160 C]) during the post vacuums, this indicates that an insufficient catalytic reaction has taken place. A catalyst check alarm will occur. On Printcon units the printer will print "CHECK DISP CATALYST"; on IRC units the "CHECK CATALYTIC BED" lamp will light up on

the alarm signal box. Repeated occurrences of this printout indicate the catalyst bed requires replacement. If the temperature setting was reached, the service pulse will not set the flip-flop circuit and no alarm will occur. The service pulse is not activated on an aborted cycle.

#### 3.4.5 Air Flush (Figures 3.2 and 3.3)

Upon reaching VS-2 after air breaking, the control turns the chamber heater control on. The control checks to see if the Disposer ready signal is still present. The vacuum pump turns on and a vacuum is pulled in the piping. The 23" Hg switch must be made within 60 seconds. When switch is made, the chamber heat is turned on and solenoids S41 (and S42 on 24x36 sterilizers) and S1 open, allowing the chamber to air flush for 20 minutes. If at any time during the air flush the Disposer ready signal is lost, solenoids S41 (and S42 on 24x36 sterilizers) and S1 will close and the sterilizer will wait for the disposer ready signal before restarting. At the end of 20 minutes, S41 (and S42 on 24x36 sterilizers) closes, allowing the chamber to air break to VS-2. S1 closes when 2" is reached. The chamber heat comes back on to maintain temperature.

#### 3.4.6 Water Dump

At this point there is some condensed moisture in the drain piping that must be removed. The vacuum pump and S7 energize and solenoid S28 opens, allowing the water to be pulled into the vacuum pump. At the end of 60 seconds, S7, S28, and the vacuum pump de-energize.

#### 3.4.7 Aeration (Figures 3.2 and 3.3)

After the air flush and water dump, the aeration phase begins. The exhaust lamp goes off and the aerate lamp comes on. The aerate timer starts the 12-hour count and the sterilizer pulls a vacuum pulse and repeats the pulse once an hour. The chamber temperature is controlled at set point, between the vacuum 1 hour pulses.

At the end of 45 minutes the control pulls a post vacuum. The combination of hold periods and post vacuum pulses continue until the aeration timer times out (the unit will pull a vacuum in the drain piping for 60 seconds prior to each post vacuum to check the integrity of the drain piping). The complete lamp comes on. The sterilizer door is unlocked, and S1 valve energizes for 60 seconds, allowing any remaining vacuum to be broken. If the sterilizer is in a post vacuum when the door is unlocked after aeration is complete, the S1 valve will remain on for 60 seconds after VS-2 is reached. At the end of 60 seconds, the Disposer is turned off.

An "OVERTEMP RESET" pushbutton, located inside the Disposer control assembly, allows qualified service personnel to reset the Disposer control when an overtemperature condition (greater than 500° F [260° C]) occurs in the catalytic bed.

### 3.5 CYCLE DESCRIPTION: MEDALLION STERILIZERS 16X16, 20X20, AND 24X36

Refer to the following drawings:

- |            |                                                |
|------------|------------------------------------------------|
| 146655-103 | Elementary diagram for Medallion 16"/20" units |
| 146655-104 | Wiring diagram for Medallion 16"/20"           |
| 146655-116 | Elementary diagram for Medallion 24x36         |
| 146655-117 | Wiring diagram for Medallion 24x36             |

The operation of the small and medium Medallion sterilizer when connected to an EO Disposer is essentially the same, with some minor differences when it comes to operation of the sterilizer vacuum pump during the water purge phase. Therefore the following description will apply to both units. Where they are different, those differences are noted.

There are also some minor differences between Envirogard and non-Envirogard equipped units. These changes are also described where applicable.

#### 3.5.1 General

An interface box assembly is mounted to the side of the Medallion sterilizer. This interface box contains a power supply, an Eaglet control board with EPROM, and an interface board. The original Medallion control, with some wiring modifications, will run the sterilizer during the conditioning and sterilize phase, but when the sterilizer enters the exhaust phase, the interface box assembly control takes over.

The Eaglet board has an EPROM which contains the exhaust cycle parameters, the interface board isolates the Eaglet board from the external loads and inputs, and the power supply supplies power to the Eaglet board.

#### 3.5.2 Cycle start

The rotary switch has a wire 837 connected to pin 45, which is hot when the rotary switch is in the off position. This 110 VAC signal is fed to the interface box assembly. When the operator starts a cycle, the rotary switch goes to the pre-vac position. This drops the 110 VAC signal to the interface box assembly. The Eaglet board, seeing the loss of 110 VAC, determines

that a cycle has started, so it then sends an output to the Disposer by wire 907. The Disposer then turns on and heats up:

#### Exhaust cycle parameters — units with EO Disposer

Summary of Cycle Differences for the various types of sterilizers (Eagle 3000 based on Rev 12)

The operation of an EO sterilizer that is connected to an EO Disposer differs from a non-Disposer equipped unit in that it turns on the Disposer and during the exhaust phase. The turn on sequence is as follows:

1. Eagle 2000 and Medallion units — Disposer ON signal occurs when cycle is started.
2. Eagle 3000 units — Disposer ON signal occurs when sterilizer starts its charge phase.

The sterilizer runs the sterilization phase in the same manner as a non-Disposer equipped unit. At the end of the exposure phase, Disposer related programming takes over. Refer to the cycle graphs.

**Disposer Ready check.** The sterilizer checks to see if the Disposer Ready signal is present from the Disposer (DISPOSER READY LED is on). If the signal is not there, the sterilizer will wait for up to 60 minutes (programmed in, not adjustable) before alarming. The alarms are as follows:

Eagle 3000 — print out "ALARM: DISPOSER NOT READY"

Eagle 2000 Printcon — print out "ALARM: DISPOSER NRDY"

Eagle 2000 IRC and Medallion units — displays "DISPOSER NOT READY" on the alarm box

The alarm can be silenced by pressing a cycle select button on the Eagle sterilizers or the silence buzzer button on Medallions.

If the Disposer Ready signal is present, the control proceeds to the Exhaust Pipe Vacuum Check.

**Exhaust Pipe Vacuum Check.** The sterilizer turns on the Vacuum signal (VACUUM LED is on), which commands the Disposer vacuum pump to turn on. A vacuum is pulled in the piping between the Disposer and sterilizer. It must reach a 23" Hg as measured by the VS23 switch within 60 seconds.

If the minimum vacuum is not met within 60 seconds, the following alarms will occur

Eagle 3000 — print out "ALARM: NO DISPOSER VACUUM"

Eagle 2000 Printcon — print out "ALARM: DISPOSER FLOW"

Eagle 2000 IRC and Medallion units — displays "DISPOSER FLOW" on the light box.

The control will try to achieve the minimum vacuum and will alarm each 60 seconds until it is reached.

If unit reaches the minimum vacuum, the control will proceed to Initial Exhaust.

**Initial Exhaust.** The sterilizer control commands the sterilizer exhaust valve S41 to open, allowing the Disposer vacuum pump to remove EO from the sterilizer chamber. The sterilizer control allows a maximum amount of time for the chamber to exhaust to a specific level, otherwise an alarm will sound. Parameters are:

Eagle 3000 — the control monitors the pressure drop for the first 40 seconds (non adjustable). The pressure in the chamber must drop a minimum of 10 mmHg or the 3000 control will alarm "ALARM: DISPOSER LINE RESTRICTED". If pressure drops more than 10 mmHg, exhaust will continue.

Then the 3000 control allows 30 minutes (adjustable in the Service Mode from 0 minutes to 59 minutes) to go from exposure pressure to 0.5 psig. If 0.5 psig is not reached, the control goes into a "TOO LONG IN EXHAUST" or TLIE alarm.

Eagle 2000 and Medallion 24x36 units — go from exposure pressure to 10" Hg within 60 minutes (not adjustable) or alarm will occur. Eagle 2000 Printcon units print "ALARM: DISP. EXHST", Eagle 2000 IRC units display "TOO LONG IN EXHAUST" on alarm box, Medallion 24x36 units display "TOO LONG IN STEP" alarm on the alarm box.

Eagle 2000 and Medallion 20x20 units — go from exposure pressure to 2" Hg within 60 minutes (not adjustable) or alarm will occur. Eagle 2000 Printcon units print "ALARM: DISP. EXHST", Eagle 2000 IRC units display "TOO LONG IN EXHAUST" on alarm box, Medallion 24x36 units display "TOO LONG IN STEP" alarm on the alarm box.

If the unit exhausts within the maximum allowable times, the control will continue to exhaust to 23" Hg.

**Continuation of first exhaust.** The sterilizer then continues to pull first post vacuum. Small sterilizers continue to have the vacuum pulled through the S41 valve, 24x36 units have a second valve, S42, which opens at 10" Hg to allow faster pull downs. The sterilizer control allows a maximum amount of time to reach VS-23 set point, otherwise an alarm will occur. Parameters are:

Eagle 3000 — go from 0.5 psig to VS-23 setting within 1 hour (adjustable in Service Mode from 00:00 to 99:99 hours:minutes) or the control prints "TOO LONG IN EVACUATION" or TLIV.

Eagle 2000 and Medallion 24x36 units — go from 10" Hg to VS-23 setting within 40 minutes (not adjustable) or Eagle 2000 Printcon prints "ALARM: DISPOSER FLOW", Eagle 2000 IRC units display "DISPOSER FLOW" on alarm box, Medallion units display "TOO LONG IN STEP" alarm on the alarm box.

Eagle 2000 and Medallion 20x20 units — go from 2" Hg to VS-23 setting within 40 minutes (not adjustable) or Eagle 2000 Printcon prints "ALARM: DISPOSER FLOW", Eagle 2000 IRC units display "DISPOSER FLOW" on alarm box, Medallion units display "TOO LONG IN STEP" alarm on the alarm box.

If unit reaches VS-23 within the maximum allowable time, the control continues to the First Air Break.

**First Air Break.** The sterilizer control turns off the Disposer vacuum pump (VACUUM LED goes off), the sterilizer S41/S42 valves deenergize, and the filtered air valve opens, allowing the chamber to break vacuum. Alarm parameters as follows:

Eagle 3000 units — The sterilizer must reach 2" Hg within 4 minutes (not adjustable) or "TOO LONG IN AIR BREAK", or TLIB, occurs.

Eagle 2000 units — current program has no too long in air break alarm.

Medallion units — The sterilizer must reach 2" Hg within 60 minutes (not adjustable) or "TOO LONG IN STEP" will be displayed on alarm box.

Once 2" Hg is reached, the control proceeds to the Second Post Vacuum.

**Second Post Vacuum.** The sterilizer control repeats the Disposer Ready check and the Exhaust Pipe Vacuum check, and if OK, turns on the Disposer vacuum pump (VACUUM LED is on), and S41 (and S42 on 24x36 units) turn on. The control allows a maximum amount of time to go from 2" to 23" Hg otherwise an alarm occurs. Alarm parameters are same as those in the Continuation of First Exhaust paragraph except that start point for timing is at 2" Hg.

Control proceeds to the Second Air Break when 23" Hg is reached.

**Second Air Break** — Operating parameters same as those for the First Air Break. Control proceeds to Air Flush when 2" Hg is reached.

**Air Flush** — The sterilizer repeats Disposer Ready check and Exhaust Pipe Vacuum check. If OK, the sterilizer filtered air valve, S41 (and S42 on 24x36 units) energize, and the Disposer vacuum pump turn on. For a 20 minute timed period (not adjustable) air is pulled through the chamber. If the Disposer Ready signal is lost at any time during air flush, the sterilizer control will allow 60 minutes for the Disposer Ready signal to return or the alarm sequence in Disposer Ready check will occur. On Eagle 3000 sterilizers, if Disposer Ready is lost any time during air flush, upon return of Disposer Ready the control will pull one post vacuum and start another 20 minute air flush.

At the end of the 20 minute air flush the control proceeds to Water Dump.

**Water Dump** — In order to remove condensed moisture from the drain piping, the control will close S41 (and S42 on 24x36 units), turn off the Disposer vacuum pump, turn on the sterilizer vacuum pump and S28 for 60 seconds (not adjustable).

**First Aerate Pulse** — Operating parameters same as the Second Post Vacuum.

**Aerate Air Break** — Operating parameters same as in the First Air Break.

**Aerate Hold Period** — Sterilizer maintains heat in chamber for 45 minutes. Vacuum level on Eagle 2000 and Medallion units stabilize at approximately 2" Hg, but may vary slightly up or down. Eagle 2000 Printcon sterilizers may print VS-2 set point more than once during hold period if the exact VS-2 setting is crossed over. Eagle 3000 units will pull down to approx 4" after reaching VS-2 set point and will maintain at 2 — 4" Hg during the hold period.

#### Special Features:

1. Check Catalyst — the sterilizer control must see a Disposer Working signal between the start of exhaust and the beginning of the 20 minute air flush, otherwise a Check Catalyst Alarm will occur.
2. Aerate Pause or Sterilize Only cycles — The only sterilizer that can be set to run a sterilize only cycle is the Eagle 3000. This is done by setting the aerate time to 00:00. Upon completion of the 20 minute air flush, the complete tone will sound and the control allows 10 minutes to unload the chamber. If the chamber is not unloaded within 10 minutes the control proceeds to Aerate.

On Eagle 2000, Medallions, and Eagle 3000's that have aerate time programmed, it is necessary to proceed with an Aerate Pause/Abort sequence. Refer to the operators manual for proper sequence.

### 3.5.3 Exhaust

The sterilizer exposure timer times out and pulses the rotary switch to the Exhaust position. The interface box assembly sees that the unit is ready to exhaust by a 110 VAC signal from the rotary switch pin 35 via wire 835. The interface box then checks for a Disposer Ready signal on wire 905. If present the interface box turns on the Disposer vacuum pump by wire 908. The control then checks the VS-23 switch (wires 823 and 825). If VS-23 actuates within 60 seconds, the interface box then opens the S41 exhaust valve (wires 803 and 828). The Disposer vacuum pump then pulls a vacuum on the chamber.

On medium sterilizers only, the control checks for actuation of the VS-10 switch (wires 885 and 825). Once VS-10 is actuated, the control then powers up S42 valve (wires 838 and 828).

This opening of S42 allows a faster exhaust rate. The S42 must delay opening until 10" Hg otherwise the Disposer catalyst will overheat. Since the same interface box is used on small sterilizers, a jumper wire is installed at the interface box VS-10 input to allow proper cycling of the control. The reason the jumper is needed is that the control must see actuation of the VS-10 input when VS-23 is reached, otherwise the control will not allow the chamber to air break.

Note that the STERILIZE light on the sterilizer control panel will remain on during the entire exhaust phase.

### 3.5.4 Air break

Once the chamber reaches VS-23 setting, the interface box turns off the Disposer vacuum pump, turns off the S41 (and S42 on medium units), and turns on the filtered air valve (wire 836).

### 3.5.5 Second post vacuum and air break

Once the chamber reaches the VS-2 setting, wires 801 and 805 signal the interface box that the chamber is at atmosphere, and the exhaust sequence above is repeated. On medium units the S42 valve now opens at the same time as the S41 valve for the remainder of the cycle.

### 3.5.6 Air flush

Once the chamber reaches the VS-2 setting, the interface box control checks for the Disposer Ready signal, and if present, turns on the Disposer vacuum pump and filtered air valve in order to perform an air flush. The timing for the 20 minute air wash is internal to the Eaglet board and is not adjustable.

### 3.5.7 Water dump

At the end of the 20 minute air flush, the interface box assembly control pulses the rotary switch to the Water Dump position via wire 850. The rotary switch wiring then powers up the sterilizer vacuum pump and S28 valve which allows water to be drawn out of the exhaust piping.

There is an input to the interface box control via wire 834 that comes from the rotary switch pin 34. This 110 VAC signal tells the interface box that Water Dump is in process and to start a 60 second timer.

The rotary switch wiring necessary to turn on the vacuum pump and S28 valve is significantly different between the small and medium units. On small units S28 is powered by pin 34 of the rotary switch and the sterilizer vacuum pump is powered by pin 48 of the rotary switch. On medium units both the S28 valve and the motor starter for the vacuum pump are powered by a connection to pin 48 of the rotary switch.

At the end of 60 seconds the interface box again pulses the rotary switch to the next position (Aerate) via wire 850.

### 3.5.8 Aerate

Once in aerate the interface box assembly control will cause one post vacuum to be pulled each hour for a 12 hour period. At the end of 12 hours the control will pull one additional post vacuum, even if it was less than 60 minutes since the last post vacuum. When this additional post vacuum is completed, the interface control will cause the rotary switch to pulse over to the next step, Complete.

While the rotary switch is in the Aerate position, the AERATING light will be lit on the control panel.

### 3.5.9 Complete

Once the control has reached the Complete position, the buzzer will sound (wire 882) for 60 seconds. The STERILE light will also come on.

### 3.5.10 Operator Unloading

On units without Envirogard, the operator must press the STOP pushbutton and the rotary switch will rotate to the OFF position, which turns off the Disposer. The load can then be removed.

On Envirogard equipped units, opening the door will cause the rotary switch to step to the OFF position, turning off the Disposer.

### 3.5.11 Short cycle

A "short cycle" has been incorporated into the Medallion program to expedite troubleshooting and "checkout". This cycle is selected by pressing and holding the "Silence Buzzer" button on the remote alarm box while powering up, and holding for 3 more seconds. The cycle times are shortened as follows:

20-MINUTE AIR WASH — 2 minutes

45-MINUTE AERATION HOLD — 2 minutes

12-HOUR AERATION — 1 hour

- a. SHORT CYCLE INDICATION: During the short cycle test mode, the Disposer Not Ready and Disposer Flow alarm lamps will blink for two seconds every five seconds as a reminder that the unit is in the short cycle mode.
- b. EXITING THE SHORT CYCLE: Turn power to the Medallion off and then on again.

## 3.6 CYCLE DESCRIPTION — SEQUENCER CONTROL

### 3.6.1 General

A sequencer control is used when there are multiple sterilizers (maximum of 3) for one Disposer. The sequencer is needed because the Disposer can handle a certain amount of EO at any one time so two sterilizers cannot exhaust to the Disposer at once.

The sterilizer interface boxes are connected to one sequencer, and then the sequencer is connected to the Disposer.

The sequencer contains an Eaglet board with EPROM, an interface board, and a power supply that supplies power to the Eaglet board. The front panel has indicator LEDs and a Auto/Manual and On/Off switch.

### 3.6.2 Cycle Operation

When a sterilizer cycle is initiated (Eagle 2000 or Medallion units) or when the sterilizer enters charge (Eagle 3000) the sequencer receives the Disposer On signal from the sterilizer, and in turn sends a signal to the Disposer to turn it on.

The sequencer programming sees this Disposer On signal from the sterilizer as the signal that the sterilizer is in cycle. The sequencer then starts a "scanning" function. During scanning the sequencer "looks" at the sterilizer(s) in cycle for 3 second every 10 seconds. On the LED display this is indicated by the STERILIZER 1, 2, or 3 LED lighting for 3 seconds every 10 seconds.

When a sterilizer is ready to exhaust, the sterilizer will send a vacuum pump turn on signal to the sequencer (VACUUM LED will be on the interface box). The sequencer sees this signal and locks on to that particular sterilizer. The sequencer will then allow that sterilizer to control the Disposer for its two post vacuums and air wash phase. At that point in the sterilizer cycle, the sterilizer control sounds a buzzer to alert the operator that the load can be removed if a sterilize only cycle is to be run. If the door is not opened, the sterilizer will progress into the aeration phase. At that point, the sequencer will drop that sterilizer from the Disposer and allow the next sterilizer to take the Disposer.


### 3.6.3 Alarms and abnormal conditions

1. Disposer not ready — if a sterilizer is in the exhaust phase and the Disposer Ready signal is lost, the LED for that particular sterilizer will flash off and on at 1 second intervals and a tone will sound in the sequencer once every second.
2. Time out condition- The sequencer counts VACUUM on signals from the sterilizer to determine where the sterilizer is at in cycle. For example if it sees two VACUUM signals, it knows the sterilizer is in the second post vacuum. The sequencer will normally allow one sterilizer to maintain control of the Disposer until it is ready to go to the aerate phase. However should a sterilizer malfunction and get stuck in some portion of the exhaust phase, the sequencer will time out and switch over to the next waiting sterilizer. The sequencer will not have any visual or audible alarms for this situation. Any alarms that would occur will be shown on the sterilizer display.

## 3.7 CYCLE DESCRIPTION — AERATORS

The aerator has two relays added to the aerator control. The first relay is powered up when the aerator is powered on. The relay contacts send a signal to the Disposer to turn ON. When the Disposer reaches operating temperature, it turns on the second relay. This second relay turns on the blower and heaters and allows the aerator to operate.

# EAGLE 3000 HI - TEMPERATURE 130°F (54°C) CYCLE WITH DISPOSER



# EAGLE 3000 LO - TEMPERATURE 100°F (38°C) CYCLE WITH DISPOSER





Figure 3-2. Cycle Graph — (Part 2 of 3)  
 Eagle 3000 Lo-Temperature  
 100°F (38°C) Cycle With Disposer

**EAGLE 3000 CHARTING CYCLE HI - TEMPERATURE 130°F (54°C) CYCLE WITH DISPOSER**  
(PREPROGRAMMED - CANNOT BE MODIFIED)



3-14/3-15  
764323-253

C- 9

**Figure 3-2. Cycle Graph — (Part 3 of 3)**  
**Eagle 3000 Charting Cycle Hi-Temperature**  
**130°F (54°C) Cycle With Disposer**  
 (Preprogrammed — Cannot Be Modified)

C- 10



**AERATE PHASE (12 HOURS) - SEE NOTE 1**






EAGLE 2000 LO - TEMPERATURE 100°F (38°C) CYCLE WITH DISPOSER

AERATE PHASE (17 HOURS) - SEE NOTE 1

PRINTCON ONLY INDICATES KEY CYCLE TRANSITION POINTS THAT ARE PRINTED DURING GAS CYCLE

- \* ALARMS  
CONDITIONING LIGHT  
A = VAC. IN 15 MIN.  
B = STEAM IN 3 MIN.  
C = GAS CHARGE IN 30 MIN.

NOTE: ALL SHADED AREAS INDICATE VALVES ENERGIZED



\* D22 - CHAMBER HEAT - SMALL (A AND B) HEATERS, MEDIUM S18 VALVE

\*\*\* D21 AUXILIARY CHAMBER HEAT SMALL (C) HEATER, MEDIUM S18 VALVE

NOTE 1 IF LONGER THAN 17 HOURS AERATE TIME IS DESIRED IT MUST BE MANUALLY TIMED

Figure 3-3. Cycle Graph — (Part 2 of 2)  
Eagle 2000 Lo-Temperature  
100°F (38°C) Cycle With Disposer

# MEDALLION GAS HI - TEMPERATURE 130°F (54°C) WITH DISPOSER





Figure 3-4. Cycle Graph — (Part 1 of 2)  
Medallion Gas Hi-Temperature  
130°F (54°C) Cycle With Disposer

**MEDALLION GAS HI - TEMPERATURE (SHORT CYCLE) 130°F (54°C) WITH DISPOSER**



**Figure 3-4. Cycle Graph — (Part 2 of 2)  
Medallion Gas Hi-Temperature  
(Short Cycle) 130°F (54°C) Cycle With Disposer**

**THIS PAGE  
INTENTIONALLY LEFT  
BLANK**

## **Section 4: General Inspection and Maintenance**

---

### **4.1 GENERAL INSPECTION AND MAINTENANCE**

The following procedures are recommended on a routine basis by the customer. A preventive maintenance guide follows this section and should be performed by a qualified service technician for optimum operation.

1. Routine inspection
  - a. Inspect Disposer door gasket seal for tears
  - b. Inspect Disposer air inlet for debris
  - c. Verify Disposer indicating lights are functioning properly
2. Quarterly
  - a. Inspect Disposer air inlet filter; replace if required. Removal of the filter is necessary for inspection as dirt collects on the inside. (Note: Disposer installations in locations where dust or dirt is a problem may require more frequent inspections and replacements of the filter.)

| Table 4-2. Preventive Maintenance Guide                                   |   |   |   |     |   |   |
|---------------------------------------------------------------------------|---|---|---|-----|---|---|
|                                                                           | 1 | 2 | 3 | 4   | 5 | 6 |
| <b>1.0 PREPARATION FOR PREVENTIVE MAINTENANCE</b>                         |   |   |   |     |   |   |
| 1.1 Discuss operation with department personnel .....                     | x | x | x | x   | x | x |
| 1.2 Review sterilizer printouts/charts for any problems .....             | x | x | x | x   | x | x |
| <b>2.0 MECHANICAL COMPONENTS</b>                                          |   |   |   |     |   |   |
| 2.1 Air filter—inspect for debris, replace each 6 months minimum .....    |   | x |   | x   |   | x |
| 2.2 Air flow switch-PS-1, PS-2—confirm setting .....                      |   |   |   |     | x |   |
| 2.3 Gas inlet solenoid valve FSV-1, rebuild .....                         |   |   |   |     | x |   |
| 2.4 Leak check piping between sterilizer and Disposer .....               | x |   |   |     |   |   |
| 2.5 Vacuum pump, rebuild .....                                            |   |   |   | A/R |   |   |
| <b>3.0 ELECTRICAL COMPONENTS</b>                                          |   |   |   |     |   |   |
| 3.1 Control relays—remove and visually inspect, replace as required ..... |   | x |   | x   |   | x |
| 3.2 Fuses—inspect for signs of overheating, replace as required .....     | x | x | x | x   | x | x |
| 3.3 Heaters—measure voltage and current .....                             |   | x |   | x   |   | x |
| 3.4 Vacuum Pump—measure voltage and current .....                         |   | x |   | x   |   | x |
| 3.5 Verify Control and Overtemp T-stat settings .....                     | x |   |   |     |   |   |
| <b>4.0 FUNCTIONAL</b>                                                     |   |   |   |     |   |   |
| 4.1 *Perform EO analysis of stack output .....                            |   |   |   | A/R |   |   |
| 4.2 Verify Analogic Controller set points .....                           |   |   | x |     |   |   |
| 4.3 Confirm operation of Disposer indicating lights .....                 | x | x | x | x   | x | x |
| <b>5.0 STERILIZER COMPONENTS</b>                                          |   |   |   |     |   |   |
| 5.1 Exhaust line solenoid—rebuild .....                                   |   |   |   |     | x |   |
| 5.2 Exhaust line strainer—inspect and clean .....                         |   |   | x |     |   | x |
| 5.3 Interconnect/Sequencer box—inspect for loose wires .....              | x | x | x | x   | x | x |
| 5.4 Confirm proper operation of indicator lights .....                    | x | x | x | x   | x | x |

\*Check with local air quality management department for frequency required for stack testing.  
A/R = As required.

Section 5: Troubleshooting

5.1 TROUBLESHOOTING CHART

Operational faults on units with an EO Disposer will not appear until the sterilizer enters the Exhaust phase. Faults that are present during Condition, Charge, and Sterilize portions of an EO sterilizer cycle do not occur with the EO Disposer, or related Disposer components.

Table 5-1. Troubleshooting Chart

| Status                                                                                | Trouble                                                                                                                                                                                                                                                                             | Cause/Correction                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. Fault — Disposer Not Ready condition                                               | <p>Eagle 3000 units — prints out "ALARM: DISPOSER NOT READY"</p> <p>Eagle 2000 Printcon — prints out "ALARM: DISPOSER NRDY"</p> <p>Eagle 2000 IRC — displays "NOT READY" on the remote alarm box</p> <p>Medallion units — displays "DISPOSER NOT READY" on the remote alarm box</p> | <p>Situation — The sterilizer control requires the Disposer to achieve Ready status within 60 minutes of starting the Exhaust phase. If the Disposer does not receive the signal, the alarm occurs.</p> <p>Causes:</p> <ol style="list-style-type: none"> <li>1. Disposer power switch OFF — this is indicated by the "POWER" LED on the sterilizer interface box being off</li> <li>2. Insufficient air flow through the Disposer — this is indicated by the "LOW AIR FLOW" LED on the sterilizer interface box being lit.</li> </ol> <p>Causes of low air flow can be any of the following:</p> <ol style="list-style-type: none"> <li>a. plugged air filter in Disposer</li> <li>b. faulty fan in Disposer</li> <li>c. faulty air pressure switch PS-1 or PS-2 in Disposer</li> <li>d. excess back pressure in Disposer exhaust duct</li> </ol> <ol style="list-style-type: none"> <li>3. Disposer has not reached minimum operating temperature (280° F). Causes of problem can be: <ol style="list-style-type: none"> <li>a. ambient air temperature too low — minimum air temp allowable is 50° F</li> <li>b. faulty incoming 3-phase power</li> <li>c. faulty heaters in Disposer</li> <li>d. faulty temperature controller or relays in Disposer</li> <li>e. fault in communication lines between Disposer and sterilizer</li> </ol> </li> <li>4. Disposer has exceeded maximum operating temperature (450° F). Causes of problem can be: <ol style="list-style-type: none"> <li>a. faulty catalytic bed</li> <li>b. faulty heater controller or relays in Disposer</li> <li>c. faulty exhaust S41 or S42 valve in sterilizer</li> </ol> </li> </ol> |
| 2. Fault — Insufficient vacuum level in exhaust line during exhaust pipe vacuum check | <p>Eagle 3000 units — prints "ALARM: NO DISPOSER VACUUM"</p> <p>Eagle 2000 Printcon — prints "ALARM: DISPOSER FLOW"</p> <p>Eagle 2000 IRC — displays "NO AIR FLOW" on the remote alarm box.</p> <p>Medallion units — displays "NO DISPOSER AIR FLOW" on the remote alarm box</p>    | <p>Situation — The sterilizer control requires that a minimum vacuum level is reached (VS-23 setting) during the exhaust pipe vacuum check within 60 seconds of starting the vacuum check. The sterilizer control continues to try to achieve the VS-23 vacuum setting even if the alarm has occurred. This alarm automatically silences if reaching the VS-23 setting.</p> <p>This alarm also occurs on Medallion sterilizers if one of the two additional conditions do not occur:</p> <ol style="list-style-type: none"> <li>a. The VS-23 switch must open back up within 10 seconds of the completion of the pipe check</li> <li>b. The 2-inch vac alarm switch must be open at the time VS-23 is actuated</li> </ol> <p><i>Note for Medallion controls only, if either of the above conditions occur, the alarm buzzer will not self clear and continues buzzing until silenced by the operator.</i></p> <p>Causes:</p> <ol style="list-style-type: none"> <li>1. Leak in exhaust pipe between sterilizer and Disposer</li> <li>2. Incorrect VS-23 setting — valve in exhaust pipe closed</li> <li>3. Faulty VS-23 switch</li> <li>4. Faulty vacuum pump in Disposer</li> <li>5. Faulty FSV-1 solenoid valve in Disposer</li> <li>6. Faulty 2" vac alarm switch (Medallion units only)</li> </ol>                                                                                                                                                                                                                                                                                                                                                       |

5-2  
784323-253

D- 9

Table 5-1. Troubleshooting Chart (continued)

| Status                                                                                                                                | Trouble                                                                                                                                                                                                                                                                     | Cause/Correction                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 3. Fault — Initial sterilizer exhaust is too slow (applies to Eagle 3000 sterilizers only)                                            | 3000 control prints "ALARM: DISPOSER LINE RESTRICTED"                                                                                                                                                                                                                       | <p>Situation — The 3000 control monitors the pressure drop during the first 40 seconds of exhaust. The pressure must drop a minimum of 10 mm Hg or the alarm will occur.</p> <p>Causes:</p> <ol style="list-style-type: none"> <li>1. Faulty S41 exhaust line solenoid valve</li> <li>2. control requires calibration</li> </ol>                                                                                                                                                                                                                                                                                                                                          |
| 4. Fault — Sterilizer takes too long to exhaust to approx. atmospheric pressure.                                                      | <p>Eagle 3000 units — prints "TOO LONG IN EXHAUST"</p> <p>Eagle 2000 Printcon — prints "ALARM: DISP. EXHST"</p> <p>Eagle 2000 IRC — displays "TOO LONG IN EXHAUST" on the remote alarm box</p> <p>Medallion units — displays "TOO LONG IN STEP" on the remote alarm box</p> | <p>Situation — The sterilizer control allows a maximum amount of time to go from exposure pressure to a set vacuum level.</p> <p>Causes:</p> <ol style="list-style-type: none"> <li>1. Faulty sterilizer S41 valve</li> <li>2. Faulty check valve on drain line of sterilizer's original plumbing</li> <li>3. Sterilizer exhaust orifice plugged</li> <li>4. Sterilizer exhaust strainer plugged</li> <li>5. Sterilizer has lost the Disposer Ready signal</li> </ol>                                                                                                                                                                                                     |
| 5. Fault — Sterilizer takes too long to exhaust to a deep vacuum setting (VS-23) during post vacuum exhaust pulses or aeration pulses | <p>Eagle 3000 units — prints "TOO LONG IN EVACUATION"</p> <p>Eagle 2000 Printcon — prints "ALARM: DISPOSER FLOW"</p> <p>Eagle 2000 IRC — displays "NO AIR FLOW" on the remote alarm box</p> <p>Medallion units — displays "TOO LONG IN STEP" on the remote alarm box</p>    | <p>Situation — The sterilizer control allows a maximum time for the sterilizer to reach the VS-23 setting of 23" Hg.</p> <p>Causes:</p> <ol style="list-style-type: none"> <li>1. Leak in exhaust pipe between sterilizer and Disposer</li> <li>2. Incorrect VS-23 setting</li> <li>3. Faulty VS-23 switch</li> <li>4. Faulty vacuum pump in Disposer</li> <li>5. Faulty S41 (or S42 on 24x36 units) exhaust valve</li> <li>6. Faulty drain line check valve on original sterilizer plumbing</li> <li>7. Water is present in exhaust line (source of water usually due to faulty sterilizer drain line check valve or vacuum pump water inlet solenoid valve).</li> </ol> |

5-3  
784323-253

D- 10

Table 5-1. Troubleshooting Chart (continued)

| Status                                                                        | Trouble                                                                                                                                                                                                                                                                                                                         | Cause/Correction                                                                                                                                                                                                                                                                            |
|-------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6. Fault — Sterilizer takes too long to air break.                            | Eagle 3000 — printer prints "TOO LONG IN AIR BREAK"<br>Eagle 2000 — no indication, printed or audible<br>Medallion — displays "TOO LONG IN STEP" on remote alarm box                                                                                                                                                            | Situation — Eagle 3000 units allow 4 minutes to reach 2 in Hg., Medallion units allow 60 minutes to reach 2 in Hg. Eagle 2000 units do not have a "TOO LONG IN AIR BREAK" alarm.<br>Causes:<br>1. Faulty filtered air valve<br>2. Plugged sterilizer filter on the filtered air valve inlet |
| 7. Fault — Disposer does not reach minimum working temperature during exhaust | Eagle 3000 — printer prints "ALARM: CHECK DISPOSER CATALYST" and "CATALYST" LED is lit on interface box<br>Eagle 2000 Printcon — printer prints "ALARM: CHECK CATALYST" and "CATALYST" LED is lit on interface box<br>Medallion — displays "CHECK CATALYTIC BED" on remote alarm box and "CATALYST" LED is lit on interface box | Situation — Disposer must reach a minimum of 320°F during sterilizer exhaust phase.<br>Causes:<br>1. Air or Refrigerant 12 is being run through sterilizer<br>Temperature rise relies on catalytic reaction with EO for temperature rise<br>2. Faulty catalyst                              |
| 8. Fault — Disposer has insufficient air flow                                 | Indication — "LOW AIR FLOW" LED is lit on sterilizer interface box                                                                                                                                                                                                                                                              | Situation — Pressure switches inside Disposer indicate insufficient air flow through the Disposer.<br>Causes:<br>1. Air filter inside Disposer is plugged<br>2. Disposer fan not operating<br>3. Back pressure in exhaust duct<br>4. Pressure switches in Disposer require adjustment       |

Table 5-1. Troubleshooting Chart (continued)

| Status                                                                 | Trouble                                                                                       | Cause/Correction                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 9. Fault — Disposer has exceeded maximum operating temperature (500°F) | "OVERTEMP" LED is lit on sterilizer interface box and "OVERTEMP" light is lit inside Disposer | Situation — Disposer goes to an overtemp condition when temperature exceeds 500° F. A latch relay in the Disposer maintains the alarm condition until manual reset button is pressed.<br>Causes:<br>1. VS-10 vacuum switch incorrectly set (applies to 24x36 IRC units only)<br>2. Faulty S42 sterilizer exhaust valve (applies to 24x36 units only)<br>3. Faulty catalyst bed<br>4. Sterilizer charged to higher than 8.5 psig<br>5. Coils on S41 and S42 valves are reversed (applies to 24 x 36 units only) |

THIS PAGE  
INTENTIONALLY LEFT  
BLANK

## Section 6: Component Replacement

### 6.1 ADJUSTMENT PROCEDURES

Tools required for this section are:

1. Digital potentiometer set for type "J" thermocouple
2. Calibrated compound gauge
3. Differential pressure gauge (P-764323-802) and gauge bit (P-764323-803) capable of reading 0 to 4 inches of water pressure
4. 50cc syringe (without needle)

#### 6.1.1 Sterilizer Components

1. VS-23 vacuum switch—Install compound gauge in plugged tee adjacent to the vacuum switch. When sterilizer is exhausting to the disposer, monitor the vacuum level in the piping. Adjust switch at 23 in Hg.
2. VS-10 vacuum switch (24x36 IRC units only)—Install compound gauge in plugged tee adjacent to the vacuum switch. When sterilizer is exhausting to the disposer, monitor the vacuum level in the piping. Adjust switch at 10 in Hg.

#### 6.1.2 Disposer Components

1. Heater Control Thermostat—Observe temperature on analogic controller. With unit stabilized at operating temperature (290° F [143° C]), compare the temperature on the thermostat dial with the digital readout. Adjust thermostat dial if they do not match.
2. Heater Overtemp Thermostat—Temporarily set the control thermostat to operate at 400° F (204° C) and the overtemp thermostat at 375° F (190° C). When temperature reaches the overtemp thermostat setting of 375° F (190° C), the heaters should turn off. Verify the dial setting of the thermostat and adjust if necessary. At the completion of this adjustment reset the control thermostat to 290° F (143° C) and the overtemp thermostat to 550° F (270° C).
3. Pressure Switches—Disposer has two pressure switches, PS1 and PS2. PS1 is to be set at 1.5 inches water, and PS2 is to be set at 1.8 inches water.

In the following procedures the syringe is used to pressurize the switch. The set point of the switch can be set by monitoring the pressure gauge and the Disposer control.

- a. PS1 measures the differential pressure across the catalyst bed. In the event of a plugged catalytic bed this switch will activate. Disconnect the 1/4 inch tube at the 90° fitting located at the top side of the heater box (See Figure 6-1). With the tube disconnected, the switch will deactivate and the Low Air Flow light in the Disposer will come on. Slide the hose onto the 1/4 inch tube. Connect the gauge, tee, and syringe. Slowly inject air into the line and watch the pressure gauge. The switch should activate at 1.5 inches. When the switch activates, the Low Air Flow light in the Disposer will turn off. Reconnect the 1/4 inch tube when completed.
- b. PS2 measures the differential between the outside of the cabinet and the inside. In the event of a plugged air filter, this switch will activate. This checks for operation of the filter element. Connect the hose to the elbow fitting at the back of the Disposer. Pressurize the line until the switch activates. The switch should activate at 1.8 inches. When the switch activates, the Low Air Flow light in the Disposer will come on.

### 6.2 COMPONENT REPLACEMENT PROCEDURES

#### 6.2.1 Prefilter (Figure 6-2)

The prefilter should be changed at least every six months. Replace filter as follows:

1. Be sure main disconnect switch is OFF.
2. Loosen door clamps and open door.
3. Remove clamp and prefilter.
4. Install new filter P/N 764323-032.
5. Replace clamp and tighten.
6. Close door and tighten door clamps.

#### 6.2.2 Catalytic Cell (Figure 6-3)

The catalytic cell must be replaced when it is no longer effective for the removal of EO, as determined by EO stack analysis. Replace as follows:

1. Be sure main disconnect switch is OFF.
2. Loosen door clamps and open door.






Figure 6-1. Disposer Pressure Switch Adjustment

6-2  
704323-253

E- 1




Figure 6-2. Prefilter Replacement

3. Loosen clamps on hose (between heater and catalytic cell transition plenum) and slide clamps and hose down.
4. Disconnect EO feed to plenum at union.
5. Check that no EO vapor is present in the duct using standard site practices.
6. Remove 3/8" bolts from perimeter of transition plenum and remove plenum.
7. Remove four 3/8" hex nuts and two clamping bars holding catalytic cell in place and remove catalytic cell.

**NOTE:** Catalytic cell weighs approximately 80 lbs. It is recommended that two people lift and remove cell or use other appropriate lifting methods.

8. Inspect interior for dirt or debris and remove.
9. Inspect new catalytic cell P764323-196. Inspect gasket for damage; repair and/or replace if necessary.
10. Install new catalytic cell. Be sure cell is seated squarely in lower housing. Install four 3/8" hex nuts and two clamping bars for catalytic cell holddown.

11. Install new gasket, if needed, around perimeter of transition plenum. Gasket and adhesive suitable for high temperature must be used.
12. Continue installation in reverse order of disassembly. Tighten all bolts and nuts evenly.

**NOTE:** Be sure all connections are tight and leak free.

#### 6.2.3. Power Switch Cable Adjustment (Figure 6-4)

1. Confirm that supplying power to the device has been turned off.
2. Move circuit breaker handle to ON position.
3. Place end of cable through hole in tab of the sliding plate on circuit breaker operator.
4. Move operating handle to its maximum ON position and replace the cable mounting threads into the slot on the fixed tab plate. Adjust the mounting nuts so that the nut on the cable rod end touches the sliding plate tab.
5. Hold operating handle in ON position and place spring over rod. Attach spring adjuster and tighten it until it touches the spring. **DO NOT OVERTIGHTEN THE SPRING ADJUSTMENT NUT.**
6. Check that the circuit breaker turns OFF and ON by moving the operating handle. If the breaker does not turn ON, loosen the cable mounting nuts at the fixed plate tab and move the cable forward until the breaker turns ON. Retighten mounting nuts.
7. Trip the circuit breaker by pushing the "trip" button.
8. Check that the breaker resets by moving the operating handle from ON to OFF and back to ON. If the breaker resets, tighten the spring adjuster one additional turn. Attach lockwasher and #10-32 nut to cable rod end and tighten to 32 in-lbs.
9. If the circuit breaker does not reset, tighten the spring adjuster one turn and repeat steps 2 and 3. Continue this until breaker resets, then tighten the adjuster one additional turn and secure with lockwasher and nut.

6-3  
704323-253

E- 2




Figure 6-3. Catalytic Cell

6-4  
704323-253

E- 3




Figure 6-4. Cable Adjustment

6-5  
704323-253

E- 4

### 6.3 DISPOSER BLOWER AIR FLOW ADJUSTMENT

Tools required for this section are:

1. Air flow meter (thermo-anemometer model 8565 used in AMSCO EO analysis or equivalent) capable of measuring between 1000 and 2000 ft/min.

#### Blower adjustment

1. Blower output — The air flow through the Disposer is factory set at  $120 \pm 10$  cfm. After installation the setting has to be confirmed and adjusted as necessary to compensate for any restriction that may be present in the exhaust duct. To adjust properly the final configuration exhaust duct must be installed. Adjust as follows:

- a. Remove air inlet cover on Disposer. If aerators are connected to air inlet, these must be disconnected also.
- b. Close cabinet door and close tightly with latches on door.
- c. Using air flow meter, measure air flow at Disposer inlet. Take four measurements approximately 1 inch in from edge (see Figure 6-6). See Table 6-1 for air flow rate to cfm conversion.




Figure 6-6. Blower Adjustment

- d. If air flow is  $120 \pm 10$  cfm, unit is within spec. and no adjustment is necessary. If not, alter the air flow rate by adjusting the baffle plate located just below the blower outlet. It will be necessary to remove the 4 bolts located around the blower outlet perimeter. Slide the baffle in or out (it is notched to allow a step adjustment). Tighten bolts and retest.

Table 6-1. Air Flow Conversion Chart

| Ft/Min | CFM |
|--------|-----|
| 1333   | 110 |
| 1345   | 111 |
| 1357   | 112 |
| 1369   | 113 |
| 1381   | 114 |
| 1393   | 115 |
| 1406   | 116 |
| 1418   | 117 |
| 1430   | 118 |
| 1442   | 119 |
| 1454   | 120 |
| 1466   | 121 |
| 1478   | 122 |
| 1490   | 123 |
| 1503   | 124 |
| 1515   | 125 |
| 1527   | 126 |
| 1539   | 127 |
| 1551   | 128 |
| 1563   | 129 |
| 1575   | 130 |

## SERVICE BULLETIN

Product: EO DISPOSER

Bulletin No.: 92020-I

Section: STERILIZERS

Supplement No.: ST-201

Approved: *[Signature]*

Date: 9/1/92

Subject: WATER TRAP KIT FOR EXHAUST LINE

**SITUATION:** EO Disposer kits are now being shipped with a water trap that is installed in the exhaust line. Kits as shipped previously did not have the water trap.

**DESCRIPTIVE CAUSE:** The Disposer vacuum pump is susceptible to damage if water is present in the exhaust line. Water can enter the exhaust line from a fault in the sterilizer or fault in the steam supply. The water trap will collect the water and prevent it from flowing to the Disposer.

**CORRECTIVE ACTION:** Order kit P755716-132 for each sterilizer connected to an EO Disposer. The kit consists of the water trap, miscellaneous pipe fittings, and a shut-off valve. Parts that need to be obtained in the field are adaptor fittings to connect the water trap's 3/4" NPT female inlet and outlet to the exhaust line. In addition, a 3/4" brass tee and a 3/4" x 1/4" bushing may be needed on some earlier Disposer installations to allow the water trap drain line to be adapted to the sterilizer drain piping.

The water trap self drains small amounts of water into the sterilizer drain piping, during the water purge phase of the cycle. If it should fill up and prevent a vacuum from being pulled in the sterilizer, then a fault exists which must be corrected. Check the following:

1. faulty sterilizer drain line check valve
2. faulty sterilizer vacuum pump water inlet solenoid valve
3. improper sterilizer steam conditioning setting
4. poor quality steam/inadequate steam trapping on steam line
5. improperly leveled sterilizer

Ref: Field Service Engineering/J. Billey



THIS CLASSIFIED INFORMATION FOR USE BY AMSCO PERSONNEL ONLY AND IS THE LATEST AVAILABLE AT DATE OF ISSUE. AMSCO ENGINEERING DIVISION RESERVES THE RIGHT TO CHANGE THIS DATA WITHOUT INCURRING OBLIGATIONS OTHER THAN ROUTINE DISTRIBUTION OF REVISED SERVICE BULLETIN. COPYRIGHT 1990 BY AMERICAN STERILIZER COMPANY, ERIE, PENNSYLVANIA, U.S.A.

Printed in U.S.A.

6-6  
764323-253

E- 5

E- 6

4-10-7 (2/89)

2. If the vacuum switch is present, then the program revision must be revision 11 or later. If the sterilizer has an earlier revision chip, order new EPROM P-93910-634. EPROMS in stock will be Revision 12 or newer and will be compatible with the VS23 switch.
3. If the interconnect kit does not have a VS23 switch, and the program is Revision 10 or earlier, then the kit can be installed as is.
4. If the interconnect kit does not have a VS23 switch, and the program is Revision 11 or later, there are two options:
  - a. Obtain a Revision 10 EPROM through Field Engineering
  - b. Install VS23 switch in the piping. This is the preferred option in that the later revision program chips have enhancements that are beneficial to both the customer and for servicing. Order P-755716-164. This kit includes the VS23 vacuum switch, wiring cable, and miscellaneous hardware needed to install the switch.

#### Existing Installations:

To allow installation of Revision 11 or newer software in an Eagle 3000 EO sterilizer that does not have VS23 in the piping, order:

- a. EPROM P-93910-634
- b. Kit P-755716-164. This kit includes the VS23 vacuum switch, wiring cable and miscellaneous hardware needed to install the switch.

## SERVICE BULLETIN

Product: EAGLE 3000 EO GAS STERILIZERS WITH DISPOSER

Bulletin No.: 92021-I

Section: STERILIZERS

Supplement No.: ST-202

Approved:

Date: 9/1/92

Subject:

COMPATIBILITY OF PROGRAM REVISION WITH DISPOSER KIT

### SITUATION

A design change has been made on the interconnect kit that allows the Eagle 3000 EO sterilizer to be connected to the EO Disposer. This design change also involves a change in the revision level of the program in the sterilizer control.

### DESCRIPTIVE CAUSE:

The design interconnect kit as originally designed used the charge line transducer to measure vacuum levels when the Disposer vacuum pump was pulling vacuum. The current design released in the 2nd quarter of 1992 uses a standard vacuum switch, labeled VS23, to perform this function. The VS23 vacuum switch, if present, is located in the exhaust piping adjacent to the S28 water dump valve.

In order for the control to function properly with the new vacuum switch, a program change was made beginning with Revision 11. Therefore, on Disposer kits that have a VS23 vacuum switch in the interconnect piping, Revision 11 or later program chips are required.

On units that do not have VS23 vacuum switch, programs Revision 10 and earlier can only be used. Therefore existing Eagle 3000 EO sterilizers with Disposer cannot update to Revision 11 or later EPROMS unless the modification described below is installed.

### CORRECTIVE ACTION:

#### New Installations:

1. When installing a Disposer kit on a 3000 EO sterilizer verify if the interconnect kit (P-146655-710 for 20x20 units, P-146655-711 for 24x36 units) has a VS23 vacuum switch (P-83866-001). Then check the program revision level on the EPROM on the control board.

continued...

Ref. Field Service Engineering/J. Bailey



THIS CLASSIFIED INFORMATION FOR USE BY AMSCO PERSONNEL ONLY AND IS THE LATEST AVAILABLE AT DATE OF ISSUE. AMSCO ENGINEERING DIVISION RESERVES THE RIGHT TO CHANGE THIS DATA WITHOUT INCURRING OBLIGATIONS OTHER THAN ROUTINE DISTRIBUTION OF REVISED SERVICE BULLETIN. COPYRIGHT 1992 BY AMERICAN STERILIZER COMPANY, ERIE, PENNSYLVANIA, U.S.A.

Printed in U.S.A.

## SERVICE BULLETIN

Product: EO DISPOSER

Bulletin No.: 92017-I

Section: STERILIZERS

Supplement No.: ST-200

Approved: *[Signature]*

Date: 6/15/92

Subject: ENVIRONMENTAL VALIDATION TESTING

**SITUATION:** EO analysis of the Disposer exhaust indicates less than advertised efficiency (99.9% destruction of EO). This EO testing is not available from AMSCO at this time, but qualified local contractors can perform this service.

**CORRECTIVE ACTION:** AMSCO has a contractual agreement with the manufacturer of the EO Disposer, Donaldson, Inc. This agreement requires that Donaldson must be consulted, should the EO Disposer not pass the initial EO test.

In the event the Disposer does not pass EO testing, contact both AMSCO Field Engineering in Erie, and Donaldson, Inc. The contact at Donaldson is the Manager of Engineering, High Purity Products (612/887-3417). As of this writing, the manager is Tim Walsh.

Do not replace any component on the EO Disposer until after making the contact with Field Engineering and Donaldson.

Ref: Field Service Engineering/J. Biley (x2962)



THIS CLASSIFIED INFORMATION FOR USE BY AMSCO PERSONNEL ONLY AND IS THE LATEST AVAILABLE AT DATE OF ISSUE. AMSCO ENGINEERING DIVISION RESERVES THE RIGHT TO CHANGE THIS DATA WITHOUT INCURRING OBLIGATIONS OTHER THAN ROUTINE DISTRIBUTION OF REVISED SERVICE BULLETIN.  
COPYRIGHT 1990 BY AMERICAN STERILIZER COMPANY, ERIE, PENNSYLVANIA, U.S.A.

Printed in U.S.A.

4-10-7 (2/89)

## SECTION 7: EXPLODED VIEWS AND PARTS LISTS

Assemblies and components of the Disposer are illustrated and identified on the following pages. The part number, the description and the quantity required for each usage is given.

### HOW TO USE THE ILLUSTRATED PARTS BREAKDOWN:

1. Determine the function and application of the part required. Turn to the List of Illustrations and select the most appropriate title. Note the illustration page number.
2. Turn to the page indicated and locate the desired part on the illustration.
3. From the illustration, obtain the index number assigned to the part desired. Refer to the same index number in the PARTS LIST accompanying the illustration to find the AMSCO PART NUMBER and description of specific information regarding the part.

7-1  
784323-253




Figure 7-1. EO Disposer Assembly (Part 1 of 3)

7-2  
764323-253

E- 11

| FIG. & INDEX NO. | PART NUMBER | SVC | DESCRIPTION                                      | UNITS PER ASSEMBLY |
|------------------|-------------|-----|--------------------------------------------------|--------------------|
| 7-1-             |             |     | ETO DISPOSER ASSEMBLY (Part 1 of 3) .....        | X                  |
| 1                |             |     | HOUSING WELDMENT .....                           | 1                  |
| 2                | P 764323    | 727 | DOOR GASKET (12-7" REQ'D.) .....                 | 1                  |
| 3                |             |     | OUTLET GASKET .....                              | 1                  |
| 4                |             |     | ADHESIVE, SILICONE RUBBER .....                  | A/R                |
| 5-12             |             |     | See Part 3 of 3                                  |                    |
| 13               |             |     | WELDMENT, LOWER UNIT .....                       | 1                  |
| 14               |             |     | FIBREX, 2" THK. ....                             | A/R                |
| 15               |             |     | FIBREX, 1" THK. ....                             | A/R                |
| 16               | P 764323    | 196 | FILTER ASSEMBLY (Catalyst Cell) .....            | 1                  |
| 17               |             |     | ROD, THREADED .....                              | 4                  |
| 18               |             |     | BAR, CLAMP .....                                 | 2                  |
| 19               |             |     | DUCT, TRANSITION ASSEMBLY .....                  | 1                  |
| 20               | P 764323    | 729 | BLANKET, TRANSITION .....                        | 1                  |
| 21               |             |     | HOUSING, HEATER .....                            | 1                  |
| 22               |             |     | CROSS BRACE, LOWER FRONT .....                   | 1                  |
| 23               |             |     | CROSS BRACE, UPPER FRONT .....                   | 1                  |
| 24               |             |     | CROSS BRACE, LOWER BACK .....                    | 1                  |
| 25               |             |     | CROSS BRACE, UPPER BACK .....                    | 1                  |
| 26               |             |     | CROSS BRACE, PUMP .....                          | 1                  |
| 27               | P 764323    | 728 | BOOT, RUBBER 4-1/2" LG. ....                     | 1                  |
| 28               |             |     | CLAMP, HOSE 6" .....                             | 2                  |
| 29               | P 764324    | 739 | HEATER ASSEMBLY, (See Figure 7-2) .....          | 1                  |
| 30               |             |     | WELDMENT, BLOWER PLATE .....                     | 1                  |
| 31               |             |     | PLATE, ORIFICE .....                             | 1                  |
| 32               |             |     | U-PLATE, ORIFICE .....                           | 1                  |
| 33               | P 764323    | 730 | GASKET, ORIFICE .....                            | 2                  |
| 34               | P 764323    | 731 | BLOWER .....                                     | 1                  |
| 35               | P 764323    | 732 | BLOWER STRIP .....                               | 2                  |
| 36               |             |     | ELECTRICAL PANEL ASSEMBLY (See Figure 7-3) ..... | 1                  |
| 37               | P 764323    | 032 | PREFILTER .....                                  | 1                  |
| 38               |             |     | CAP, INLET .....                                 | 1                  |
| 39               |             |     | CLAMP, HOSE 4" .....                             | 1                  |
| 40               |             |     | INTERFACE, OUTLET .....                          | 1                  |
| 41               |             |     | CLAMP ASSEMBLY .....                             | 1                  |
| 42               |             |     | ETO PLUMBING (See Figure 7-4) .....              | 1                  |
| 43               |             |     | NOT USED                                         |                    |
| 44               | P 764323    | 248 | PUMP .....                                       | 1                  |
|                  | P 764323    | 249 | • KIT, Repair .....                              | 1                  |
|                  | P 764324    | 713 | • MOTOR .....                                    | 1                  |
|                  | P 764324    | 714 | • • CAPACITOR .....                              | 1                  |
| 45               |             |     | PLATE, PUMP .....                                | 1                  |
| 46               |             |     | PAD, PUMP .....                                  | 4                  |
| 47               |             |     | SNAP BUSHING, SUPPORT ANGLE .....                | 1                  |
| 48               |             |     | LABEL, DONALDSON .....                           | 1                  |
| 49               |             |     | LABEL, CAUTION .....                             | 1                  |
| 50               |             |     | LABEL, CAUTION, HOT SURFACE .....                | 1                  |
| 51               |             |     | LABEL, ETO INLET CONNECTION .....                | 1                  |
| 52               |             |     | LABEL, EXHAUST CONNECTION .....                  | 1                  |
| 53               |             |     | LABEL, READ AND UNDERSTAND .....                 | 1                  |
| 54               |             |     | LUG, LIFTING .....                               | 2                  |

(Continued on Part 2 of 3)

7-3  
764323-253

E- 12






Figure 7-1. EO Disposer Assembly (Part 2 of 3)

| FIG. & INDEX NO.           | PART NUMBER | SVC | DESCRIPTION                           | UNITS PER ASSEMBLY |
|----------------------------|-------------|-----|---------------------------------------|--------------------|
| 7-1-                       |             |     | ETO DISPOSER ASSEMBLY (Part 2 of 3)   | X                  |
| 55                         |             |     | BOLT, 3/8-16 HEX HD. x 1              | 24                 |
| 56                         |             |     | NUT, 3/8-16 HEX                       | 28                 |
| 57                         |             |     | WASHER, FLAT, 3/8                     | 38                 |
| 58                         |             |     | LOCKWASHER, 3/8                       | 36                 |
| 59                         |             |     | BOLT, 5/16-18 HEX HD. x 1             | 4                  |
| 60                         |             |     | NUT, 5/16-18 HEX, LOCKING, NYLOK      | 4                  |
| 61                         |             |     | WASHER, FLAT, 5/16                    | 8                  |
| 62                         |             |     | LOCKWASHER, 5/16                      | 10                 |
| 63                         |             |     | BOLT, 5/16-18 HEX HD. x 3/4           | 6                  |
| 64                         |             |     | BOLT, 1/4-20 HEX HD. x 1              | 31                 |
| 65                         |             |     | NUT, 1/4-20 HEX                       | 3                  |
| 66                         |             |     | WASHER, FLAT, 1/4                     | 34                 |
| 67                         |             |     | WASHER, LOCK, 1/4                     | 31                 |
| 68                         |             |     | NUT, HEX, 1/2-13                      | 2                  |
| 69                         |             |     | WASHER, LOCK, 1/2                     | 2                  |
| 70                         |             |     | SCREW, NO. 10-32 x 1/2 THREAD FORMING | 1                  |
| 71                         |             |     | SCREW, NO. 8-32 x 1/2 ROUND HD        | 2                  |
| 72                         |             |     | LOCKWASHER, #8                        | 2                  |
| 73                         |             |     | NUT, #8-32                            | 2                  |
| 74                         |             |     | PLUG, PIPE, 1/4 NPT STD., SQ. HEAD    | 1                  |
| 75                         |             |     | BOLT, 5/16" x 1 SELF-LOCKING, HEX HD. | 16                 |
| 76                         |             |     | LABEL, AMSCO                          | 1                  |
| 77                         |             |     | LABEL, WARNING, DOOR                  | 1                  |
| 78                         |             |     | LABEL, WARNING, AIR CLEANER           | 1                  |
| 79                         |             |     | CAP, PLUG, ELEC.                      | 2                  |
| 80                         |             |     | CAP, PLUG, ETO INLET                  | 1                  |
| 81                         |             |     | CAP, PLUG, INLET                      | 1                  |
| 82                         |             |     | CAP, PLUG, OUTLET                     | 1                  |
| 83                         |             |     | BUSHING, SNAP, TOP OF HEATER          | 1                  |
| 84                         |             |     | BUSHING, SNAP, 1/2" K.O. ELEC.        | 3                  |
| 85                         |             |     | BUSHING, SNAP, BOTTOM OF HEATER       | 1                  |
| 86                         |             |     | LABEL, FRENCH                         | 1                  |
| 87                         |             |     | LABEL, ELEC. DIAGRAM                  | 1                  |
| 88                         |             |     | LABEL, CSA                            | 1                  |
| 89                         |             |     | LABEL, ETL                            | 1                  |
| (Continued on Part 3 of 3) |             |     |                                       |                    |



THERMOCOUPLE RECORDER (OPTIONAL)



REAR VIEW

Figure 7-1. EO Disposer Assembly (Part 3 of 3)

| FIG. & INDEX NO. | PART NUMBER | S V C | DESCRIPTION                                    | UNITS PER ASSEMBLY |
|------------------|-------------|-------|------------------------------------------------|--------------------|
| 7-1-             |             |       | ETO DISPOSER ASSEMBLY (Part 3 of 3) .....      | X                  |
| 5                |             |       | PIPE, 1/2 NPT .....                            | 1                  |
| 6                |             |       | PIPE, 1/4 NPT .....                            | 1                  |
| 7                |             |       | COUPLING, REDUCING 1/2 x 1/4 NPT .....         | 1                  |
| 8                |             |       | BUSHING, REDUCING 1/4 x 1/8 NPT .....          | 1                  |
| 9                |             |       | CONNECTOR, MALE 1/8 NPT x 1/4 ODT .....        | 2                  |
| 10               |             |       | CONNECTOR, FEMALE 1/4 NPT x 1/4 ODT .....      | 1                  |
| 11               |             |       | ELBOW, MALE 1/4 NPT x 1/4 ODT .....            | 1                  |
| 12               |             |       | TUBE, SOFT COPPER 1/4" O.D. (73" REQ'D.) ..... | 1                  |
| 90               | P 764324    | 130   | RECORDER, Thermocouple (Optional) .....        | 1                  |
| 91               | P 764324    | 131   | • PAPER, Box .....                             | A/R                |
| 92               | R 4024      | 123   | WIRE, Thermocouple Lead .....                  | 1                  |






Figure 7-2. Heater Assembly, EO Disposer

7-8  
784323-253

F- 3

| FIG. &<br>INDEX<br>NO. | PART<br>NUMBER |        |     | S<br>V<br>C | DESCRIPTION                  | UNITS PER<br>ASSEMBLY |  |  |  |
|------------------------|----------------|--------|-----|-------------|------------------------------|-----------------------|--|--|--|
| 7-2-                   | P              | 784324 | 739 |             | HEATER ASSEMBLY .....        | X                     |  |  |  |
| 1                      | P              | 784323 | 882 |             | THERMOSTAT, ELECTRONIC ..... | 2                     |  |  |  |
| 2                      | P              | 784323 | 883 |             | HEATER ELEMENT .....         | 11                    |  |  |  |
| 3                      | P              | 784323 | 884 |             | HEATER ELEMENT, RTD .....    | 1                     |  |  |  |
| 4                      | P              | 784323 | 885 |             | RTD, HIGH LIMIT .....        | 1                     |  |  |  |
| 5                      | P              | 784323 | 886 |             | RTD, PROCESS CONTROL .....   | 1                     |  |  |  |
| 6                      | P              | 784323 | 887 |             | TERMINAL BLOCK SECTION ..... | 5                     |  |  |  |
| 7                      | P              | 784323 | 888 |             | TERMINAL END .....           | 2                     |  |  |  |
| 8                      | P              | 784323 | 889 |             | FUSE BLOCK .....             | 1                     |  |  |  |
| 9                      | P              | 784323 | 890 |             | FUSE, 1/4 AMP .....          | 1                     |  |  |  |
| 10                     | P              | 784323 | 891 |             | CONTACTOR .....              | 3                     |  |  |  |
| 11                     | P              | 784323 | 892 |             | CLAMP .....                  | 1                     |  |  |  |

7-9  
784323-253

F- 4






Figure 7-3. Electrical Assembly, EO Disposer

| FIG. & INDEX NO. | PART NUMBER  | SVC | DESCRIPTION                                                                               | UNITS PER ASSEMBLY |
|------------------|--------------|-----|-------------------------------------------------------------------------------------------|--------------------|
| 7-3-             |              |     | ELECTRICAL ASSEMBLY .....                                                                 | X                  |
| 1                | P 764323 009 |     | HOUSING, Pilot Light, Full Voltage, (PL1-4) .....                                         | 4                  |
| 2                | P 764323 733 |     | LAMP, 24 Volt .....                                                                       | 4                  |
| 3                | P 764323 010 |     | LENS, Red, (PL1) .....                                                                    | 1                  |
| 4                | P 764323 734 |     | LENS, Green, (PL2,3) .....                                                                | 2                  |
| 5                | P 764323 735 |     | LENS, Amber, (PL4) .....                                                                  | 1                  |
| 6                | P 764323 736 |     | CIRCUIT BREAKER, Main Disconnect, 40 Amp, 3-Pole, 240 Volt ..                             | 1                  |
| 7                |              |     | • OPERATING MECHANISM .....                                                               | 1                  |
| 8                |              |     | • HANDLE, Flange Operating .....                                                          | 1                  |
| 9                |              |     | • CABLE, 48" Long .....                                                                   | 1                  |
| 10               |              |     | • LINE LUGS .....                                                                         | 3                  |
| 11               | P 764323 737 |     | STARTER CONTACTOR, 2-Pole, 15 Amp, 120 Volt Coll, (M1, M2) .....                          | 2                  |
| 12               | P 764323 007 |     | SOCKET .....                                                                              | 2                  |
| 13               | P 764323 738 |     | TRANSFORMER, 208/380 Volt Primary, 115 VAC Secondary, 500 VA, (T1) .....                  | 1                  |
| 14               | P 764323 002 |     | BLOCK, Fuse, 2-Pole, 250 Volt (FU1,2,3) .....                                             | 3                  |
| 15               | P 764323 014 |     | FUSE, 8 Amp, 250 Volt, (FU1,2) .....                                                      | 4                  |
| 16               | P 764323 015 |     | FUSE, 3 Amp, 250 Volt, (FU3) .....                                                        | 2                  |
| 17               | P 764323 739 |     | BLOCK, Fuse, 3-Pole, 250 Volt (FU3,5) .....                                               | 1                  |
| 18               | P 764323 016 |     | FUSE, 4 Amp, 250 Volt, (FU4) .....                                                        | 1                  |
| 19               | P 764323 017 |     | FUSE, 1 Amp, 250 Volt, (FU7) .....                                                        | 1                  |
| 20               | P 764323 740 |     | FUSE, 1/2 Amp, 250 Volt, (FU5) .....                                                      | 1                  |
| 21               | P 764323 006 |     | TRANSFORMER, 120 Volt Primary, 24 VAC Secondary, 40 VA, (T2) .....                        | 1                  |
| 22               | P 764323 012 |     | RELAY, 4-Pole, Double Throw, 24 VDC, (CR2,3,4) .....                                      | 3                  |
| 23               | P 764323 008 |     | SOCKET .....                                                                              | 3                  |
| 24               | P 764323 741 |     | TIMER, 30 Min. Fixed Delay, 120 VAC Coll, (TD1) .....                                     | 1                  |
| 25               | P 764323 012 |     | RELAY, 24 VAC Coll, 4-Pole, Double Throw, (CR1,CR6,CR7,CR8) .....                         | 4                  |
| 26               | P 764323 008 |     | SOCKET .....                                                                              | 4                  |
| 27               | P 764323 742 |     | LATCH RELAY, Double Pole, Double Throw, Magnetic Latch Type, 2VDC Colls, (CR5) .....      | 1                  |
| 28               | P 764323 007 |     | SOCKET .....                                                                              | 1                  |
| 29               | P 764323 034 |     | TERMINAL BLOCK, 600 V, 30 Amp, Retained, (TB1,TB3) .....                                  | 2                  |
| 30               | P 764323 035 |     | END BARRIER .....                                                                         | 2                  |
| 31               | P 764323 036 |     | END CLAMPS .....                                                                          | 2                  |
| 32               | P 764323 743 |     | POWER SUPPLY, 24 DC, 115 VAC Input, 1.2 Amp, (PU1) .....                                  | 1                  |
| 33               | P 764323 039 |     | DIGITAL TEMP. INDICATOR, 24 VDC Power, 4 Programmable Open Collector Outputs, (TI1) ..... | 1                  |
| 34               | P 764323 744 |     | PRESSURE SWITCH, Differential Type, 0-2" H2O, (PS1,PS2) .....                             | 2                  |
| 35               | P 764323 745 |     | FUSE HOLDER, Panel Mount, (FU6) .....                                                     | 1                  |
| 36               | P 764323 746 |     | FUSE, .75 Amp, 120 VDC .....                                                              | 1                  |
| 37               | P 764323 747 |     | SWITCH, Door Bypass, Push Momentary, Pull Maintained, (LS2) ..                            | 1                  |
| 38               | P 764323 748 |     | SWITCH, Door, Single Pole, Double Throw Momentary, (LS1) .....                            | 1                  |
| 39               | P 764323 749 |     | PUSH BUTTON, Illuminated, 24 VDC, Full Voltage, (ILPB1) .....                             | 1                  |
| 40               | P 764323 872 |     | THERMOCOUPLE, Type "J" Dual, 18" Length, Exposed Junction, (TE1) .....                    | 1                  |
| 41               | P 764323 751 |     | TRANSFORMER, Buck-Boost, 120/240 Primary, 12/24 Secondary, 150 KVA, (T3) .....            | 1                  |



**SIDE VIEW**  
(Vacuum Plumbing)



**FRONT VIEW**  
(Pressure Plumbing)

Figure 7-4. Piping Assembly, EO Disposer

| FIG. & INDEX NO. | PART NUMBER  | SVC | DESCRIPTION                                                                                   | UNITS PER ASSEMBLY |
|------------------|--------------|-----|-----------------------------------------------------------------------------------------------|--------------------|
| 7-4-             |              |     | PIPING ASSEMBLY .....                                                                         | X                  |
| 1                |              |     | BAR, Entrance .....                                                                           | 1                  |
| 2                | P 764323 005 |     | VALVE, Solenoid, 3/4 NPT .....                                                                | 1                  |
| P 764070 002     |              |     | • COIL, 120V, 60 Hz .....                                                                     | 1                  |
| P 764324 031     |              |     | • REPAIR KIT .....                                                                            | 1                  |
| 3                |              |     | NIPPLE, Close, 1/4 NPT .....                                                                  | 1                  |
| 4                |              |     | BUSHING, Reducing, 1/2 x 1/4 NPT .....                                                        | 1                  |
| 5                |              |     | NIPPLE, 1/2 NPT x 1-1/2 LG. ....                                                              | 1                  |
| 6                |              |     | NIPPLE, 1/2 NPT x 2 LG. ....                                                                  | 1                  |
| 7                |              |     | ELBOW, Street, 1/2 NPT x 90° .....                                                            | 2                  |
| 8                |              |     | ELBOW, 1/2 NPT x 90° .....                                                                    | 2                  |
| 9                |              |     | UNION, Standard, 1/2 NPT .....                                                                | 2                  |
| 10               |              |     | NIPPLE, 3/4 NPT x 2 LG. ....                                                                  | 1                  |
| 11               |              |     | BUSHING, Reducing, 3/4 x 1/2 NPT .....                                                        | 2                  |
| 12               |              |     | TUBE, Copper, 3/8 ODT x 2-5/8 LG. ....                                                        | 2                  |
| 13               |              |     | TUBE, Copper, 3/8 ODT x 9-3/8 LG. ....                                                        | 1                  |
| 14               |              |     | ELBOW, 3/8 ODT x 90° .....                                                                    | 1                  |
| 15               |              |     | TEE, 1/2 NPT x 1/2 ODT .....                                                                  | 1                  |
| 16               |              |     | ELBOW, 1/4 NPT x 1/2 ODT .....                                                                | 2                  |
| 17               | P 764324 030 |     | HOSE, 1/2" Corrugated, W/1/2 NPT Hex Male Fitting Welded Both Ends x 11-1/2" O/A Length ..... | 2                  |

THIS PAGE  
INTENTIONALLY LEFT  
BLANK

## REFERENCE DRAWINGS

### DISPOSER

|                                       |        |          |
|---------------------------------------|--------|----------|
| Electrical Schematic Disposer 100 CFM | 1 of 3 | P19-3610 |
|                                       | 2 of 3 |          |
|                                       | 3 of 3 |          |

|                                        |        |           |
|----------------------------------------|--------|-----------|
| Wiring Connections Sterilizer/Disposer | 1 of 2 | 93910-145 |
|                                        | 2 of 2 |           |

### MEDALLION STERILIZERS

|                                                          |        |  |
|----------------------------------------------------------|--------|--|
| Small Medallion Cryotherm Piping Schematic with Disposer | 1 of 1 |  |
|----------------------------------------------------------|--------|--|

|                                                           |        |  |
|-----------------------------------------------------------|--------|--|
| Medium Medallion Cryotherm Piping Schematic with Disposer | 1 of 1 |  |
|-----------------------------------------------------------|--------|--|

|                                   |        |            |
|-----------------------------------|--------|------------|
| ETO Disposer Kit, Small Medallion | 1 of 3 | 146655-109 |
|                                   | 2 of 3 |            |
|                                   | 3 of 3 |            |

|                                             |        |            |
|---------------------------------------------|--------|------------|
| Valve Assembly - Small Medallion w/Disposer | 1 of 1 | 136807-659 |
|---------------------------------------------|--------|------------|

|                                    |        |            |
|------------------------------------|--------|------------|
| EO Disposer Kit - Medium Medallion | 1 of 3 | 146655-115 |
|                                    | 2 of 3 |            |
|                                    | 3 of 3 |            |

|                                              |        |            |
|----------------------------------------------|--------|------------|
| Valve Assembly - Medium Medallion w/Disposer | 1 of 1 | 136807-660 |
|----------------------------------------------|--------|------------|

|                                                |        |            |
|------------------------------------------------|--------|------------|
| Elementary Diagram - All Cryotherm Sterilizers | 1 of 1 | 146655-103 |
|------------------------------------------------|--------|------------|

|                                            |        |            |
|--------------------------------------------|--------|------------|
| Wiring Diagram - All Cryotherm Sterilizers | 1 of 1 | 146655-104 |
|--------------------------------------------|--------|------------|

|                                              |        |            |
|----------------------------------------------|--------|------------|
| Elementary Diagram Medallion Gas Sterilizers | 1 of 1 | 146655-116 |
|----------------------------------------------|--------|------------|

|                                                                                 |        |            |
|---------------------------------------------------------------------------------|--------|------------|
| Wiring Diagram - All 24 x 36 Medallion Sterilizers with Disposer and Envirogard | 1 of 1 | 146655-117 |
|---------------------------------------------------------------------------------|--------|------------|

|                                                             |        |            |
|-------------------------------------------------------------|--------|------------|
| Electrical Schematic Medallion-Disposer/Sequencer Interface | 1 of 2 | 146657-254 |
|                                                             | 2 of 2 |            |

|                                       |        |            |
|---------------------------------------|--------|------------|
| Interface Box Assembly (Before 11/91) | 1 of 2 | 146655-106 |
|                                       | 2 of 2 |            |

|                                      |        |            |
|--------------------------------------|--------|------------|
| Interface Box Assembly (After 11/91) | 1 of 3 | 146655-718 |
|                                      | 2 of 3 |            |
|                                      | 3 of 3 |            |

# **EAGLE 2000 STERILIZERS**

|                                                        |                                      |            |
|--------------------------------------------------------|--------------------------------------|------------|
| Piping Schematic Eagle 20 x 20 w/IRC                   | 1 of 1                               |            |
| Piping Schematic Eagle 20 x 20 w/Printcon              | 1 of 1                               |            |
| Piping Schematic Eagle 24 x 36 w/IRC                   | 1 of 1                               |            |
| Piping Schematic Eagle 24 x 36 w/Printcon              | 1 of 1                               |            |
| ETO Disposer Kit Assembly 20 x 20 Eagle IRC 2000       | 1 of 3<br>2 of 3<br>3 of 3           | 146655-119 |
| ETO Disposer Kit Assembly 20" Printcon                 | 1 of 3<br>2 of 3<br>3 of 3           | 146653-868 |
| Valve Assembly - 20" Gas Print. or IRC with Disposer   | 1 of 1                               | 136806-983 |
| ETO Disposer Kit Assembly 24 x 36 Eagle IRC Eagle 2000 | 1 of 4<br>2 of 4<br>3 of 4<br>4 of 4 | 146655-120 |
| ETO Interface Kit Assembly 24 x 36 Printcon            | 1 of 4<br>2 of 4<br>3 of 4<br>4 of 4 | 146655-118 |
| Valve Assembly 24 x 36 Eagle 2000 w/Disposer           | 1 of 1                               | 136807-661 |
| Electrical Schematic Eagle IRC w/Disposer              | 1 of 3<br>2 of 3<br>3 of 3           | 146655-121 |
| Electrical Schematic                                   | 1 of 3<br>2 of 3<br>3 of 3           | 146653-873 |
| IRC - Disposer Sequencer Interface                     | 1 of 1                               | 146657-256 |
| Printcon - Disposer/Sequencer Interface                | 1 of 1                               | 146657-255 |
| Interface Box Assembly Prnt. & Rnd Chrt w/Disposer     | 1 of 1                               | 136807-997 |

# **EAGLE 3000 STERILIZERS**

|                                                  |                            |            |
|--------------------------------------------------|----------------------------|------------|
| Piping Schematic 20 x 20 Gas-Eagle 3000          | 1 of 1                     | 146655-712 |
| Piping Schematic Medium Gas-Eagle 3000           | 1 of 1                     | 146655-713 |
| ETO Disposer Kit Assembly 20" Eagle 3000 Gas     | 1 of 3<br>2 of 3<br>3 of 3 | 146655-710 |
| Valve Assembly - 20" Eagle 3000 w/Disposer       | 1 of 1                     | 136807-904 |
| ETO Disposer Kit Assembly 24 x 36 Eagle 3000     | 1 of 3<br>2 of 3<br>3 of 3 | 146655-711 |
| Valve Assembly 24 x 36 Eagle 3000 w/Disposer     | 1 of 1                     | 136807-905 |
| EO Sterilizer/ Aerator Machine Schematic         | 1 of 2<br>2 of 2           | 146655-583 |
| Eagle 3000 Disposer/Sequencer Interface          | 1 of 1                     | 146657-275 |
| Interface Box Assembly Eagle 3000/Gas w/Disposer | 1 of 1                     | 136807-908 |

# **AERATORS & MISCELLANEOUS**

|                                               |                            |                |
|-----------------------------------------------|----------------------------|----------------|
| Wiring Connections Aerator w/Disposer         | 1 of 3<br>2 of 3<br>3 of 3 | 136807-670     |
| Schematic Modified Aerator w/Disposer         | 1 of 1                     | 136807-671     |
| Electrical Schematic Eagle Disposer/Sequencer | 1 of 3<br>2 of 3<br>3 of 3 | 146657-295     |
| Wiring Connections Sequencer Control          | 1 of 2<br>2 of 2           | 136807-669     |
| Interface Box Assembly                        | 1 of 2<br>2 of 2           | 146655-107     |
| PC Board Assembly, Disposer Interface         | 1 of 3<br>2 of 3<br>3 of 3 | 146655-715     |
| Control Board, Eaglet II                      | 1 of 2<br>2 of 2           | 146653-874/875 |

# **Service Kits**

Service Repair Kit for Pump

1 of 1

764323-249

Disposer Water Trap Kit

1 of 2

755716-132

2 of 2


VS23 Installation Kit

1 of 2


755716-164

2 of 2


THIS PAGE  
INTENTIONALLY LEFT  
BLANK



ABATOR  
HOUSING ASSEMBLY  
SCALE: 1/8



LEFT SIDE SUBPANEL  
SCALE: 1/4



RIGHT SIDE SUBPANEL  
SCALE: 1/4

# BILL OF MATERIAL

| ITEM NO. | QTY. | DESCRIPTION                                                                        | CODE       |
|----------|------|------------------------------------------------------------------------------------|------------|
| 1        | 4    | Pilot Light, Full Voltage Housing                                                  | PL-1-4     |
| 2        | 1    | Diablo 481-5710-01-303                                                             | PL-2       |
| 3        | 1    | Lamp, 24 Volt - Sylvania 81819                                                     | PL-1,3     |
| 4        | 1    | Lens, Red - Diablo 826-1191-300                                                    | PL-4       |
| 5        | 1    | Lens, Green - Diablo 826-1192-300                                                  |            |
| 6        | 1    | Lens, Amber - Diablo 826-1193-300                                                  |            |
| 7        | 1    | Main Disconnect Circuit Breaker, 40 Amp., 3-Pole, 240 Volt                         | MAIN DISC  |
| 8        | 1    | G.E. 8TCB132049H                                                                   |            |
| 9        | 1    | Operating Mechanism - Siemens 8FNDCE80                                             |            |
| 10       | 1    | Flange Operating Handle 8FNDH                                                      |            |
| 11       | 1    | Cable, 48" Long - Siemens 8FNDCE-048                                               |            |
| 12       | 1    | Line Lugs                                                                          |            |
| 13       | 1    | LTE Electrical - BLU-1100                                                          |            |
| 14       | 2    | Starter Contactor, 2-Pole, 15 Amp. 120 Volt Coil                                   | H-1, H-2   |
| 15       | 1    | Porter & Brunfield 8KLP-11A18-120                                                  |            |
| 16       | 2    | Socket - 1/2" 8KSL711                                                              |            |
| 17       | 1    | Transformer, 200/380 Volt Primary, 115 VAC Secondary, 500 VA                       | T-1        |
| 18       | 1    | G.E. 89T58899                                                                      |            |
| 19       | 3    | Fuse Block, 2-Pole, 250 Volt                                                       | FU-1,2,3   |
| 20       | 4    | Marathon 83FA2SP                                                                   |            |
| 21       | 3    | Fuse, 8 Amp., 250 Volt - Buss 8FBN-R8                                              | FU-1,2     |
| 22       | 2    | Fuse, 3 Amp., 250 Volt - Buss 8FBN-R3                                              | FU-3       |
| 23       | 1    | Fuse Block, 3-Pole, 250 Volt                                                       | FU-3,5     |
| 24       | 1    | Marathon 83FA2SP                                                                   |            |
| 25       | 1    | Fuse, 4 Amp., 250 Volt - Buss 8FBN-R4                                              | FU-4       |
| 26       | 1    | Fuse, 1 Amp., 250 Volt - Buss 8FBN-R1                                              | FU-7       |
| 27       | 1    | Fuse, 1/2 Amp., 250 Volt - Buss 8FBN-R1/2                                          | FU-5       |
| 28       | 1    | Transformer, 120 Volt Primary, 24 VAC Secondary, 48 VA                             | T-2        |
| 29       | 1    | Honeywell 8ATT20-1009                                                              |            |
| 30       | 3    | Relay, 4-Pole, Double Throw, 24 VDC                                                | CR-2,3,4   |
| 31       | 2    | Porter & Brunfield 8KMU17B11-24                                                    |            |
| 32       | 1    | Socket - TE 8KSL715                                                                |            |
| 33       | 1    | Timer, 30 Minute Fixed Delay, 120 VAC Coil                                         | TB-1       |
| 34       | 1    | Syracuse Electronics 8SPABA13-30H                                                  |            |
| 35       | 4    | Relay, 24 VAC Coil, 4-Pole Double Throw                                            | CR1, CR7   |
| 36       | 4    | Porter & Brunfield 8KMU17A11-24 VAC                                                | CR6, CR8   |
| 37       | 1    | Socket - TE 8KSL715                                                                |            |
| 38       | 1    | Latch Relay, Double Pole, Double Throw, Magnetic Latch Type, 24VDC Coils           | CR5        |
| 39       | 1    | Porter & Brunfield 8KMU17B11-24                                                    |            |
| 40       | 1    | Socket - TE 8KSL711                                                                |            |
| 41       | 45   | Terminal Block, 600 Volt, 30 Amp. Retained                                         | TB-1, TB-3 |
| 42       | 1    | Square D 89880-GR6                                                                 |            |
| 43       | 2    | End Barrier - Square D 89880-GR6                                                   |            |
| 44       | 1    | End Clamps - Square D 89880-GR10                                                   |            |
| 45       | 1    | Power Supply, 24 DC, 115 VAC Input                                                 | PU-1       |
| 46       | 1    | Condor 8H24-1.2                                                                    |            |
| 47       | 1    | Digital Temperature Indicator, 24 VDC Power, 4 Programmable Open Collector Outputs | TI-1       |
| 48       | 1    | Analogic 81TG-2600                                                                 |            |
| 49       | 1    | Pressure Switch, Differential Type                                                 | PS-1       |
| 50       | 1    | 8" 100                                                                             |            |
| 51       | 1    | Dwyer 81823-2                                                                      |            |
| 52       | 1    | Fuse Holder, Panel Mount                                                           | FU-4       |
| 53       | 1    | Buss 8HOP                                                                          |            |
| 54       | 1    | Fuse, .75 Amp., 120 VDC                                                            |            |
| 55       | 1    | Buss 8HOB-3/4                                                                      |            |
| 56       | 1    | Door Bypass Switch, Push Momentary, Pull Maintained                                | LS-2       |
| 57       | 1    | Micro Switch 81BH-430                                                              |            |
| 58       | 1    | Door Switch, Single Pole, Double Throw Momentary                                   | LS-1       |
| 59       | 1    | Micro Switch 81BH1                                                                 |            |
| 60       | 1    | Illuminated Pushbutton, 24 VDC, Full Voltage                                       | ILPB-1     |
| 61       | 1    | G.E. 8KCB18A11A18                                                                  |            |
| 62       | 1    | Type "J" Terminal Board                                                            | TB-2       |
| 63       | 1    | Pyramation 824-228-82                                                              |            |
| 64       | 1    | Thermocouple Type "J" Buss, 18" Length                                             | TE-1       |
| 65       | 1    | Capacitor Junction Pyramation 8J2918181884                                         |            |
| 66       | 1    | Transformer, Buck - Boost, 120/240 Primary 12/24 secondary, .150 KVA               | T3         |
| 67       | 1    | G.E. 8 9T518185                                                                    |            |


| REV | ZONE | REVISION DESCRIPTION    | APPROVAL/DATE |
|-----|------|-------------------------|---------------|
| 1   |      | SHOP AS BUILT           | 7/24/89       |
| 2   |      | ADDED BOOST TRANSFORMER | 7/31/89       |
| 3   |      | RELEASED TO PRODUCTION  | 9/18/89       |

|                  |  |                     |  |                                          |  |
|------------------|--|---------------------|--|------------------------------------------|--|
| MATERIAL         |  | RELEASE NO. 89A5742 |  | TITLE 100 CFM ABATOR ELECTRICAL ASSEMBLY |  |
| DATE 9-18-89     |  | DATE 9-18-89        |  | NUMBER                                   |  |
| REF DATE 6-20-89 |  | REF DATE 6-20-89    |  | P19-3610                                 |  |
| DATE 9/20/89     |  | DATE 9/20/89        |  | SCALE AS NOTED SHEET 1 OF 3              |  |
| DATE 9/20/89     |  | DATE 9/20/89        |  | REVISION 3                               |  |
| DATE 9/20/89     |  | DATE 9/20/89        |  | SECTIONAL SIZE OUTLINE SIZE              |  |
| DATE 9/20/89     |  | DATE 9/20/89        |  | SECTION CONTROL DWG LOCATION             |  |
| DATE 9/20/89     |  | DATE 9/20/89        |  | DWG LOCATION                             |  |
| DATE 9/20/89     |  | DATE 9/20/89        |  | DWG LOCATION                             |  |

G- 1

G- 2

| CADD | REV | DESCRIPTION                   | DATE        |
|------|-----|-------------------------------|-------------|
| CADD | 1   | SHOP AS BUILT                 | 7/24/88     |
| CADD | 2   | ADDED BOOSTER TRANSFORMER     | 7/24/88     |
| CADD | 3   | ADDED SHIELDING AND GROUNDING | 8/2/88      |
| CADD | 4   | RELEASED TO PRODUCTION        | NOV 9, 1988 |





| REV | ZONE | REVISION DESCRIPTION    | APPROVAL/DATE |
|-----|------|-------------------------|---------------|
| 1   |      | SHOP AS BUILT           | 7/25/89       |
| 2   |      | ADDED BOOST TRANSFORMER | 7/31/89       |
| 3   |      | RELEASED TO PRODUCTION  | 9/16/89       |

# MFC ANALOGIC

PART # ITC2600-OX-0  
SUPPLIER: PROCESS MEASUREMENT CO.  
MPLS. MINN.

INPUT THERMOCOUPLE TYPE J OR K (CIRCLE 1)

CONTROL SECTION SETUP PARAMETERS:

| OUTPUT # | ACTION | SETPOINT | DEADBAND |
|----------|--------|----------|----------|
| 1        | HIGH   | 280      | 10       |
| 2        | HIGH   | 330      | 10       |
| 3        | LOW    | 440      | 10       |
| 4        | HIGH   | 500      | 10       |

NOTE: IF HIGH ACTION IS SET, THEN THE RELAYS ENERGIZE (CLOSE) WHEN THE INPUT INCREASES ABOVE THE SETPOINT; THEY WILL DE-ENERGIZE (OPEN) BELOW THE SETPOINT MINUS THE DEADBAND VALUE.

IF LOW ACTION IS SET, THEN THE RELAYS ENERGIZE (CLOSE) WHEN THE INPUT DECLINES BELOW THE SETPOINT; THEY WILL DE-ENERGIZE ABOVE THE SETPOINT PLUS THE DEADBAND VALUE.

|                        |  |                                       |  |                                                    |  |                                               |  |
|------------------------|--|---------------------------------------|--|----------------------------------------------------|--|-----------------------------------------------|--|
| MATERIAL               |  | RELEASE NO. 89A5742                   |  | TITLE 100 CFM ABATOR PROGRAM SETPOINT FOR ITC 2600 |  | 07/31/89 12:43 R3401A03                       |  |
| REFERENCE              |  | DATE 9-18-89                          |  | EFFECT DATE 6-20-89                                |  | NUMBER P19-3610                               |  |
| NEXT ASSY. 61X00-6500  |  | DATE 8/20/89                          |  | CNC                                                |  | SCALE 1/4" = 1"                               |  |
| PREVIOUS DRAWING NO.   |  | APPROVED 9/2/89                       |  | DATE 9-21-89                                       |  | SHEET 3 OF 3                                  |  |
| DO NOT SCALE           |  | TOLERANCES UNLESS OTHERWISE SPECIFIED |  | 3 PL. 3                                            |  | SECTIONAL SIZE OUTLINE SIZE                   |  |
| THIRD ANGLE PROJECTION |  | DESIGN CONTROL                        |  | DWG LOCATION                                       |  | VERSION 3                                     |  |
|                        |  |                                       |  |                                                    |  | Donaldson Industrial Group<br>Minneapolis, MN |  |


G- 5

G- 6

|      |  |                    |  |              |  |          |  |                       |  |            |  |                        |  |                       |  |                        |  |
|------|--|--------------------|--|--------------|--|----------|--|-----------------------|--|------------|--|------------------------|--|-----------------------|--|------------------------|--|
| REF. |  | DWG. NO. 93910-145 |  | SHEET 2 OF 2 |  | QUANTITY |  | PART NUMBER 93910-145 |  | ITEM NO. 1 |  | PART NAME WIRING CONN. |  | DESCRIPTION, MATERIAL |  | CAD DWG.               |  |
|      |  |                    |  |              |  |          |  |                       |  |            |  |                        |  | REV NO. 7             |  | REVISION DATE 11-01-91 |  |


  

INTERFACE BOX  
SEE SHT. 2 OF 2




STERILIZER  
SIDE VIEW  
20X20 SHOWN

NOTE: MAXIMUM CABLE  
LENGTH 200'



DISPOSER  
SIDE VIEW

TERMINAL STRIP  
SEE VIEW B-B



DISPOSER  
FRONT VIEW  
DOOR OPEN

| FUNCTION       | WIRE # | TERM STRIP |
|----------------|--------|------------|
| POWER          | 27     | TB#27      |
| POWER RETURN   | 26     | TB#26      |
| VACUUM ON      | 38     | TB#38      |
| READY          | 32     | TB#32      |
| DISPOSER ON    | 37     | TB#37      |
| WORKING        | 31     | TB#31      |
| AIRFLOW        | 34     | TB#34      |
| OVERTEMP       | 28     | TB#28      |
| OVERTEMP RESET | 10     | TB#10      |
| OVERTEMP RESET | 12     | TB#12      |

CABLE SPEC:  
10 CONDUCTOR JACKETED CABLE, 22  
OR 24 AWG. STRANDED, BELDEN #8456  
OR EQUAL. CHECK STATE AND/OR LOCAL  
CODES.

7 11-01-91 910211

|                        |  |                                                      |  |
|------------------------|--|------------------------------------------------------|--|
| TOLERANCE STANDARD     |  | FRACTIONAL $\pm 1/64$ DECIMAL $\pm .005$             |  |
| UNLESS OTHERWISE NOTED |  | ANGULARITY $\pm 1^\circ$ MACH. SURF. $\pm .25^\circ$ |  |

|   |          |        |   |          |        |
|---|----------|--------|---|----------|--------|
| 3 | 11-28-89 | RA #82 | 6 | 07-17-90 | 900200 |
| 2 | 11-06-89 | RA #71 | 6 | 07-17-90 | RA#145 |
| 1 | 09-06-89 | RA #54 | 5 | 06-13-90 | RA#133 |
| 0 | 6-20-89  | PER 25 | 4 | 05-08-90 | RA 122 |

SCALE  $\frac{1}{16}$   
1



AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA

THIS DRAWING, SPECIFICATIONS, AND ALL INFORMATION THEREIN ARE  
THE PROPERTY OF AMERICAN STERILIZER COMPANY AND ARE TO BE  
MAINTAINED CONFIDENTIAL. THEY ARE NOT TO BE COPIED OR USED IN ANY  
MANNER DETRIMENTAL TO THE INTEREST OF AMERICAN STERILIZER  
COMPANY AND SHALL BE RETURNED UPON REQUEST.

AMSCO

TITLE WIRING, CONNECTIONS  
STERILIZER/DISPOSER

FIRST MADE FOR:


|             |                  |              |
|-------------|------------------|--------------|
| DWN JMT     | CKD JMT-GENG LML | MFG          |
| DATE 1-9-89 | DATE 6-20-89     | DATE 6-20-89 |

|                    |  |              |  |
|--------------------|--|--------------|--|
| DWG. NO. 93910-145 |  | SHEET 2 OF 2 |  |
|--------------------|--|--------------|--|



|                |              |          |             |          |              |                       |         |               |
|----------------|--------------|----------|-------------|----------|--------------|-----------------------|---------|---------------|
| REF. 93910-145 | 2 SHEET OF 2 | QUANTITY | PART NUMBER | ITEM NO. | PART NAME    | DESCRIPTION, MATERIAL | REV NO. | REVISION DATE |
|                |              |          | 93910-145   | 1        | WIRING CONN. |                       | 7       | 11-01-91      |

REF.-EAGLE 3000 &  
EAGLE 2000  
INTERFACE BOX




REF.-EAGLE 3000 &  
EAGLE 2000

FROM STERILIZER / 10



INTERFACE  
CONTROL  
BOX



VIEW A-A

**VIEW B-B**  
DISPOSER  
TERMINAL STRIP

MEDALLION STERILIZER

TOLERANCE STANDARD FRACTIONAL  $\pm 1/64$  DECIMAL  $\pm .005$   
UNLESS OTHERWISE NOTED ANGULARITY  $\pm 1^\circ$  MACH. SURF.  $\pm .025$

| NO. | DATE | E.C.A. NUMBER | NO. | DATE | E.C.A. NUMBER |
|-----|------|---------------|-----|------|---------------|
|     |      |               |     |      |               |
|     |      |               |     |      |               |
|     |      |               |     |      |               |
|     |      |               |     |      |               |
|     |      |               |     |      |               |



AMSCO

AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA

THIS DRAWING, SPECIFICATIONS AND ALL INFORMATION THEREIN ARE  
THE PROPERTY OF AMERICAN STERILIZER COMPANY AND ARE TO BE  
MAINTAINED CONFIDENTIAL. ARE NOT TO BE COPIED OR USED IN ANY  
WAY DETRIMENTAL TO THE INTEREST OF AMERICAN STERILIZER  
COMPANY AND SHALL BE RETURNED UPON REQUEST.


DWN JMT CKD JMT-GENG LML MFG  
DATE 1-5-89 DATE 6-20-89 DATE 6-20-89 DATE

TITLE  
WIRING, CONNECTIONS  
STERILIZER/DISPOSER


FIRST MADE FOR:

DWG. NO. 93910-145

2 SHEET OF 2



**Small Medallion Cryotherm  
Piping Schematic  
With Disposer**




# RAMSCO

## EO DISPOSER SYSTEM

### INTERFACE EAGLE 3000 & 7000 MEDIA AERATION CABINET

2000-0000-0000






| QUANTITY | PART NUMBER | PART NAME                 | DESCRIPTION, MATERIAL           | REVISION   |
|----------|-------------|---------------------------|---------------------------------|------------|
| 1        | 146655-109  | DISPOSER KIT              | 20" MEDALLION                   | 1903-06-92 |
| 1        | 36807-659   | 2 VALVE ASS'Y             | S41, S28                        |            |
| 1        | 46655-103   | 3 SCHEMATIC, ELEC         | REF. DWG.                       |            |
| 1        | 46655-718   | 4 BOX, INTERFACE          | REF. DWG.                       |            |
| 1        | 29361-106   | 5 INSTALL. INST.          | REF. DWG.                       |            |
| 1        | 83866-001   | 6 SWITCH, VACUUM          |                                 |            |
| 1        | 93910-136   | 7 LIGHT BOX ASS'Y         |                                 |            |
| 1        | 93910-138   | 8 BRKT. PIPE 'SPT.        |                                 |            |
| 2        | 18398-045   | 9 U BOLT                  | 1/4" (W/NUTS & CLAMP)           |            |
| 1        | 36807-915   | 10 CABLE                  |                                 |            |
| 1        | 29361-132   | 11 BOX, MODIFIED          |                                 |            |
| 1        | 29361-131   | 12 COVER, BOX             |                                 |            |
| 1        | 29361-116   | 13 FLAT ELL               |                                 |            |
| 2        | 29361-115   | 14 WALL TRACK             | 5' LG.                          |            |
| 2        | 29361-117   | 15 SPLICE WALL TRACK      |                                 |            |
| 6        | 9316-041    | 16 SCREW, RD. HD.         | #10-32 X 3/4" LG.               |            |
| 6        | 20844-061   | 17 WASHER, FLAT           | 203 ID X .438 OD X .031 THK     |            |
| 6        | 129359-532  | 18 NUT, KEPS              | #10-32                          |            |
| 14       | 18538-091   | 19 NUT, WIRE              |                                 |            |
| 3        | 29326-001   | 20 MOUNT, HARNESS         |                                 |            |
| 6        | 84104-001   | 21 TIE, WIRE              |                                 |            |
| 3        | 34218-091   | 22 FITTING, STR.          | 1/4 OD X 1/4" NPT               |            |
| 1        | 4906-091    | 23 TEE                    | 1/4" OD X 1/4" LG. SEE NOTE #2  |            |
| 2        | 28917-091   | 24 NIPPLE                 | 1/4" NPT X 1" LG.               |            |
| 1        | 129359-533  | 25 NUT, KEPS              | 1/4-20                          |            |
| 1        | 150473-295  | 26 WASHER, FLAT           | .281 ID X .625 OD X .063 THK    |            |
| 1        |             | 27 PRINT                  | REF. P/N 146655-109 SEE NOTE #3 |            |
| 1        | 146655-104  | 28 WIRING DIA.            |                                 |            |
| 6        | 91354-045   | 29 SQUEEZE CLAMP          |                                 |            |
| 4        | 3952-041    | 30 SCREW, RD. HD.         | 1/4-20 X 1 1/2" LG.             |            |
| 1        | 150828-475  | 31 VALVE, SAFETY          | 1/2 NPT. SEE NOTE #5            |            |
| 1        | 93910-118   | 32 LABEL, PHASE INDICATOR |                                 |            |
| 1        | 29013-091   | 33 NIPPLE, BRASS          | 3/8 NPT X 1" LG.                |            |
| 1        | 84443-001   | 34 BUSHING, REDUCING      | 3/4 NPT X 3/8 NPT               |            |
| 1        | 129361-125  | 35 EQUIPMENT, MANUAL      | REF. DWG.                       |            |
| 2        | 12529-061   | 36 SCREW, RD. HD.         | #6-32 X 1/4" LG.                |            |
| 1        | 49087-044   | 37 PLUG, BUTTON           |                                 |            |
| 1        | 93910-145   | 38 WIRING CONNECTIONS     | DISPOSER - REF. DWG.            |            |
| 1        | 44500-091   | 39 ELL, 90°               | SOLDER JOINT FITTING            |            |
| 1        | 836-091     | 40 BUSHING, REDUCING      | 3/4 NPT X 1/2 NPT               |            |
| 1        | 150822-588  | 41 LENS, AERATING         |                                 |            |
| 3        | 150778-001  | 42 VARIATOR               |                                 |            |
| 1        | 150777-001  | 43 CAPACITOR              | 0.1 MFD. 500V                   |            |
| 1        | 45589-091   | 44 COUPLING               | 1/4 NPT                         |            |
| 1        | 150828-262  | 45 VALVE, CHECK           |                                 |            |
| 1        | 90331-091   | 46 CONN. STR.             | 3/8                             |            |
| 1        | 24748-091   | 47 BUSHING                | 3/8                             |            |
| 25       | 129361-897  | 48 CONTACT PTN            | 24-20 AWG                       |            |
| 1        | 939-091     | 49 SOLDER BUSHING         | 1/8 NPT X 1/4 NPT               |            |
| 1        | 76053-091   | 50 RED BUSHING            | 1/2 NPT X 1/4 NPT               |            |
| 1        | 90525-091   | 51 GAUGE                  | 1/8 NPT                         |            |
| 1        | 129356-508  | 52 TEE SOLDER             | 7/8 ODT X 7/8 ODT X 1/2 NPT     |            |
| 1        | 93910-642   | 53 VALVE, BALL            | 7/8 ODT                         |            |
| 1        | 764317-524  | 54 TAPE, TEFLON           |                                 |            |
| 1        | 42510-091   | 55 FTG ELL, COMP          | 1/4 ODT X 1/8 NPT               |            |
| 1        | 93910-653   | 56 TRAP, WATER            | 3/4 NPT                         |            |
| 1        | 22472-091   | 57 BUSHING, RED.          | 3/4 NPT X 1/4 NPT               |            |
| 1        | 45407-091   | 58 ELL, COMP.             | 1/4 ODT X 1/4 NPT               |            |
| 1        | 129362-980  | 59 GUIDELINES, DISPOSER   |                                 |            |
| 1        | 146657-254  | 60 SCHEMATIC              |                                 |            |
| 1        | 150822-048  | 61 VALVE, BALL            |                                 |            |

# NOTES:

- USE PIPE SEALANT PER AMSCO DWG. 150822-020 ON ALL THREADED PIPE JOINTS.
- COPPER TUBING .030 WALL X 1/4" OD, SOFT ANNEALED, ASTM B280.
- INCLUDE WITH EACH KIT ONE PRINT, P/N 146655-109.
- MOUNT LIGHT BOX ASS'Y TO BEST SUIT. EITHER USE WALL TRACK COMPONENTS (ITEM #13, #14, #15) OR USE CLAMP (ITEM #30) FOR ENTRY INTO BACK OF LIGHT BOX ASS'Y.
- REPLACE EXISTING SAFETY VALVE WITH ITEM #32.




AMERICAN STERILIZER COMPANY  
ETIO DISPOSER KIT  
SMALL MEDALLION  
FIRST MADE FOR  
146655-109 11 3




A- 3


A- 4


















B- 1


B- 2
















CONTINUED ON SHEET 2 OF 2


CAD DWG.

|                                              |      |     |     |                                                                |      |     |     |
|----------------------------------------------|------|-----|-----|----------------------------------------------------------------|------|-----|-----|
| TOLERANCE STANDARD<br>UNLESS OTHERWISE NOTED |      |     |     | FRACTIONAL 11/64 DECIMAL : 005<br>ANGULARITY: 1° MACH SURF 125 |      |     |     |
|                                              |      |     |     |                                                                |      |     |     |
|                                              |      |     |     |                                                                |      |     |     |
|                                              |      |     |     |                                                                |      |     |     |
|                                              |      |     |     |                                                                |      |     |     |
|                                              |      |     |     |                                                                |      |     |     |
| 0 03-26-92 320052                            |      |     |     |                                                                |      |     |     |
| BY                                           | DATE | REV | QTY | BY                                                             | DATE | REV | QTY |





|                                                                                       |                                                                                                                                                                                                                                                                                      |                                            |
|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------|
|  | AMERICAN STERILIZER COMPANY<br>EPIC, PENNSYLVANIA<br>We design and manufacture all types of sterilizers for<br>the medical, dental, pharmaceutical, and food processing<br>industries. We have a complete line of sterilizers for<br>your selection. Write for our literature today. | TITLE<br>INTERFACE BOX ASS'Y.<br>MEDALLION |
|                                                                                       | FIRST MADE FOR: 146655-109, 146655-115<br>Q. 146655-106                                                                                                                                                                                                                              | 1 of 2                                     |



146655-106  
146655-718 146655-718

NOTES:

1. WIRE #22 GA STRANDED COPPER, 19/34 STRANDING, TEFLON INSULATION, U.L. (STYLE 1213, 300 VOLT, 105°C) AND CSA (APPLIANCE WIRE, .010 THK INSULATION, 150 VOLT, 150°C) LABELS REQ'D.
2. WIRE #22 GA STRANDED COPPER, 7/30 STRANDING, 300V, 80°C, PVC INSULATION.
3. WIRE 18 GA STRANDED COPPER, 16/30 STRANDING, 600V, 105°C, THERMOPLASTIC INSULATION.
4. CONNECTION METHOD "A":  ITEM #41, STRIP 3/16"  
CONNECTION METHOD "B":  ITEM #42, STRIP 1/4"  
CONNECTION METHOD "C":  ITEM #43, STRIP 1/4"  
CONNECTION METHOD "D": SOLDER - STRIP 1/4" AND TIN. SEE NOTE #7.  
CONNECTION METHOD "E": STRIP 5/16" AND TIN  
CONNECTION METHOD "F": NO STRIP  
CONNECTION METHOD "G":  ITEM #44, STRIP 1/4"
5. CAREFULLY APPLY COVER LABEL, ITEM #25, BEFORE WIRING SWITCH, ITEM #15, AND L. LIGHTS, ITEM #14.
6. BLACK INSULATION SLEEVING, BRAIDED GLASS FIBER/VINYL COVERED NEMA VS-1, GRADE B, SIZE 14 (.064/.072 ID).
7. ALL SOLDER MATERIAL SHALL BE IN ACCORDANCE WITH AMSCO DWG. 52595.
8. LABEL COMPONENTS AS SHOWN WITH BRADY TYPE MARKERS, BLACK LETTERS ON YELLOW BACKGROUND, OR EQUIVALENT, LABELS REQ'D: TB2, TB3, TB4, TB5, OVERTEMP RESET, JS4.
9. LABEL WIRES 4" LONG OR LESS ONCE IN CENTER, LABEL LONGER WIRES 1" FROM EACH END USING BRADY TYPE MARKERS, BLACK LETTERS ON WHITE BACKGROUND.
10. PLACE TIE WRAPS, ITEM #47, WHERE NECESSARY TO MAINTAIN A NEAT ASS'Y.
11. DO NOT TIE WRAP WIRES #4 & #5 TO ANY OTHER WIRING. USE HARNESS MOUNTS, ITEM #46, TO SECURE TO BOX.
12. PULL EXCESS WIRE LENGTH, IF ANY, OF CABLE ASSEMBLIES, ITEMS 23 & 24, OUT OF BOX PRIOR TO TIE-WRAPPING.
13. DRESS WIRES AWAY FROM ALL FUSES. FUSE ACCESS MUST NOT BE OBSTRUCTED.
14. ITEMS 72, 73 REFERRED TO AS "GROUND CABLE" IN WIRING CHART.
15. SHUNTS FOR JUMPERS JP1 & JP2 ON INTERFACE BOARD, ITEM 7, MAY BE IN EITHER POSITION.
16. SHUNT FOR JUMPER JP1 ON EAGLET BOARD, ITEM 9, MUST BE IN "WDOG" POSITION.
17. CHART FOR JUMPER BAR ITEM #18
  - TB5-1 TO TB5-2
  - TB5-2 TO TB5-3
  - TB5-3 TO TB5-4
  - TB5-4 TO TB5-5
  - TB5-5 TO TB5-6
  - TB5-6 TO TB5-7
18. APPLY TIE MARKER, ITEM 65, TO UNTERMINATED END OF CONDUIT, LABEL MARKER GROUND.
19. LABEL "OVERTEMP RESET" BELOW PUSHBUTTON ON OUTSIDE OF BOX, PER NOTE 8.
20. LABEL "JS4" BELOW CONNECTOR ON OUTSIDE OF BOX, PER NOTE 8.

| TOLERANCE STANDARD     | FRACTIONAL 1/164 | DECIMAL .000 |
|------------------------|------------------|--------------|
| UNLESS OTHERWISE NOTED | ANGULARITY       | MACH. SURF.  |
| 2103-18-92820052       |                  |              |
| 1112-06-919102111      |                  |              |
| 1112-06-919102111      |                  |              |
| 0111-04-919102111      |                  |              |

| QUANTITY | PART NUMBER | PART NAME           | DESCRIPTION, MATERIAL         | UNIT       |
|----------|-------------|---------------------|-------------------------------|------------|
| 1        | 146655-718  | INTERFACE BOX       | MEDALLION                     | 2 03-18-92 |
| 1        | 146655-718  | BOX MODIFICATION    |                               |            |
| 1        | 93910-110   | PACKET MODIFICATION |                               |            |
| 1        | 93821-001   | LINE FILTER         |                               |            |
| 1        | 146655-719  | HARNESS ASS'Y       | POWER SUPPLY                  |            |
| 1        | 44155-091   | GUARD BUTTON        | RESET                         |            |
| 1        | 146655-715  | PC BOARD ASS'Y      | DISPOSER INTERFACE            |            |
| 1        | 146655-875  | PC BOARD ASS'Y      | EAGLET II                     |            |
| 1        | 146655-720  | HARNESS ASS'Y       | P54-                          |            |
| 1        | 93902-900   | TERMINAL BOARD      | 7 POS. TB5                    |            |
| 1        | 93908-525   | FUSEBLOCK           | 3 POS. F2, F3, F4             |            |
| 2        | 150822-678  | FUSEBLOCK           | SINGLE F1, F5                 |            |
| 5        | 20341-091   | FUSE                | 0.5A 250V, F1, F2, F3, F4, F5 |            |
| 13       | 129361-013  | LED, SNAP-IN        |                               |            |
| 1        | 93909-993   | SWITCH, ROCKER      | AUTO/MAN SW1                  |            |
| 13       | 150778-001  | VARIABLE            |                               |            |
| 1        | 171         | VINYL INSL TUBE     | 1/16 X .48" LG. SEE NOTE #6   |            |
| 6        | 84123-001   | BAR, JUMPER         |                               |            |
| 19       |             |                     |                               |            |
| 4        | 129361-680  | STAND-OFF, M/F      | #6-32 X 7/16" LG.             |            |
| 4        | 129361-082  | STAND-OFF M/F       | #6-32 X 1 3/4" LG.            |            |
| 1        | 93910-097   | CABLE ASS'Y         | J2/P2                         |            |
| 1        | 136807-651  | CABLE ASS'Y         | CONTROL BOX                   |            |
| 1        | 136807-653  | CABLE ASS'Y         | ROTARY SOLENOID               |            |
| 1        | 93910-112   | LABEL               | LED/SWITCH ID SEE NOTE #5     |            |
| 1        | 129361-019  | PLUG, 6 POS.        | P1                            |            |
| 1        | 136807-917  | HARNESS ASS'Y       | P3                            |            |
| 28       |             |                     |                               |            |
| 29       |             |                     |                               |            |
| 1        | 129361-086  | PLUG, 12 POS.       | J1                            |            |
| 1        | 56396-486   | PLUG, 34 POS.       | P7                            |            |
| 34       |             |                     |                               |            |
| 1        |             | WIRE, BLACK         | 7'-0" LG. SEE NOTE #1         |            |
| 1        |             | WIRE, WHITE         | 0'-9" LG. SEE NOTE #2         |            |
| 1        |             | WIRE, BLACK         | 3'-9" LG. SEE NOTE #2         |            |
| 1        |             | WIRE, WHITE         | 0'-10" LG. SEE NOTE #3        |            |
| 1        |             | WIRE, BLACK         | 1'-2" LG. SEE NOTE #3         |            |
| 39       |             |                     |                               |            |
| 1        |             | WIRE, GREEN         | 6'-6" LG. SEE NOTE #3         |            |
| 31       | 150822-810  | CONTACT, SOCKET     | .025 SO POST 24/20 AWG        |            |
| 34       | 129359-024  | REC., FAST-ON       | .020 X .187 TAB 22/18 AWG     |            |
| 17       | 150476-410  | REC., FAST-ON       | .032 X .250 TAB 22/18 AWG     |            |
| 3        | 129359-454  | TERM. RING          | #18 AWG X #10 STUD            |            |
| 2        | 32118-091   | ADAPTER, OK DISC    |                               |            |
| 8        | 129326-001  | HARNESS MOUNT       |                               |            |
| 45       | 84104-001   | TIE WRAP            | SEE NOTE #10                  |            |
| 12       | 93908-039   | SCREW, SEMS         | #6-32 X 3/8" LG               |            |
| 12       | 93908-039   | SCREW, SEMS         | #6-32 X 1/4" LG               |            |
| 12       | 93908-039   | SCREW, SEMS         | #8-32 X 1/2" LG               |            |
| 5        | 129359-531  | NUT, KEPS           | #8-32                         |            |
| 4        | 93908-036   | SCREW, SEMS         | #10-32 X 3/8" LG              |            |
| 5        | 129359-532  | NUT, KEPS           | #10-32                        |            |
| 55       |             |                     |                               |            |
| 5        | 92615-001   | LABEL               | 1/2 AMP F1, F2, F3, F4, F5    |            |
| 3        | 50527-061   | SCREW, SOCKET       | #8-32 X 1/2" LG.              |            |
| 3        | 17796-091   | WASHER, FLAT        | .203 ID X .438 OD X .031 THK. |            |
| 1        | 93910-116   | PROGRAMMED EPROM    | UM                            |            |
| 1        | 129361-108  | CONTACT, RECEPTACLE | J7                            |            |
| 1        | 129361-109  | HOUSING, RECEPTACLE |                               |            |
| 63       |             |                     |                               |            |
| 12       | 84121-002   | NUT, KEPS           | #6-32                         |            |
| 2        | 150822-003  | WASHER, TIE         |                               |            |
| 69       |             |                     |                               |            |
| 1        | 129360-814  | SCREW, GREEN END    | #10-32 X 1/2" LG              |            |
| 1        | 129359-025  | LABEL               |                               |            |
| 69       |             |                     |                               |            |
| 70       |             |                     |                               |            |
| 71       |             |                     |                               |            |
| 1        | 18001-091   | CONN. ELL           | 3/8                           |            |
| 1        | 73608-091   | CONDUIT             | 3/8 X 52" LG                  |            |





AMERICAN STERILIZER COMPANY  
EPITE, PENNSYLVANIA  
TITLE: INTERFACE BOX ASS'Y.  
MEDALLION  
FIRST MADE FOR: 146655-109, 146655-115  
146655-718 146655-718














**Eagle 20 x 20  
Piping Schematic  
(Units With Printcon)  
With Disposer**





REF 146653-868, 146655-120  
 146655-119



NOTES:

1. USE PIPE SEALANT PER AMSCO DNG. 150822-020 ON ALL THREADED PIPE JOINTS.

2. COPPER TUBING .030 WALL X 1/4" OD. SOFT ANNEALED, ACTM B280.

3. INCLUDE WITH EACH KIT ONE PRINT, P/N 146655-119.


4. LABEL VACUUM SWITCH 'VS23', ITEM 6, WITH MARKER, BLACK LETTERS ON YELLOW BACKGROUND.

5. WIRE SPEC. #22 AWG STRANDED COPPER, 19/34 STRANDING, 300 VOLT, 105°C, TEFLON INSULATION, UL/CSA LABELS REQ'D, 48" LONG.

| QUANTITY | PART NUMBER | ITEM | PART NAME             | DESCRIPTION, MATERIAL          | REMARKS    |
|----------|-------------|------|-----------------------|--------------------------------|------------|
| 1        | 146655-119  | 1    | DISPOSER KIT          | 20" ROUND CHART 20X20          | 7 03-09-92 |
| 1        | 136806-983  | 2    | VALVE ASS'Y           | S41                            |            |
| 1        | 146655-121  | 3    | SCHEMATIC, ELEC.      | REF. DNG.                      |            |
| 1        | 136806-997  | 4    | BOX INTERFACE         |                                |            |
| 1        | 129361-107  | 5    | INSTALL. INST.        | REF. DNG.                      |            |
| 1        | 83866-001   | 6    | SWITCH, VAC           | VS23                           |            |
| 1        | 93909-980   | 7    | BRKT, VAC SW          |                                |            |
| 1        | 150822-295  | 8    | SUPPORT, PIPE         | 1/4" NPT X 3 1/2" LG.          |            |
| 1        | 93909-979   | 9    | BRKT, PIPE SPT.       |                                |            |
| 1        | 39590-010   | 10   | CLAMP, PIPE           | 1/4" NPT                       |            |
| 1        | 136807-667  | 11   | CABLE, AC             |                                |            |
| 1        | 136807-668  | 12   | CABLE, DC             |                                |            |
| 1        | 93910-135   | 13   | LIGHT BOX ASS'Y       |                                |            |
| 4        | 3985-041    | 14   | SCREW, RD. HD.        | #6-32 X 3/4" LG.               |            |
| 4        | 19675-041   | 15   | LOCKWASHER            | #6                             |            |
| 1        | 129084-001  | 16   | JUNCTION BOX          |                                |            |
| 2        | 91354-045   | 17   | CLAMP, CORD           |                                |            |
| 2        | 93829-001   | 18   | HOUSING, PLUG         | P1A, P2A                       |            |
| 1        | 19          |      |                       |                                |            |
| 1        | 20          |      |                       |                                |            |
| 4        | 129361-716  | 21   | SCREW, FLG LOCK       | 1/4-20 X 1" LG.                |            |
| 4        | 129360-015  | 22   | NUT, FLG LOCK         | 1/4-20                         |            |
| 1        | 4908-091    | 23   | TEE                   | 1/4" NPT                       |            |
| 8        | 8538-091    | 24   | NUT, WIRE             |                                |            |
| 1        | 129361-022  | 25   | PLUG, 10 POS          | P5                             |            |
| 1        | 129361-023  | 26   | PLUG, 21 POS          | P3                             |            |
| 1        | 93910-652   | 27   | ASSEMBLY              | PROGRAMMED CHIPS               |            |
| 1        | 93910-653   | 28   | TRAP, WATER           |                                |            |
| 1        | 22472-091   | 29   | RED. BUSHING          | 3/4 NPT X 1/4 NPT              |            |
| 1        | 45407-091   | 30   | ELL. COMP.            | 1/4 NPT X 1/4 CDT              |            |
| 2        | 34218-091   | 31   | FITTING, STR.         | 1/4 OD X 1/4" NPT              |            |
| 1        | 32          |      |                       |                                |            |
| 3        | 28917-091   | 33   | NIPPLE, BRASS         | 1/4" OD X 1/4" LG. SEE NOTE #2 |            |
| 2        | 150828-689  | 34   | NUT, LOCK             | 1/4" NPT X 1" LG.              |            |
| 2        | 52148-045   | 35   | LOCKWASHER            | 9/16" ID                       |            |
| 1        | 36          |      |                       |                                |            |
| 2        | 129361-115  | 37   | WALLTRAK              | REF PRINT P/N 146655-119       |            |
| 2        | 129361-116  | 38   | FLAT ELBOW            | 5" LG.                         |            |
| 2        | 129361-117  | 39   | SPLICE                | WALLTRAK                       |            |
| 1        | 40          |      |                       |                                |            |
| 1        | 41          |      |                       |                                |            |
| 2        | 150822-003  | 42   | TIE MARKER            | LABEL "P1A" AND "P2A"          |            |
| 1        | 150828-475  | 43   | VALVE, SAFETY         | 1/2" NPT                       |            |
| 1        | 129361-124  | 44   | EQUIPMENT MANUAL      | REF. DNG.                      |            |
| 1        | 93910-145   | 45   | WIRING CONNECTIONS    | DISPOSER - REF. DNG.           |            |
| 1        | 44500-091   | 46   | ELL. 90°              | SOLDER JOINT FITTING           |            |
| 1        | 84205-001   | 47   | HOUSING, PIN          | P3A                            |            |
| 2        | 84192-001   | 48   | JACKSCREW             | FEMALE, P3A                    |            |
| 2        | 3439-091    | 49   | PIPE PLUG             | 1/8 NPT                        |            |
| 1        | 146581-004  | 50   | MOTHERBOARD ASSEMBLY  |                                |            |
| 1        | 146590-058  | 51   | PC BOARD              |                                |            |
| 1        | 146645-226  | 52   | CONTROL, PRIMARY      |                                |            |
| 1        | 53          |      |                       |                                |            |
| 1        | 129352-535  | 54   | HOUSING               | #22 AWG BLACK, SEE NOTE #6     |            |
| 2        | 84187-003   | 55   | CONNECTOR, 3 POS. P13 |                                |            |
| 2        | 84207-002   | 56   | SOCKET                | CONTACT                        |            |
| 2        | 84207-002   | 57   | PIN, CONTACT          |                                |            |
| 1        | 83528-001   | 58   | I.D. LOCKING TAB      |                                |            |
| 2        | 84223-001   | 59   | CARD GUIDE            |                                |            |
| 2        | 84223-001   | 60   | LOCKING TAB           |                                |            |
| 1        | 150828-252  | 61   | VALVE, CHECK          |                                |            |
| 1        | 45585-091   | 62   | COUPLING              | 1/4 NPT                        |            |
| 1        | 939-091     | 63   | BUSHING, RED.         | 1/4 NPT X 1/8 NPT              |            |
| 1        | 76053-091   | 64   | BUSHING, RED.         | 1/2 NPT X 1/4 NPT              |            |
| 1        | 90525-091   | 65   | GAUGE                 | 1/8 NPT                        |            |
| 1        | 129356-508  | 66   | TEE, SOLDER           | 3/4 X 3/4 X 1/2 NPT            |            |
| 1        | 93910-642   | 67   | VALVE, BALL           | 3/4 SOLDER                     |            |
| 1        | 764317-524  | 68   | TAPE, TEFLON          |                                |            |
| 1        | 129362-980  | 69   | GUIDELINES, DISP.     |                                |            |
| 1        | 146657-256  | 70   | SCHEMATIC             |                                |            |
| 1        | 150822-048  | 71   | VALVE, BALL           |                                |            |


| TOLERANCE STANDARD                                 | FRACTIONAL 1/64 DECIMAL 0.005      |
|----------------------------------------------------|------------------------------------|
| UNLESS OTHERWISE NOTED                             | ANGULARITY 1/16" MACH. SURF. 0.002 |
| 3 102-12-90 RA 107 6 102-10-92 920034              |                                    |
| 2 9-06-89 RA 47 5 11-11-91 910211                  |                                    |
| 1 6-20-89 RA 29 4 10-08-90 900200                  |                                    |
| 3 6-7-89 RA 15 4 10-08-90 RA 170 7 03-09-92 920052 |                                    |

|                                                                 |                                                                                                                 |
|-----------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| <p>AMERICAN STEEL TIE COMPANY<br/>       ERIE, PENNSYLVANIA</p> | <p>TITLE ETO DISPOSER KIT ASS'Y<br/>       20X20 IRC EAGLE 2000</p> <p>FIRST MADE FOR<br/>       146655-119</p> |
|-----------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|



|      |            |     |
|------|------------|-----|
| REF. | 1          | 2   |
| NO.  | 146655-119 | 393 |

| 8        |              | 9        |                | 10                    |      | 11            |  |
|----------|--------------|----------|----------------|-----------------------|------|---------------|--|
| QUANTITY | PART NUMBER  | ITEM NO. | PART NAME      | DESCRIPTION, MATERIAL | REV. | REVISION DATE |  |
|          | X 146655-119 |          | DISPOSABLE KIT | 20" ROUND CHART 2000  | 7    | 03-09-92      |  |




7. LOCATE ITEMS 62, 63, 64, 65, & 66 AS CLOSE TO STERILIZER AS POSSIBLE.
8. LOCATE COMPOUND GAUGE, ITEM 64, BETWEEN BALL VALVE, ITEM 66, & DISPOSER.
9. BALL VALVE (ITEM 66) IS DESIGNED TO BE SOFT SOLDERED. APPLY HEAT WITH THE FLAME DIRECTED AWAY FROM THE CENTER OF VALVE BODY. SOLDER THE BALL VALVE TO THE TEFLON SEATS. SOLDER VALVES ONLY IN THE FULLY CLOSED POSITION.
10. MAKE DIM 'A' AS SHORT AS POSSIBLE TO REDUCE WATER FORMATION. SLOPE LINE TOWARDS STERILIZER.

[illegible][illegible]











ALF.  
DWG. NO. 136806-983  
SHEET 1

| QUANTITY | PART NUMBER | REV. NO. | PART NAME          | DESCRIPTION, MATERIAL   | REV. NO. | REVISION DATE |
|----------|-------------|----------|--------------------|-------------------------|----------|---------------|
| 1        | 136806-983  | 1        | VALVE ASSY         | 20" PRINTICON/IRC GAS   | 4        | 11-01-91      |
| 1        | 83228-002   | 2        | VALVE, SOL.        | 3/4" NPT                |          |               |
| 1        | 150828-844  | 3        | VALVE, CHECK       | 3/4" NPT                |          |               |
| 1        | 129360-765  | 4        | ORIFICE PL. CLMT.  | 3/4" NPT X 7/64" DIA ID |          |               |
| 1        | 4932-091    | 5        | TEE                | 3/4" NPT                |          |               |
| 1        | 84327-001   | 6        | TEE, REDUCING      | 3/4" X 3/4" X 1/4"      |          |               |
| 1        | 89986-001   | 7        | ELL. UNION         | 3/4" NPT                |          |               |
| 1        | 3443-091    | 8        | PLUG, PIPE         | 3/4" NPT                |          |               |
| 2        | 3439-091    | 9        | PLUG, PIPE         | 1/8" NPT                |          |               |
| 1        | 29293-091   | 10       | NIPPLE, BRASS      | 3/4" NPT X 1 3/4" LG.   |          |               |
| 5        | 29291-091   | 11       | NIPPLE, BRASS      | 3/4" NPT X 1 1/4" LG.   |          |               |
| 1        | 29307-091   | 12       | NIPPLE, BRASS      | 3/4" NPT X 5 1/4" LG.   |          |               |
| 1        | 28899-091   | 13       | NIPPLE, BRASS      | 1/8" NPT X 1" LG.       |          |               |
| 1        | 1345-091    | 14       | CROSS              | 3/4" NPT                |          |               |
| 1        | 22472-091   | 15       | BUSHING, RED.      | 3/4" NPT X 1/4" NPT     |          |               |
| 1        | 150822-934  | 16       | VALVE, SOLENOID    | 1/4" NPT                |          |               |
| 1        | 939-091     | 17       | BUSHING, RED.      | 1/4" NPT X 1/8" NPT     |          |               |
| 1        | 77774-091   | 18       | STRAINER, NIPPLE   | 1/2" NPT                |          |               |
| 1        | 47708-091   | 19       | STRAINER           | 3/4" NPT                |          |               |
| 1        | 34218-091   | 20       | FITTING, STR. COMP | 1/4" ODT X 1/4" NPT     |          |               |

NOTES:

1. USE PIPE SEALANT PER AMSCO DWG. 150822-020 ON ALL THREADED PIPE JOINTS.
2. LABEL VALVES AS SHOWN WITH BRADY TYPE MARKERS, BLACK LETTERS ON YELLOW BACKGROUND.






| TOLERANCE STANDARD     | FRACTIONAL 1/64 | DECIMAL ±.005    |
|------------------------|-----------------|------------------|
| UNLESS OTHERWISE NOTED | ANGULARITY ±1°  | MACH. SURF. ±25° |
| 3                      | 07-19-90        | 900129           |
| 2                      | 11-28-89        | RA #84           |
| 1                      | 06-15-89        | RA #24           |
| 0                      | 04-28-89        | PER #2           |
| 3                      | 07-19-90        | 900200           |



AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA  
THEY SPECIFY SPECIFICATIONS AND ALL INFORMATION THEREIN ARE THE PROPERTY OF AMERICAN STERILIZER COMPANY AND TO BE KEPT CONFIDENTIAL. NOT TO BE COPIED OR USED IN ANY MANNER WITHOUT THE WRITTEN PERMISSION OF AMERICAN STERILIZER COMPANY AND SHALL BE RETURNED UPON REQUEST.

TITLE VALVE ASSEMBLY - 20" GAS  
PRINT. OR IRC W/DISPOSER  
FIRST MADE FOR: 146653-863  
DWG. NO. 136806-983  
SHEET 1

REF 146655-119, 146655-118  
 146655-120



# NOTES:

1. USE PIPE SEALANT PER AMSCO DWG. 150822-020 ON ALL THREADED PIPE JOINTS.
2. COPPER TUBING .030 WALL X 1/4" OD. SOFT ANNEALED, ASTM B280.
3. INCLUDE WITH EACH KIT ONE PRINT, P/N 146655-120.
4. LABEL VACUUM SWITCHES 'VS10' AND 'VS23', ITEM 6, WITH MARKERS, BLACK LETTERS ON YELLOW BACKGROUND.
6. WIRE SPEC: #22 AWG STRANDED COPPER, 19/34 STRANDING, 300 VOLT, 105°C., TEFLON INSULATION, UL/CSA LABELS REQ'D, 48" LONG.

| QUANTITY | PART NUMBER | ITEM NO. | PART NAME            | DESCRIPTION, MATERIAL          | REMARKS  |
|----------|-------------|----------|----------------------|--------------------------------|----------|
| X        | 146655-120  | 1        | DISPOSER KIT         | 24 X 36 IRC EAGLE 2000         | CC-17-92 |
| 1        | 136807-661  | 2        | VALVE ASS'Y          | S41, S42, S28                  |          |
| 1        | 146655-121  | 3        | SCHEMATIC, ELEC.     | REF. DWG.                      |          |
| 1        | 136806-997  | 4        | BOX, INTERFACE       |                                |          |
| 1        | 129361-107  | 5        | INSTALL. INST.       | REF. DWG.                      |          |
| 2        | 83866-001   | 6        | SWITCH, VAC          | VS10, VS23                     |          |
| 2        | 93909-980   | 7        | BRKT, VAC SW         |                                |          |
| 1        | 150822-296  | 8        | SUPPORT, PIPE        | 1/4" NPT X 5 1/4" LG.          |          |
| 1        | 93910-143   | 9        | BRKT, PIPE SPT.      |                                |          |
| 1        | 93990-010   | 10       | CLAMP, PIPE          | 1/4" NPT                       |          |
| 1        | 136807-665  | 11       | CABLE ASS'Y, AC      |                                |          |
| 1        | 136807-666  | 12       | CABLE ASS'Y, DC      |                                |          |
| 1        | 93910-135   | 13       | LIGHT BOX ASS'Y      |                                |          |
| 8        | 3984-041    | 14       | SCREW, RD.HD.        | #6-32 X 3/8" LG.               |          |
| 8        | 19675-041   | 15       | LOCKWASHER           | #6                             |          |
| 1        | 93617-001   | 16       | JUNCTION BOX         |                                |          |
| 2        | 91354-045   | 17       | CLAMP, CORD          |                                |          |
| 2        | 93829-001   | 18       | HOUSING, PLUG        | P1A, P2A                       |          |
| 4        | 9374-041    | 19       | SCREW, RD.HD.        | #10-32 X 3/8" LG.              |          |
| 4        | 81682-006   | 20       | LOCKWASHER           | #10                            |          |
| 8        | 3949-042    | 21       | SCREW, RD.HD.        | 1/4-20 X 3/4" LG.              |          |
| 8        | 29357-533   | 22       | NUT, KEPS            | 1/4-20                         |          |
| 8        | 150473-295  | 23       | WASHER, FLAT         | .281 ID X .625 OD X .063 THK   |          |
| 12       | 18538-091   | 24       | NUT, WIRE            |                                |          |
| 1        | 129361-022  | 25       | PLUG, 10 POS         | P5                             |          |
| 1        | 129361-023  | 26       | PLUG, 21 POS         | P3                             |          |
| 1        | 93910-652   | 27       | ASSEMBLY             | PROGRAMMED CHIPS               |          |
| 1        | 33910-653   | 28       | TRAP, WATER          |                                |          |
| 1        | 28899-091   | 29       | NIPPLE, BRASS        | 1/8 NPT X 1"                   |          |
| 1        | 45407-091   | 30       | SOLEIL, COMP         | 1/4 NPT X 1/4 COT              |          |
| 2        | 34218-091   | 31       | FITTING, STR.        | 1/4 OD X 1/4" NPT              |          |
| 1        | 33319-091   | 32       | TUBING, COPPER       | 1/4" OD X 156" LG. SEE NOTE #2 |          |
| 2        | 150828-689  | 33       | NIPPLE, BRASS        | 3/4 NPT X 17 1/2" LG.          |          |
| 2        | 52148-045   | 34       | NUT, LOCK            | 1/4-18 NPSL                    |          |
| 2        | 52148-045   | 35       | LOCKWASHER           | 9/16" ID                       |          |
| 2        | 129361-115  | 36       | WALLTRAK             | REF PRINT P/N 146655-120       |          |
| 2        | 129361-116  | 37       | FLAT ELBOW           | WALLTRAK                       |          |
| 2        | 129361-117  | 38       | SPLICE               | WALLTRAK                       |          |
| 2        | 150822-003  | 39       | TIE MARKER           | LABEL P1A, P2A                 |          |
| 1        | 21950-091   | 40       | ELBOW, COMP          | 1/4 COT X 1/8" NPT             |          |
| 1        | 42581-091   | 41       | TEE, COMP            | 1/4 COT X 1/4" NPT             |          |
| 1        | 150828-475  | 42       | VALVE, SAFETY        | 1/2" NPT                       |          |
| 1        | 129361-124  | 43       | EQUIPMENT MANUAL     | REF. DWG.                      |          |
| 1        | 93910-145   | 44       | WIRING CONNECTIONS   | DISPOSER - REF. DWG.           |          |
| 1        | 836-042     | 45       | BUSHING, RED.        | 3/4 X 1/2                      |          |
| 1        | 44500-091   | 46       | ELL, 90°             | SOLDER JOINT FITTING           |          |
| 1        | 84205-001   | 47       | HOUSING, PIN         | P3A                            |          |
| 2        | 84192-001   | 48       | JACKSCREW            | FEMALE, P3A                    |          |
| 2        | 3419-091    | 49       | PIPE PLUG            | 1/8 NPT                        |          |
| 1        | 146581-004  | 50       | MOTHERBOARD ASSEMBLY |                                |          |
| 1        | 146590-058  | 51       | PC BOARD             |                                |          |
| 1        | 146645-226  | 52       | CONTROL, PRIMARY     |                                |          |
| 1        | 54          | 53       | WIRE                 | #22 AWG BLACK, SEE NOTE #6     |          |
| 1        | 129352-535  | 54       | HOUSING              | CONNECTOR, 3 POS, P13          |          |
| 2        | 84187-003   | 55       | SOCKET               | CONTACT                        |          |
| 2        | 84207-002   | 56       | PIN, CONTACT         |                                |          |
| 1        | 83528-001   | 57       | I.D. LOCKING TAB     |                                |          |
| 2        | 84225-001   | 58       | CARD GUIDE           |                                |          |
| 3        | 84226-001   | 59       | LOCKING TAB          |                                |          |
| 1        | 150828-262  | 60       | VALVE, CHECK         |                                |          |
| 1        | 45589-091   | 61       | COPPLING             | 1/4 NPT                        |          |
| 1        | 939-091     | 62       | BUSHING, RED.        | 1/4 NPT X 1/8 NPT              |          |
| 1        | 76053-091   | 63       | BUSHING, RED.        | 1/2 NPT X 1/4 NPT              |          |
| 1        | 90525-091   | 64       | GAUGE                | 1/8 NPT                        |          |
| 1        | 129356-508  | 65       | TEE, SOLDER          | 3/4 X 3/4 X 1/2 NPT            |          |
| 1        | 93910-642   | 66       | VALVE, BALL          | 3/4 SOLDER                     |          |
| 1        | 764317-524  | 67       | TAPE, TEFLON         |                                |          |
| 1        | 129362-980  | 68       | GUIDELINES, DISP.    |                                |          |
| 1        | 146657-256  | 69       | SCHEMATIC            |                                |          |
| 1        | 150822-048  | 70       | VALVE, BALL          |                                |          |
| 1        | 28917-091   | 71       | NIPPLE, BRASS        |                                |          |

| TOLERANCE STANDARD | FRACTIONAL 11/64 | DECIMAL .005 | ANGULARITY 11/64 | MACH. SURF. 11/64 |
|--------------------|------------------|--------------|------------------|-------------------|
| 1                  | 1-17-90          | PA 99        | 6                | 11-11-91          |
| 2                  | 9-05-89          | PA 48        | 5                | 10-08-90          |
| 3                  | 20-89            | PA 30        | 5                | 10-08-90          |
| 4                  | 6-7-89           | PER 19       | 4                | 02-12-90          |
| 5                  |                  |              |                  |                   |
| 6                  |                  |              |                  |                   |
| 7                  |                  |              |                  |                   |
| 8                  |                  |              |                  |                   |
| 9                  |                  |              |                  |                   |
| 10                 |                  |              |                  |                   |
| 11                 |                  |              |                  |                   |
| 12                 |                  |              |                  |                   |
| 13                 |                  |              |                  |                   |
| 14                 |                  |              |                  |                   |
| 15                 |                  |              |                  |                   |
| 16                 |                  |              |                  |                   |
| 17                 |                  |              |                  |                   |
| 18                 |                  |              |                  |                   |
| 19                 |                  |              |                  |                   |
| 20                 |                  |              |                  |                   |
| 21                 |                  |              |                  |                   |
| 22                 |                  |              |                  |                   |
| 23                 |                  |              |                  |                   |
| 24                 |                  |              |                  |                   |
| 25                 |                  |              |                  |                   |
| 26                 |                  |              |                  |                   |
| 27                 |                  |              |                  |                   |
| 28                 |                  |              |                  |                   |
| 29                 |                  |              |                  |                   |
| 30                 |                  |              |                  |                   |
| 31                 |                  |              |                  |                   |
| 32                 |                  |              |                  |                   |
| 33                 |                  |              |                  |                   |
| 34                 |                  |              |                  |                   |
| 35                 |                  |              |                  |                   |
| 36                 |                  |              |                  |                   |
| 37                 |                  |              |                  |                   |
| 38                 |                  |              |                  |                   |
| 39                 |                  |              |                  |                   |
| 40                 |                  |              |                  |                   |
| 41                 |                  |              |                  |                   |
| 42                 |                  |              |                  |                   |
| 43                 |                  |              |                  |                   |
| 44                 |                  |              |                  |                   |
| 45                 |                  |              |                  |                   |
| 46                 |                  |              |                  |                   |
| 47                 |                  |              |                  |                   |
| 48                 |                  |              |                  |                   |
| 49                 |                  |              |                  |                   |
| 50                 |                  |              |                  |                   |
| 51                 |                  |              |                  |                   |
| 52                 |                  |              |                  |                   |
| 53                 |                  |              |                  |                   |
| 54                 |                  |              |                  |                   |
| 55                 |                  |              |                  |                   |
| 56                 |                  |              |                  |                   |
| 57                 |                  |              |                  |                   |
| 58                 |                  |              |                  |                   |
| 59                 |                  |              |                  |                   |
| 60                 |                  |              |                  |                   |
| 61                 |                  |              |                  |                   |
| 62                 |                  |              |                  |                   |
| 63                 |                  |              |                  |                   |
| 64                 |                  |              |                  |                   |
| 65                 |                  |              |                  |                   |
| 66                 |                  |              |                  |                   |
| 67                 |                  |              |                  |                   |
| 68                 |                  |              |                  |                   |
| 69                 |                  |              |                  |                   |
| 70                 |                  |              |                  |                   |
| 71                 |                  |              |                  |                   |
| 72                 |                  |              |                  |                   |
| 73                 |                  |              |                  |                   |
| 74                 |                  |              |                  |                   |
| 75                 |                  |              |                  |                   |
| 76                 |                  |              |                  |                   |
| 77                 |                  |              |                  |                   |
| 78                 |                  |              |                  |                   |
| 79                 |                  |              |                  |                   |
| 80                 |                  |              |                  |                   |
| 81                 |                  |              |                  |                   |
| 82                 |                  |              |                  |                   |
| 83                 |                  |              |                  |                   |
| 84                 |                  |              |                  |                   |
| 85                 |                  |              |                  |                   |
| 86                 |                  |              |                  |                   |
| 87                 |                  |              |                  |                   |
| 88                 |                  |              |                  |                   |
| 89                 |                  |              |                  |                   |
| 90                 |                  |              |                  |                   |
| 91                 |                  |              |                  |                   |
| 92                 |                  |              |                  |                   |
| 93                 |                  |              |                  |                   |
| 94                 |                  |              |                  |                   |
| 95                 |                  |              |                  |                   |
| 96                 |                  |              |                  |                   |
| 97                 |                  |              |                  |                   |
| 98                 |                  |              |                  |                   |
| 99                 |                  |              |                  |                   |
| 100                |                  |              |                  |                   |



SCALE: 1/2:1


AMERICAN STERILIZER COMPANY  
 EST. 1924  
 146655-120

TITLE  
 ETO DISPOSER KIT ASSY  
 24X36 IRC EAGLE 2000


FIRST MADE FOR:  
 146655-120
















F- 1

F- 2

























|                        |                               |
|------------------------|-------------------------------|
| TOLERANCE STANDARD     | FRACTIONAL 1/64 DECIMAL .005  |
| UNLESS OTHERWISE NOTED | ANGULARITY .010 MACH SURF 125 |
| DATE                   | 03-18-92                      |
| REV                    | 0                             |
| QTY                    | 1                             |
| PART NUMBER            | 146657-256                    |
| PART NAME              | SCHEMATIC                     |
| DESCRIPTION, MATERIAL  | IRC INTERFACE                 |
| REV                    | 0                             |
| DATE                   | 03-18-92                      |

**AMERICAN STERILIZER COMPANY**  
 Erie, Pennsylvania

**IRC-DISPOSER SEQUENCER INTERFACE**

DATE: 03-18-92 BY: J. J. J. FOR: J. J. J.

QTY: 1 PART NO: 146657-256


CAD DWG.






SCALE  $\frac{3/4}{1}$

**VIEW C-C**  
**SCALE 1:1**





| QUANTITY | PART NUMBER | ITEM NO | PART NAME     | DESCRIPTION, MATERIAL | REV | DATE     |
|----------|-------------|---------|---------------|-----------------------|-----|----------|
|          | 146655-712  | 1       | PPG SCHEMATIC | 20 X 20 GAS           | 3   | 03-13-92 |

|                        |                                   |
|------------------------|-----------------------------------|
| TOLERANCE STANDARD     | FRACTIONAL: ±1/64 DECIMAL: ±0.005 |
| UNLESS OTHERWISE NOTED | ANGULARITY: ±1° MACH SURF: ±2°    |
| 1103-13-92020052       |                                   |
| 2103-11-91910055       |                                   |
| 1108-13-90RA#325       |                                   |
| 0106-19-90PER#66       |                                   |

AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA

**TITLE** PIPING SCHEMATIC  
20 X 20 GAS-EAGLE 3000

**FIGS** MADE FOR  
146655-712

G- 13

G- 14

CAD DWG



AMSCO

EO DISPOSER SYSTEM

INTERFACE: EAGLE 3000 & 7000, MEDALLION, AERATION CABINET


P-75-123-251

9/8/82

2 of 3







| QTY | PART NUMBER | ITEM | PART NAME     | DESCRIPTION     | MATERIAL | ST       | DATE |
|-----|-------------|------|---------------|-----------------|----------|----------|------|
| 1   | 146655-713  | 1    | PPG SCHEMATIC | MEDIUM 3000 GAS | 3        | 03-13-92 |      |


  


| TOLERANCE STANDARD | FRACTIONAL | DECIMAL | ANGULARITY | MACH SURF |
|--------------------|------------|---------|------------|-----------|
| 3103-1-92B2005     |            |         |            |           |
| 2103-1-91B1005     |            |         |            |           |
| 1103-09-90B1325    |            |         |            |           |
| 0103-19-90B1325    |            |         |            |           |

AMERICAN STERILIZER COMPANY  
CRIE, PENNSYLVANIA

**PIPING SCHEMATIC**  
MEDIUM GAS-EAGLE 3000


ESTD MADE FOR  
146655-713





A- 5

A- 6



A- 7

A- 8


REF 136806-983

DWG NO 136807-904

SHEET OF 1

## NOTE:

1. USE PIPE SEALANT PER AMSCO DWG. 150822-020 ON ALL THREADED PIPE JOINTS.
2. LABEL VALVES AS SHOWN WITH BRADY TYPE MARKERS, BLACK LETTERS ON YELLOW BACKGROUND.

TOLERANCE STANDARD FRACTIONAL  $\pm 1/64$  DECIMAL  $\pm .005$   
UNLESS OTHERWISE NOTED ANGULARITY  $\pm 1^\circ$  MACH. SURF  $\pm .25^\circ$ 

|   |          |        |   |          |        |
|---|----------|--------|---|----------|--------|
| 3 | 12-10-90 | RA#386 |   |          |        |
| 2 | 09-13-90 | RA#346 | 5 | 03-02-92 | 920049 |
| 1 | 06-27-90 | RA#285 | 4 | 11-18-91 | 910211 |
| 0 | 06-18-90 | PER#63 | 3 | 12-10-90 | 910043 |


© AMSCO 1989


A-9

AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA  
THIS COMPANY'S SPECIFICATIONS AND ALL INFORMATION THEREON ARE THE PROPERTY OF AMERICAN STERILIZER COMPANY, AND ARE TO BE KEPT CONFIDENTIAL, AND NOT TO BE COPIED OR USED IN ANY MANNER WITHOUT THE WRITTEN PERMISSION OF AMERICAN STERILIZER COMPANY.  
DWG NO 136807-904  
DATE 11-18-91  
DATE 6-15-90  
DATE 6-18-90  
DATE 6-18-90TITLE VALVE ASSY-20" GAS  
EAGLE 3000 W/DISPOSER  
FIRST MADE FOR 146655-710  
DWG NO 136807-904  
SHEET OF 1SCALE  $\frac{1}{2}$ 

A-10

CAD DWG.





A- 13

A- 14





REF 136807-661  
DWG NO 136807-905

2

3

4


5

6


7

NOTES:

1. USE PIPE SEALANT PER AMSCO DWG. 150822-020 ON ALL THREADED PIPE JOINTS.
2. LABEL VALVES AS SHOWN WITH BRADY TYPE MARKERS. BLACK LETTERS ON YELLOW BACKGROUND.



VIEW A-A






| TOLERANCE STANDARD     |          | FRACTIONAL $\pm 1/64$    | DECIMAL $\pm .005$    |
|------------------------|----------|--------------------------|-----------------------|
| UNLESS OTHERWISE NOTED |          | ANGULARITY $\pm 1^\circ$ | MACH. SURF $\pm .005$ |
| 2                      | 03-02-92 | 920049                   |                       |
| 1                      | 09-13-90 | 910043                   |                       |
| 1                      | 09-13-90 | RA#347                   |                       |
| 0                      | 06-18-90 | PER#63                   | 3 06-01-92 920081     |

|  |                                                                                                                                                                                                                                                                                                    |  |                                                              |  |
|--|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------------------------------------------------------------|--|
|  | <b>AMERICAN STERILIZER COMPANY</b><br>ERIE, PENNSYLVANIA                                                                                                                                                                                                                                           |  | <b>TITLE</b><br>VALVE ASS'Y 24 X 36<br>EAGLE 3000 W/DISPOSER |  |
|  | <small>FOR SERVICE, SPECIFICATIONS AND ALL INFORMATION THEREON ARE THE PROPERTY OF AMERICAN STERILIZER COMPANY. AND TO BE KEPT CONFIDENTIAL. AND NOT TO BE COPIED OR USED IN ANY MANNER WITHOUT THE WRITTEN PERMISSION OF AMERICAN STERILIZER COMPANY. AND SHALL BE RETURNED UPON REQUEST.</small> |  | <b>FIRST MADE FOR</b> 146655-711                             |  |
|  | <b>DATE</b> PD <b>DATE</b> PD <b>DATE</b> MEC <b>DATE</b> MFG <b>DATE</b> SDN                                                                                                                                                                                                                      |  | <b>DWG NO</b> 136807-905                                     |  |
|  | <b>DATE</b> 03-11-90 <b>DATE</b> 06-10-90 <b>DATE</b> 08-10-90 <b>DATE</b> 08-10-90                                                                                                                                                                                                                |  | <b>SHEET</b> 1 OF 1                                          |  |



B- 3

B- 4


CAD DWG.




NOTE:  
N.O.E. SCHEMATIC  
SEE DWG. NO. 1368



DETAIL 'A'  
PUMP MOTOR CONNECTIONS





|                                                                                                                  |                                 |
|------------------------------------------------------------------------------------------------------------------|---------------------------------|
| TOLERANCE STANDARD                                                                                               | FRACTIONAL 1/164 DECIMAL + .005 |
| UNLESS OTHERWISE NOTED                                                                                           | ANGULARITY .125 MACH SURF 125   |
| <div style="display: flex; justify-content: space-between;"> <span>03-12-92 020052</span> <span>10</span> </div> |                                 |

**AMERICAN STERILIZER COMPANY**  
ERIE, PENNSYLVANIA

THIS DRAWING IS THE PROPERTY OF AMERICAN STERILIZER COMPANY. IT IS TO BE USED FOR THE MANUFACTURE OF EQUIPMENT SPECIFIED HEREON. IT IS NOT TO BE REPRODUCED OR COPIED IN ANY MANNER WITHOUT THE WRITTEN PERMISSION OF AMERICAN STERILIZER COMPANY.

DATE: 10/27/92 BY: JMD/ELC/PLM/PLM  
DATE: 10/27/92 BY: JMD/ELC/PLM/PLM


**TITLE** EAGLE 3000 DISPOSER/SEQUENCER INTERFACE

**FIRST MADE FOR**


**DWG NO** 146657-275

**SHEET** 1















|                     |              |          |             |          |           |                        |
|---------------------|--------------|----------|-------------|----------|-----------|------------------------|
| REF.                | 2            | 3        | 4           | 5        | 6         | 7                      |
| DWG. NO. 136807-671 | SHEET 1 OF 1 | QUANTITY | PART NUMBER | ITEM NO. | PART NAME | DESCRIPTION, MATERIAL  |
|                     |              | X        | 136807-671  | 1        | SCHEMATIC |                        |
|                     |              |          |             |          |           | REVISION DATE 09-06-89 |



|                                     |                                |
|-------------------------------------|--------------------------------|
| TOLERANCE STANDARD                  | FRACTIONAL ±1/64 DECIMAL ±.005 |
| UNLESS OTHERWISE NOTED              | ANGULARITY ±1° MACH. SURF. 125 |
| 1 09-06-89 RA #51                   |                                |
| 0 6-20-89 PER 25                    |                                |
| NO. DATE E.C. NO. NO. DATE E.C. NO. |                                |
| DESIGNER                            | REVISIONS                      |



AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA

TITLE SCHEMATIC, MODIFIED  
AERATOR W/DISPOSER

FIRST MADE FOR:

DWG. NO. 136807-671

SHEET 1 OF 1


C- 5

C- 6

CAD DWG.







FOR FACTORY TEST  
NOT APPLICABLE  
FOR FIELD  
INSTALLED  
SEQUENCERS

INTERFACE  
BOARD

P6-4

| TOLERANCE STANDARD     |     | FRACTIONAL 11/64 |     | DECIMAL 0.005 |     |
|------------------------|-----|------------------|-----|---------------|-----|
| UNLESS OTHERWISE NOTED |     | ANGULARITY 11°   |     | MACH SURF 11° |     |
| 1                      | 2   | 3                | 4   | 5             | 6   |
| 7                      | 8   | 9                | 10  | 11            | 12  |
| 13                     | 14  | 15               | 16  | 17            | 18  |
| 19                     | 20  | 21               | 22  | 23            | 24  |
| 25                     | 26  | 27               | 28  | 29            | 30  |
| 31                     | 32  | 33               | 34  | 35            | 36  |
| 37                     | 38  | 39               | 40  | 41            | 42  |
| 43                     | 44  | 45               | 46  | 47            | 48  |
| 49                     | 50  | 51               | 52  | 53            | 54  |
| 55                     | 56  | 57               | 58  | 59            | 60  |
| 61                     | 62  | 63               | 64  | 65            | 66  |
| 67                     | 68  | 69               | 70  | 71            | 72  |
| 73                     | 74  | 75               | 76  | 77            | 78  |
| 79                     | 80  | 81               | 82  | 83            | 84  |
| 85                     | 86  | 87               | 88  | 89            | 90  |
| 91                     | 92  | 93               | 94  | 95            | 96  |
| 97                     | 98  | 99               | 100 | 101           | 102 |
| 103                    | 104 | 105              | 106 | 107           | 108 |
| 109                    | 110 | 111              | 112 | 113           | 114 |
| 115                    | 116 | 117              | 118 | 119           | 120 |
| 121                    | 122 | 123              | 124 | 125           | 126 |
| 127                    | 128 |                  |     |               |     |



AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA

FOR FACTORY TEST NOT APPLICABLE FOR FIELD INSTALLED SEQUENCERS

DATE: 03-26-92



|                |                                                  |
|----------------|--------------------------------------------------|
| TITLE          | SCHEMATIC, ELECTRICAL - EAGLE DISPOSER/SEQUENCER |
| FIRST MADE FOR | 146657-295                                       |
| DATE           | 03-26-92                                         |
| REV            | 3                                                |

C- 11


C- 12

CAD DWG

|          |            |   |            |          |             |          |                   |
|----------|------------|---|------------|----------|-------------|----------|-------------------|
| REF. NO. | 136807-669 | 2 | 3          | 4        | 5           | 6        | 7                 |
| DWG. NO. | 136807-669 | 1 | SHEET OF 2 | QUANTITY | PART NUMBER | ITEM NO. | PART NAME         |
|          |            |   |            | X        | 136807-669  | 1        | WIRING            |
|          |            |   |            |          |             |          | SEQUENCER CONTROL |
|          |            |   |            |          |             |          | REV. NO. 4        |
|          |            |   |            |          |             |          | DATE 11-01-91     |



# SEQUENCER CONNECTIONS TO STERILIZER




- REFERENCE SEQUENCER KIT. P/N 129361-130
- A - CORD CLAMP
- B - WIRE NUT
- NOTES:
1. WHEN A SEQUENCER CONTROL ASSEMBLY IS PART OF THE INSTALLATION, THE STERILIZER (S) (AND AERATOR(S), IF APPLICABLE) DO NOT CONNECT DIRECTLY TO THE DISPOSER IN THIS CASE. THE STERILIZERS CONNECT TO THE SEQUENCER ASSEMBLY, AND THE SEQUENCER ASSEMBLY CONNECTS TO THE DISPOSER.
  2. 24 VAC CONTROL CIRCUITRY IS UTILIZED ON ALL INTERCONNECTING WIRING. MULTICONDUCTOR 22 OR 24 AWG STRANDED CABLE IS RECOMMENDED. OBTAIN 10 CONDUCTOR CABLE FOR THE CONNECTION FROM THE DISPOSER TO THE SEQUENCER. 10 CONDUCTOR FOR CONNECTION FROM THE STERILIZERS TO THE SEQUENCER. AND 4 CONDUCTOR FOR THE CONNECTION OF AERATOR(S) TO THE SEQUENCER.
  3. THE SEQUENCER BOX IS NORMALLY MOUNTED IN THE ACCESS AREA OF THE STERILIZERS TO AID IN SERVICING. WALL-MOUNTED IN THE ACCESS AREA IS ACCEPTABLE.
  4. BEFORE WIRING THE STERILIZERS TO THE SEQUENCER BOX, THEY MUST BE IDENTIFIED #1, #2, #3 AS DESIRED IN ORDER TO BE PROPERLY WIRED TO THE SEQUENCER ASSEMBLY.

## POWER INPUT CONNECTIONS TO SEQUENCER CONTROL

|                        |                                  |
|------------------------|----------------------------------|
| TOLERANCE STANDARD     | FRACTIONAL : 1/64 DECIMAL : .005 |
| UNLESS OTHERWISE NOTED | ANGULARITY : 1° MACH. SURF. 125° |
| 3 09-04-90 RA#163      |                                  |
| 2 07-17-90 RA#149      |                                  |
| 1 09-06-89 RA #52      | 4 11-01-91 910211                |
| 0 6-20-89 PER 25       | 3 07-04-90 900200                |
| REV. NO.               | REV. NO.                         |
| DATE                   | DATE                             |


|                                                   |                                                  |
|---------------------------------------------------|--------------------------------------------------|
| AMERICAN STERILIZER COMPANY<br>ERIE, PENNSYLVANIA | TITLE<br>WIRING CONNECTIONS<br>SEQUENCER CONTROL |
| FIRST MADE FOR:                                   | DWG. NO. 136807-669                              |
| DATE 11-01-91                                     | DATE 11-01-91                                    |
| DATE 11-01-91                                     | DATE 11-01-91                                    |

















146653-874/875 11 2



| QUANTITY | PART NUMBER | PART NAME            | DESCRIPTION    | MATERIAL | DATE     |
|----------|-------------|----------------------|----------------|----------|----------|
| 1        | 146653-874  | CONTROL BD.          | 901 STERILIZER |          | 11-06-89 |
| 1        | 146653-875  | CONTROL BD.          | E10 DISPOSER   |          | 11-06-89 |
| 1        | 136806-595  | BASE BOARD           |                |          |          |
| 1        | 129360-191  | MICROPROCESSOR       | 411 80C31      |          |          |
| 1        | 129360-427  | AL. TANK 1/2" X 1/2" | 02 62555       |          |          |
| 1        | 129360-184  | AL. TANK 1/2" X 1/2" | 03 7361383     |          |          |
| 1        | 129360-227  | 700V SUPPLY SUP      | 05 651239      |          |          |
| 1        | 129360-549  | RAM AND CLOCK        | 06 4456102     |          |          |
| 1        | 93902-368   | 700V CONVERTER       | 07 460804      |          |          |
| 1        | 93902-368   | 700V CONVERTER       | 08 4456102     |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 09 741554      |          |          |
| 1        | 129046-001  | 700V SUPPLY          | 10 741554      |          |          |
| 2        | 129360-050  | 1/2" X 1/2" X 1/2"   | 11 741554      |          |          |
| 1        | 129360-075  | 1/2" X 1/2" X 1/2"   | 12 741554      |          |          |
| 4        | 129360-770  | 1/2" X 1/2" X 1/2"   | 13 741554      |          |          |
| 1        | 93902-368   | 700V CONVERTER       | 14 4456102     |          |          |
| 1        | 93902-368   | 700V CONVERTER       | 15 4456102     |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 16 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 17 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 18 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 19 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 20 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 21 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 22 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 23 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 24 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 25 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 26 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 27 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 28 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 29 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 30 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 31 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 32 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 33 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 34 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 35 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 36 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 37 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 38 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 39 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 40 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 41 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 42 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 43 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 44 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 45 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 46 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 47 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 48 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 49 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 50 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 51 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 52 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 53 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 54 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 55 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 56 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 57 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 58 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 59 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 60 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 61 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 62 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 63 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 64 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 65 741554      |          |          |
| 1        | 129360-551  | 700V CONVERTER       | 66 741554      |          |          |

NOTES  
 1. STAMP LATEST REVISION NUMBER USING WHITE LETTERING  
 2. TOP SEALING TAP SHOULD BE REMOVED AFTER WELDING OPERATION. (THIS NOTE SHOULD APPEAR ON ALL ASSEMBLIES CALLING OUT THIS PART.)  
 3. INSD SCHEMATIC FILE NAME: 146653-874


| TOLERANCE STANDARD     | FRACTIONAL: 1/64 DECIMAL: .005 |
|------------------------|--------------------------------|
| UNLESS OTHERWISE NOTED | ANGULARITY: MACH SURF: 12      |
| 2 11-06-89             | RA #67                         |
| 1 5-23-89              | RA #17                         |
| 0 3-15-89              | PER #6                         |

CONTROL BOARD  
EAGLE II

146653-874/875 11 2

146653-874/875 2 2

| QUANTITY | PART NUMBER    | PART NAME   | DESCRIPTION, MATERIAL | REV | DATE     |
|----------|----------------|-------------|-----------------------|-----|----------|
|          | 146653-874/875 | CONTROL BD. |                       | 2   | 11-06-89 |



| TOLERANCE STANDARD     | FRACTIONAL -11/64 | DECIMAL -0.005 |
|------------------------|-------------------|----------------|
| UNLESS OTHERWISE NOTED | ANGULARITY: 11    | MACH SURF: 25  |

| SPRNG CAPACITOR | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 | 31 | 32 | 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | 41 | 42 | 43 | 44 | 45 | 46 | 47 | 48 | 49 | 50 | 51 | 52 | 53 | 54 | 55 | 56 | 57 | 58 | 59 | 60 | 61 | 62 | 63 | 64 | 65 | 66 | 67 | 68 | 69 | 70 | 71 | 72 | 73 | 74 | 75 | 76 | 77 | 78 | 79 | 80 | 81 | 82 | 83 | 84 | 85 | 86 | 87 | 88 | 89 | 90 | 91 | 92 | 93 | 94 | 95 | 96 | 97 | 98 | 99 | 100 |
|-----------------|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|-----|
|-----------------|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|-----|

AMERICAN STERILIZED COMPANY  
COLE, PENNSYLVANIA

TITLE CONTROL BOARD  
EAGLET II

FIRST MADE FOR  
146653-874/875 2 2

E-1

E-2

CAD DWG.

THIS PAGE  
INTENTIONALLY LEFT  
BLANK

|                  |                |          |          |               |
|------------------|----------------|----------|----------|---------------|
| SHEET<br>OF<br>1 | DRAWING NUMBER | MATERIAL | REV. NO. | REVISION DATE |
|                  |                |          | 0        | 10-5-89       |


  

KIT, SERVICE REPAIR

GAST PUMP MODEL NUMBER  
IS 4 VCF-10-M400X

MFR  
GAST MFR CORP  
BENTON HARBOR, MICH.

AMSCO PART NUMBER  
FOR COMPLETE PUMP  
IS P-764323-248



| Ref. No. | Description      | Part Qty. | P.T. N2. |
|----------|------------------|-----------|----------|
| • 2      | Felt             | 2         | 8344A    |
| • 5      | Head Gasket      | 2         | AF520A   |
| • 6      | Valve, Outlet    | 2         | AF545    |
| • 8      | Valve, Inlet     | 2         | AF544    |
| • 9      | Gasket, Cylinder | 2         | AF521    |
| • 11     | Piston Ring      | 4         | AF541    |
| • 12     | Piston Seal      | 4         | AF540    |
| • 14     | Rider Ring       | 2         | AF595    |
| • 18     | Manifold Sleeve  | 2         | AF567A   |
|          | Service Kit      | 1         | K263     |

• Denotes parts included in service kit.

|     |         |               |           |  |  |  |  |
|-----|---------|---------------|-----------|--|--|--|--|
| NO. | DATE    | E.C.A. NUMBER | REVISIONS |  |  |  |  |
| 0   | 10-5-89 |               |           |  |  |  |  |

AMSCO

PRINTED IN U.S.A. © AMSCO 19

AMERICAN STERILIZER COMPANY

ERIE, PENNSYLVANIA

THIS DRAWING, SPECIFICATIONS, AND ALL INFORMATION THEREIN ARE THE PROPERTY OF AMERICAN STERILIZER COMPANY. ARE TO BE MAINTAINED CONFIDENTIAL. ARE NOT TO BE COPIED OR USED IN ANY WAY DETRIMENTAL TO THE INTEREST OF AMERICAN STERILIZER COMPANY AND SHALL BE RETURNED ON REQUEST.

TITLE: KIT, SERVICE REPAIR

FOR PUMP P-764323-248

TOLERANCES UNLESS OTHERWISE SPECIFIED: MACHINED SURF. 125  
FRACTIONAL DIMS. ± 1/64 DECIMAL DIMS. ± .005 ANGULARITY ± 1°

FIRST MADE FOR:

DWN. *RAM* CRD.  ENG. *RL* MFC.

DATE *10-89* DATE  DATE *10-5-89* DWT

DWC. *P-764323-249* NO. *4*

SHEET *1* OF *1*

| DWG. NO. | 755716-132 | 1 SHEET OF 2 | QUANTITY | PART NUMBER | ITEM NO. | PART NAME      | DESCRIPTION, MATERIAL | REV NO. | REVISION DATE |
|----------|------------|--------------|----------|-------------|----------|----------------|-----------------------|---------|---------------|
|          |            |              | X        | 755716-132  | 1        | KIT, DISPOSER  | WATER TRAP            | 0       | 6-9-92        |
|          |            |              | 1        | 150822-048  | 2        | VALVE, BALL    |                       |         |               |
|          |            |              | 1        | 93910-653   | 3        | TRAP, WATER    |                       |         |               |
|          |            |              | 1        | 22472-091   | 4        | BUSHING, RED.  | 3/4NPT X 1/4NPT       |         |               |
|          |            |              | 1        | 45407-091   | 5        | ELBOW, COMP.   | 1/4NPT X 1/4ODT       |         |               |
|          |            |              | 1        | 34218-091   | 6        | FITTING, STR.  | 1/4OD X 1/4"NPT       |         |               |
|          |            |              | 12       | R915-103    | 7        | TUBING, COPPER | 1/4OD X 1/4" LG       |         |               |
|          |            |              | 1        | 28917-091   | 8        | NIPPLE, BRASS  | 1/4NPT X 1" LG        |         |               |

THE PURPOSE OF THE WATER TRAP IS TO PREVENT WATER FROM ENTERING THE EXHAUST PIPING TO THE DISPOSER. EXCESS WATER IN THE EXHAUST PIPING CAN CAUSE DISPOSER PUMP FAILURE.

IF EXCESS AMOUNTS OF WATER ARE FOUND AFTER INSTALLATION OF THE TRAP, LOOK AT THE FOLLOWING AREAS:

1. FAULTY DRAIN LINE CHECK VALVE
2. FAULTY VACUUM PUMP WATER INLET SOLENOID VALVE
3. IMPROPER STEAM CONDITION SETTING
4. POOR QUALITY STEAM
5. IMPROPERLY LEVELED STERILIZER

THE WATER TRAP CONSISTS OF THE TRAP, MISCELLANEOUS PIPE FITTINGS, AND A SHUT OFF VALVE. THE WATER TRAP USES A FLOATING BALL THAT WILL CHOKE THE AIR FLOW TO THE DISPOSER PUMP WHEN THE TRAP FILLS WITH WATER. THE DISPOSER WILL THEN NO LONGER BE ABLE TO PULL VACUUM FROM THE STERILIZER AND AN ALARM ON THE STERILIZER WILL OCCUR (EITHER A DISPOSER FLOW ALARM OR A TOO LONG IN STEP ALARM).

THE 1/4" LINE ON THE BOTTOM OF THE TRAP ALLOWS THE TRAP TO SELF DRAIN SMALL AMOUNTS OF COLLECTED WATER DURING THE STERILIZER WATER DUMP PHASE. THE SHUTOFF VALVE AT THE BOTTOM OF THE TRAP IS TO BE LEFT OPEN DURING NORMAL OPERATION. IT IS TO BE CLOSED ONLY WHEN IT IS NECESSARY TO PERFORM A VACUUM LEAK TEST ON THE EXHAUST PIPING.

TOLERANCE STANDARD | FRACTIONAL  $\pm 1/64$  DECIMAL  $\pm .005$   
UNLESS OTHERWISE NOTED | ANGULARITY  $\pm 1^\circ$  MACH. SURF.  $\sqrt{25}$

|           |        |               |           |      |               |
|-----------|--------|---------------|-----------|------|---------------|
| 0         | 6-9-92 |               |           |      |               |
| NO.       | DATE   | E.C.A. NUMBER | NO.       | DATE | E.C.A. NUMBER |
| REVISIONS |        |               | REVISIONS |      |               |



AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA

THIS DRAWING, SPECIFICATIONS, AND ALL INFORMATION THEREIN ARE THE PROPERTY OF AMERICAN STERILIZER COMPANY. ARE TO BE MAINTAINED CONFIDENTIAL. ARE NOT TO BE COPIED OR USED IN ANY WAY DETRIMENTAL TO THE INTEREST OF AMERICAN STERILIZER COMPANY AND SHALL BE RETURNED UPON REQUEST.

AMSCO

DWN BMS CKD JB ENG MFG  
DATE 6-9-92 DATE 6-9-92 DATE DATE

TITLE  
KIT, DISPOSER  
WATER TRAP

FIRST MADE FOR:


DWG. NO. 755716-132

1 SHEET OF 2


DWG. NO. 755716-132

2 SHEET 2  
OF

| QUANTITY |  |  |  | PART NUMBER | ITEM NO. | PART NAME     | DESCRIPTION, MATERIAL | REV NO | REVISION DATE |
|----------|--|--|--|-------------|----------|---------------|-----------------------|--------|---------------|
|          |  |  |  | 755716-132  | 1        | KIT, DISPOSER | WATER TRAP            | 0      | 6-9-92        |



MEDIUM STERILIZER TYPICAL



SMALL STERILIZER TYPICAL

TOLERANCE STANDARD FRACTIONAL  $\pm 1/64$  DECIMAL  $\pm .005$   
UNLESS OTHERWISE NOTED ANGULARITY  $\pm 1^\circ$  MACH. SURF.  $\sqrt{25}$

|           |      |              |           |      |              |
|-----------|------|--------------|-----------|------|--------------|
|           |      |              |           |      |              |
|           |      |              |           |      |              |
|           |      |              |           |      |              |
|           |      |              |           |      |              |
|           |      |              |           |      |              |
| NO        | DATE | E C A NUMBER | NO        | DATE | E C A NUMBER |
| REVISIONS |      |              | REVISIONS |      |              |



AMERICAN STERILIZER COMPANY  
ERIE, PENNSYLVANIA

THIS DRAWING, SPECIFICATIONS, AND ALL INFORMATION THEREIN ARE THE PROPERTY OF AMERICAN STERILIZER COMPANY. ARE TO BE MAINTAINED CONFIDENTIAL. ARE NOT TO BE COPIED OR USED IN ANY WAY DETRIMENTAL TO THE INTEREST OF AMERICAN STERILIZER COMPANY AND SHALL BE RETURNED UPON REQUEST.

**AMSCO**

|             |      |      |      |     |
|-------------|------|------|------|-----|
| DWN         | BMS  | CKD  | ENG  | MFG |
| DATE 6-9-92 | DATE | DATE | DATE |     |

TITLE KIT, DISPOSER  
WATER TRAP

FIRST MADE FOR:

DWG. NO. 755716-132

2 SHEET OF 2




# INSTALLATION INSTRUCTIONS:

THE VS 23 SWITCH IS TO BE INSTALLED AT THE POINT WHERE THE 3/4" ODT EXHAUST LINE CONNECTS WITH THE STERILIZERS PIPING. SEE FIGURE


1. INSTALL VS23 SWITCH IN EXHAUST PIPING AS SHOWN. NOTE THAT SWITCH IS TO BE LOCATED AT THE POINT WHERE THE STERILIZERS EXHAUST LINE CONNECTS TO THE 3/4" ODT COPPER TUBING.

SWITCH MUST BE INSTALLED AS VERTICAL AS POSSIBLE.

2. INSTALL CABLE AT VS23 AND P53 PLUG AS SHOWN.
3. STERILIZER MUST HAVE REV 11 PROGRAM CHIP OR LATER.
4. ADJUST VS23 SO THAT SWITCH ACTUATES WHEN CHAMBER IS AT 23" HG.




E- 9



VIEW A-A  
VACUUM SWITCH

## BILL OF MAT'L

| ITEM | QTY | PART NO.   | PART NAME                | DESCRIPTION      |
|------|-----|------------|--------------------------|------------------|
| 1    | 1   | 84323-001  | REDUCING TEE             | 3/4 X 3/4 X 1/4  |
| 2    | 1   | 939-091    | BUSHING                  | 1/4 X 1/8        |
| 3    | 1   | 29291-091  | NIPPLE                   | 3/4 X 1" 1/4 LG. |
| 4    | 1   | 28899-091  | NIPPLE                   | 1/8 X 1" LG.     |
| 5    | 1   | 83866-001  | VACCUM SWITCH            |                  |
| 6    | 1   | 93910-654  | CABLE                    |                  |
| 7    | 1   | 93909-980  | BRACKET                  |                  |
| 8    | 1   | 76262-091  | SNAP RING                |                  |
| 9    | 2   | 3984-041   | SCREW                    | *6-32 X 3/8" LG. |
| 10   | 2   | 19675-041  | LOCKWASHER               |                  |
| 11   | 4   | 129326-001 | HARNESS MOUNT            |                  |
| 12   | 6   | 84104-001  | WIRE TIE                 |                  |
| 13   | W-  | 755716-164 | INSTALLATION INSTRUCTION |                  |




VIEW B-B

P53  
SHOWN ENLARGED FOR CLARITY


KIT, VS23 INSTALLATION

P/N 755716-164  
REV 0 DATE 07/23/92  
PAGE 1 OF 2

E- 10



24 X 36



20 X 20

KIT. VS23 INSTALLATION

P/N 755716-164  
REV 0 DATE 07/23/92  
PAGE 2 OF 2



AMSCO

## EAGLE SERIES STAGE 2

INTERFACE: EAGLE 3000 & 2000, MEDALLION, AERATION CABINET  
P-764323-254

9/92

3 of 3

