

Introducing...

Steri-Pak™ LF

BOWIE-DICK TEST SHEETS

Large size Bowie-Dick Test Sheets with new blue to black ink technology

- **IMPROVED INK CHEMISTRY AIDS INTERPRETATION**

Vivid blue color remains if air leak is present.

- **HIGH STANDARDS**

Improved chemistry and improved pack design meets EN and BS standards for performance. ISO compliant.


- **PERMANENT COLOR CHANGE**

Color change from blue to black is not reversible. Indicator sheet will not fade in sunlight. Can be used as permanent record.

- **ENVIRONMENTALLY SAFE**

Improved chemistry is "Lead-Free" and Heavy Metals Free.

FAILURE

Site	Operator	Date	Department	CAT. NO. BD110
Steri-Pak™ LF LEAD FREE BOWIE-DICK TEST SHEET				
				
Supervisor	Machine No.	Cycle No.	Result	LOT 990201 134°C-3½ Min.

SIZE= 280 CM X 215 CM

STERITEC...Over 30 years of sterilization monitoring experience.