

Introducing...

The GREEN CARD



Bowie-Dick Test

Why do we call it the "Green Card"?

- A color change to **GREEN** indicates a Bowie-Dick Test PASS. If any purple remains, a FAIL is indicated.
- It is ecologically smart to save the forests. No disposable paper is used and this also saves precious storage space.
- Special holder is **GREEN**.
- It will save you a lot of **GREEN** due to its small size.
- Our competitors will be **GREEN** with envy.
- Meets AAMI Standard / International version meets EN 867-4 Standard
- Not useful for immigration purposes.

From the people who gave you the first small disposable Bowie-Dick Pack (See US Patent No. 4,579,715 Sept., 1983).

CARD



ACTUAL SIZE