

Introducing...

Steri-Pak™ LF

Disposable Bowie-Dick Test Pack with new blue to black ink technology

- **IMPROVED INK CHEMISTRY AIDS INTERPRETATION**

Vivid blue color remains if air leak is present.

- **HIGH STANDARDS**

Improved chemistry and improved pack design meets EN Standards (BD 108 S-T) and U.S. Standards (BD 111). ISO compliant.

- **CONTINUOUS INK DESIGN FROM CENTER TO EDGE**

Make air leaks easier to see.

- **PERMANENT COLOR CHANGE**


Can be used as permanent record.

- **ENVIRONMENTALLY SAFE**

Improved chemistry is "Lead-Free" and Toxic Heavy Metals Free.

- **EXTERNAL INDICATOR**

Verifies processing.



FAILURE

Site	Operator	Date	Department	
Steri-Pak™ LF LEAD FREE BOWIE-DICK TEST				
Supervisor	Machine No.	Cycle No.	Result	LOT 000000 134°C-3½ Min.