


Sterilizer Main Components

SECTION VI PARTS LIST



MA242701

Used on Units with Serial Numbers:
CP1000 thru CP1061 CR1000 thru CR1154 CS1000 thru CS1975
RB1000 thru RB7986 DR1000 thru DR1014

Item	Part No.	Description	Qty.	Item	Part No.	Description	Qty.
1		Panel (Refer to "Front Panel Components" Elsewhere)	Ref	9	MZY100027	Nut	1
2	MZY100019	Screw	1	10	MZY100028	Chamber Leg Bracket	2
3	MZY100020	Door Stop	1	11	MZY100029	Screw	1
4		Door Assembly (Refer to "Door Assembly" Elsewhere)	Ref	12	MZY100030	Tank Cover	1
5	MZY100021	Torque Nut	1	13		Tank Assembly (Refer to "Tank Assembly" Elsewhere)	Ref
6	MZY100022	Screw	7	14	MZY100031	Clip	4
7		Base Assembly (Refer to "Base Components" Elsewhere)	Ref	15		Vent Valve Components (Refer to "Valve Assembly" Elsewhere)	Ref
8	MZY100023	Shell Assembly [M7-001 & M7-004] (Refer to "Chamber Components {Side or Bottom}" Elsewhere)	1	16	MZY100032	Lever	1
	MZY100024	Shell Assembly [M7-005] (Refer to "Chamber Components {Side or Bottom}" Elsewhere)	1	17	MZY100021	Nut	1
	MZY100025	Shell Assembly [M7-002] (Refer to "Chamber Components {Side or Bottom}" Elsewhere)	1	18	MZY100034	Screw	1
	MZY100026	Shell Assembly [M7-003] (Refer to "Chamber Components {Side or Bottom}" Elsewhere)	1	19	MZY100035	Insulation	1
				20	MZY100036	Banding Strap (shown)	2
					MZY100037	Banding Strap End	1
				21	MZY100038	Banding Strap Clip	2
				22	MZY100039	Tubing	1
				23	MZY100040	Tubing	1
				24	MZY100041	Screw	4

Always Specify Model & Serial Number