


Sterilizer Main Components

SECTION VI PARTS LIST



MA242702

Used on Units with Serial Numbers:RB17993, CR1686, DR1268, CP1150, CS7754, FM1426 and FN11230 thru Present

Item	Part No.	Description	Qty.	Item	Part No.	Description	Qty.
1		Panel (Refer to "Front Panel Components" Elsewhere)	Ref	9	MZY100027	Nut	1
2	MZY100019	Screw	1	10	MZY100033	Chamber Leg Bracket	1
3	MZY100020	Door Stop	1	11	MZY100029	Screw	1
4		Door Assembly (Refer to "Door Assembly" Elsewhere)	Ref	12	MZY100030	Tank Cover	1
5	MZY100021	Torque Nut	1	13		Tank Assembly (Refer to "Tank Assembly" Elsewhere)	Ref
6	MZY100007	Screw	5	14	MZY100046	Screw	4
7		Base Assembly (Refer to "Base Components" Elsewhere)	Ref	15		Vent Valve Components (Refer to "Valve Assembly" Elsewhere)	Ref
8	MZY100023	Shell Assm. [M7-001,M7-004 & M7-005] (Refer to "Chamber Components (Side or Lower)" Elsewhere)	1	16	MZY100032	Lever	1
	MZY100043	Shell Assembly [M7-002 & M7-006] (Refer to "Chamber Components (Side or Lower)" Elsewhere)	1	17	MZY100021	Nut	1
	MZY100044	Shell Assembly [M7-003] (Refer to "Chamber Components (Side or Lower)" Elsewhere)	1	18	MZY100034	Screw	1
				19	MZY100035	Insulation	1
				20	MZY100036	Banding Strap (shown)	2
					MZY100037	Banding Strap End	1
				21	MZY100038	Banding Strap Clip	2
				22	MZY100039	Tubing	1
				23	MZY100033	Tubing	1

Always Specify Model & Serial Number