


Front Panel Components

SECTION VI PARTS LIST



MA241102

**Used on Units with Serial Numbers CP1000 thru CP1069,
CR1000 thru CR1125, CS1000 thru CS1694, RB1000 thru RB7499**

Item	Part No.	Description	Qty.	Item	Part No.	Description	Qty.
1	MZY100047 MZY100048	Ritter Panel (M7-001) Midmark Panel (M7-002, M7-003, M7-004)	1 1	8	MZY100054	Timer Knob	1
2	MZY100049	Flexible Shaft Assembly (Includes Items 3 and 4)	1	9		Nut	Ref
3	MZY100050	• Lockwasher	1	10	MZY100055 MZY100056	Pilot Light (M7-001, M7-003, M7-004) Pilot Light (M7-002)	1 1
4	MZY100051	• Nut	1	11	MZY100057	Fill / Vent Knob Assembly (Includes Item 12)	1
5	MZY100052	Thermometer	1	12	MZY100058	• Setscrew	1
6	MZY100053	Control Knob (Includes Item 7)	1	13	MZY100059 MZY100060	Timer (M7-001, M7-003, M7-004) Timer (M7-002)	1 1
7	•	• Setscrew	1				

Always Specify Model & Serial Number