


Chamber Components - Lower

SECTION VI PARTS LIST



MA241002

**Used on Units with Serial Numbers RB14565, CP1092, DR1015
CR1260, CS5297, FM1000 and FN1000 thru Present**

Item	Part No.	Description	Qty.	Item	Part No.	Description	Qty.
1		Chamber (Refer to "Chamber Components - Side" Elsewhere)	Ref	10	MZY100021	Nut	1
2	MZY100139	Heater (M7-001, M7-004, M7-005)	1	11	MZY100149	Spacer	1
	MZY100161	Heater (M7-002, M7-006)	1		MZY100151	Rocker Assembly (Includes Item 6, Item 15, Item 14)	1
	MZY100162	[Includes Items 3]]	1	12	MZY100152	Screw	1
3	MZY100142	Heater (M7-003) [Includes Items 3]	1	13	MZY100163	Rocker Bracket	1
4	MZY100143	Nut	2	14	MZY100153	Locknut	2
5	MZY100144	Gasket	2	15	MZY100041	Screw	1
6	MZY100166	Thermostat (Primary)	2	16	MZY100155	Bracket	1
7	MZY100146	Bracket	1	17	MZY100164	Reset Button	1
8	MZY100147	Washer	1	18	MZY100165	Rod	1
9	MZY100148	Lockwasher	2	19		Front Panel (Refer to "Front Panel Components" Elsewhere)	Ref
		Spacer	1	20	MZY100167	Thermostat (Secondary)(M7-001 & 004)	1

Always Specify Model & Serial Number