


Tank Assembly

SECTION VI PARTS LIST



MA240804

Used on Units with Serial Numbers:
CP1000 thru CP1061 CR1000 thru CR1113 CS1000 thru CS1569
RB1000 thru RB7193

Item	Part No.	Description	Qty.	Item	Part No.	Description	Qty.
	MZY100168	Tank Assembly (Incls. Items 1 thru 19) .	1	12	MZY100178	• Compression Nut	1
1	MZY100169	• Tank	1	13	MZY100175	• Male Elbow	1
2	MZY100170	• Safety Valve	1	14	MZY100179	• Connector Drilled	1
		(Apply Teflon Tape #014-0026-00)	1	15	•	• Compression Nut	1
3	MZY100171	• Pipe Nipple	1	16	MZY100180	Bracket	1
4	MZY100172	• Pipe "T"	1	17	MZY100181	Thermometer Nut	1
5	MZY100173	• Bellows Assembly	1	18	MZY100039	Tube	Ref
6	MZY100174	• Bellows	1	19		Tube (Refer to "Sterilizer Main Components" Item 23)	6-5
7	MZY100175	• Elbow	1	20	MZY100182	Drain (Includes O-Ring)	1
8	MZY100136	• Compression Nut	1		MZY100183	• O-ring	1
9	MZY100176	• Tube	1	21	MZY100052	Thermometer (Includes Item 22)	1
10	MZY100177	• Straight "T"	1	22	MZY100184	• Thermometer Sleeve	1
11	MZY100137	• Male Elbow	1				

Always Specify Model & Serial Number